

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-006689

Location: Emmons County
S.O. (LINT)

A. General Setup and Checks:

1. Diagnostic Tests Pass and Instrument in Ready Mode
2. Breath tube heated
3. Date, time and location code (Level 2,E). Re-set if necessary.
4. Print test (Level 1,P). Sign and attach test record.
5. Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
 - a. Display: 911 psi Regulator: 900 psi

B. Tests (Sign and attach test records):

- Configure simulator for the following tests (Level 1,S).
1. Wet Calibration Check
 - a. Low AC (Level 1,C): Known Value 17160 AC: 0.020 AC
Sim. Ser #: DR3375 Lot #: 717609 Exp. Date: 6.8.19
 - b. High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC
Sim. Ser #: DR7111 Lot #: 17350 Exp. Date: 10.11.19

Configure dry gas standard for the following tests (Level 1,S).

2. Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone
Sim. Ser #: D7089 Lot #: ICS 5 Exp. Date: NA
 Display reads "Interferent Detect"
3. RFI Check (CMS Mode) Display reads "RFI Detect"
4. Dry Calibration Check (Level 1,C): Known Value 0.080 AC
Gas Cyl Lot #: 19817080A4 Cyl #: 1 Exp. Date: 9/5/2019

Test 1 <u>0.083</u>	Test 4 <u>0.082</u>	Test 7 <u>0.083</u>
Test 2 <u>0.083</u>	Test 5 <u>0.083</u>	Test 8 <u>0.083</u>
Test 3 <u>0.083</u>	Test 6 <u>0.083</u>	Test 9 <u>0.083</u>
Average <u>0.083</u>		

C. Remarks/Maintenance: ANNUAL INSPECTION, REPLACED
INTERNAL PRINTER ROLL

Instrument is acceptable to be used in the field.

Charles E. Ed
Breath Analyst Signature
NA

4.29.2019
Date
NA

Reviewed by _____

Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006689
Location = LINT 8164.14.00 09/16
04/29/2019 13:30

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMNPO

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNOPQRSTU

24 characters per line
abcdefghijklmnopqrstuvw
ABCDEFGHIJKLMNOPQRSTUWX

32 characters per line
abcdefghijklmnopqrstuvwxy123456
ABCDEFGHIJKLMNOPQRSTUWXYZ123456

42 characters per line
abcdefghijklmnopqrstuvwxy1234567890abcdef
ABCDEFGHIJKLMNOPQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup
Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On MemFull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 19817080A4
Standard Cyl #? 1
Standard Expiration? 09/05/2019
Oper No? 666666

Flow Cal. Date: 09/19/2017
Slope 625
Intercept -591019

IR Calibration Date: 09/19/2017
3um 9um

0th Coef(*100): -14072 -18220
1st Coef(*100): 273930 133081
2nd Coef(*100): 2533 1211
H2O adj(mg/l)*10K): 588 599

***** Test End *****



Operator Signature
CHARLES EDER

Remarks: PRINT
TEST

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006689
Location = LINT 8164.14.00 09/16
04/29/2019 13:48

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:48
02 Std. Sol.	0.021	13:49
03 Room Air	0.000	13:49
04 Std. Sol.	0.022	13:50
05 Room Air	0.000	13:51
06 Std. Sol.	0.022	13:51
07 Room Air	0.000	13:52

08 Sim Temp = 34.0°C

Simul Ser No = DR3375
Std Sol No = 17160
County = 15 Oper No. = 666666


Operator Signature

CHARLES EDER

Remarks: *Low AC*
0.020 AC

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006689
Location = LINT 8164.14.00 09/16
04/29/2019 13:54

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:55
02 Std. Sol.	0.302	13:56
03 Room Air	0.000	13:56
04 Std. Sol.	0.302	13:57
05 Room Air	0.000	13:57
06 Std. Sol.	0.304	13:58
07 Room Air	0.000	13:59

08 Sim Temp = 34.0°C

Simul Ser No = DR7111
Std Sol No = 17350
County = 15 Oper No. = 666666


Operator Signature

CHARLES EDER

Remarks: *HIGH AC*
0.300 AC

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006689
Location = LINT 8164.14.00 09/16
04/29/2019 13:39

Test	AC	Time
01 Room Air	0.000	13:40
02 *Subject Test	INT*	13:41
03 Room Air	0.000	13:41

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/D15821456
Lot No = 19817080A4
Cyl No = 1
Expiration Date = 09/05/2019
County = 15 Oper No. = 666666


Operator Signature
CHARLES EDER

Remarks: *Interferent*
Check
0.10 WAC + 0.05% ~~X~~ CEE
Form 106-18000
ACETONE

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006689
Location = LINT 8164.14.00 09/16
04/29/2019 13:36

Test	AC	Time
01 Diagnostic	OK	13:37
02 Room Air	0.000	13:37
03 *Subject Test	RFI*	13:38
04 Room Air	0.000	13:38

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/D15821456
Lot No = 19817080A4
Cyl No = 1
Expiration Date = 09/05/2019
County = 15 Oper No. = 666666

I followed the Approved Method
and the instructions displayed
by the Intoxilyzer in conducting
this test.


Operator Signature
CHARLES EDER

Remarks: *RFI CHECK*

Form 106-18000


Intoxilyzer Test Record and Checklist
 NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
 North Dakota Model 8000 SN 80-006689
 Location = LINT 8164.14.00 09/16
 04/29/2019 13:31

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:31
02 Std. Gas	0.083	13:32
03 Room Air	0.000	13:32
04 Std. Gas	0.083	13:33
05 Room Air	0.000	13:33
06 Std. Gas	0.083	13:34
07 Room Air	0.000	13:34

Lot No = 19817080A4
 Cyl No = 1
 Exp Date = 09/05/2019
 County = 15 Oper No. = 666666


 Operator Signature
 CHARLES EDER

Remarks: CALIBRATION
 CHECK

Form 106-I8000
 0.080 AC


Intoxilyzer Test Record and Checklist
 NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
 North Dakota Model 8000 SN 80-006689
 Location = LINT 8164.14.00 09/16
 04/29/2019 13:42

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:43
02 Std. Gas	0.082	13:43
03 Room Air	0.000	13:44
04 Std. Gas	0.083	13:44
05 Room Air	0.000	13:45
06 Std. Gas	0.083	13:45
07 Room Air	0.000	13:46

Lot No = 19817080A4
 Cyl No = 1
 Exp Date = 09/05/2019
 County = 15 Oper No. = 666666


 Operator Signature
 CHARLES EDER

Remarks: CALIBRATION
 CHECK
 0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
 NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
 North Dakota Model 8000 SN 80-006689
 Location = LINT 8164.14.00 09/16
 04/29/2019 14:01

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:02
02 Std. Gas	0.083	14:02
03 Room Air	0.000	14:03
04 Std. Gas	0.083	14:03
05 Room Air	0.000	14:04
06 Std. Gas	0.083	14:04
07 Room Air	0.000	14:04

Lot No = 19817080A4
 Cyl No = 1
 Exp Date = 09/05/2019
 County = 15 Oper No. = 666666


 Operator Signature
 CHARLES EDER

Remarks: CALIBRATION
 CHECK
 0.080 AC

Form 106-I8000