

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-006679

Location: VALLEY CITY P.D.  
(VALL)

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 860 psi Regulator: 850 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check
  - a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.020 AC  
Sim. Ser #: DR5131 Lot #: 201810D Exp. Date: 10.24.20
  - b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.300 AC  
Sim. Ser #: DR7111 Lot #: 17350 Exp. Date: 10.11.19

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR7089 Lot #: JCS 5 Exp. Date: NA  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 19817080A4 Cyl #: 6 Exp. Date: 9.5.19

Test 1	<u>0.080</u>	Test 4	<u>0.081</u>	Test 7	<u>0.080</u>
Test 2	<u>0.081</u>	Test 5	<u>0.080</u>	Test 8	<u>0.080</u>
Test 3	<u>0.079</u>	Test 6	<u>0.081</u>	Test 9	<u>0.080</u>
Average	<u>0.080</u>				

C. Remarks/Maintenance: ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

Charles E. Eden  
Breath Analyst Signature

5.15.2019  
Date

NA  
Reviewed by

NA  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000       SN 80-006679  
Location = VALL                8164.14.00 09/16  
05/15/2019                       19:16

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?         No  
Display Volume?             No  
Disable On Memfull?         Yes  
# of Print Copies?           1  
Select Std (D/W/I)?         Dry  
Standard Value?             0.080  
Standard Lot #?             19817080A4  
Standard Cyl #?             6  
Standard Expiration?        09/05/2019  
Oper No?                    666666

Flow Cal. Date:             09/06/2016  
Slope                        701  
Intercept                   -682926

IR Calibration Date:        09/06/2016  
                              3um                    9um

-----  
0th Coef(\*100):            -18983            -23165  
1st Coef(\*100):            280986            134481  
2nd Coef(\*100):            2443              1290  
H2O adj(mg/l\*10k):        569                553

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006679  
Location = VALL      8164.14.00 09/16  
05/15/2019      19:31

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:31
02 Std. Sol.	0.018	19:32
03 Room Air	0.000	19:33
04 Std. Sol.	0.019	19:33
05 Room Air	0.000	19:34
06 Std. Sol.	0.019	19:35
07 Room Air	0.000	19:35

08 Sim Temp = 34.0°C

Simul Ser No = DR5131  
Std Sol No = 201810D  
County = 02      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

LOW AC  
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006679  
Location = VALL      8164.14.00 09/16  
05/15/2019      19:40

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:41
02 Std. Sol.	0.293	19:42
03 Room Air	0.000	19:42
04 Std. Sol.	0.294	19:43
05 Room Air	0.000	19:44
06 Std. Sol.	0.295	19:44
07 Room Air	0.000	19:45

08 Sim Temp = 34.0°C

Simul Ser No = DR7111  
Std Sol No = 17350  
County = 02      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.300 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006679  
Location = VALL      8164.14.00 09/16  
05/15/2019      19:21

Test	AC	Time
01 Room Air	0.000	19:22
02 *Subject Test	INT*	19:22
03 Room Air	0.000	19:23

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 19817080A4  
Cyl No = 6  
Expiration Date = 09/05/2019  
County = 02      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10 AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006679  
Location = VALL                    8164.14.00 09/16  
05/15/2019                            19:24

Test	AC	Time
01 Diagnostic	OK	19:25
02 Room Air	0.000	19:25
03 *Subject Test	RFI*	19:25
04 Room Air	0.000	19:26

\*Invalid Test  
  Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 19817080A4  
Cyl No = 6  
Expiration Date = 09/05/2019  
County = 02                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: **RFI CHECK**

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006679  
Location = VALL      8164.14.00 09/16  
05/15/2019      19:16

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:17
02 Std. Gas	0.080	19:18
03 Room Air	0.000	19:18
04 Std. Gas	0.081	19:18
05 Room Air	0.000	19:19
06 Std. Gas	0.079	19:19
07 Room Air	0.000	19:20

Lot No = 19817080A4  
Cyl No = 6  
Exp Date = 09/05/2019  
County = 02      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *CALIBRATION CHECK*  
*0.080 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006679  
Location = VALL      8164.14.00 09/16  
05/15/2019      19:26

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:27
02 Std. Gas	0.081	19:27
03 Room Air	0.000	19:27
04 Std. Gas	0.080	19:28
05 Room Air	0.000	19:28
06 Std. Gas	0.081	19:29
07 Room Air	0.000	19:29

Lot No = 19817080A4  
Cyl No = 6  
Exp Date = 09/05/2019  
County = 02      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *CALIBRATION CHECK*  
*0.080 AC*

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006679  
Location = VALL      8164.14.00 09/16  
05/15/2019      19:47

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:47
02 Std. Gas	0.080	19:48
03 Room Air	0.000	19:48
04 Std. Gas	0.080	19:48
05 Room Air	0.000	19:49
06 Std. Gas	0.080	19:49
07 Room Air	0.000	19:50

Lot No = 19817080A4  
Cyl No = 6  
Exp Date = 09/05/2019  
County = 02      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000