

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 006685	Instrument Location Garri	son		
Reason for Install/Repair				
Other (Specify)				
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
🔀 6. Run Tests:				
🕱 A. Print Test (Level1, Function P).				
☑ B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
7. Repair and/or Maintenance Performed (if any):				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5)	0496, Form 120-G) and Place it by the Ir	ntoxilyzer® for Use.		
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form		
Field Inspector Signature		Date _ / /		
Nichard Orfuson		7/12/2019		
Crime Laboratory Use Only		/ /		
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from				
the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the				
County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.				
Reviewed/Certified By		Certified Date		

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006685
Location = GARR 8164.14.00 09/16
07/12/2019 09:57

********** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34917080A3 Standard Cyl #? Standard Expiration? 02/05/2020 Oper No? 118707

Flow Cal. Date: 09/06/2016 Slope 705 Intercept -687370

IR Calibration Date: 09/06/2016 3um 9um

 0th Coef(*100):
 -14155
 -22915

 1st Coef(*100):
 271436
 135296

 2nd Coef(*100):
 2016
 1251

 H2O adj(mg/l*10k):
 372
 505

******* Printer Test End *******

Operator Signature RICHARD JOHNSON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GARR
07/12/2019

Alcohol Analyzer SN 80-006685 8164.14.00 09/16 09:57

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:58
02 Std. Gas	0.078	09:58
03 Room Air	0.000	09:59
04 Std. Gas	0.078	09:59
05 Room Air	0.000	10:00
06 Std. Gas	0.079	10:00
07 Room Air	0.000	10:01

Lot No = 34917080A3

Cyl No = 11

Exp Date = 02/05/2020

County = 28

Oper No. = 118707

Operator Signature RICHARD JOHNSON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006685
Location = GARR 8164.14.00 09/16
07/12/2019 10:02

Test		AC	Time
01 Room Air	Test	0.000	10:04
02 *Subject		RFI*	10:05
03 Room Air		0.000	10:05

*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST T Sub DOB = 01/01/2000 Sub Sex = Male

Sub Sex = Male Weight = 100
Test = OTH Cit = 12345

Dr. Lic. = ND/TES001234

Lot No = 34917080A3Cyl No = 11

Expiration Date = 02/05/2020

County = 28 Oper No. = 118707

Operator Signature RICHARD JOHNSON

Remarks:

Form 106-I8000