



# INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program

SFN 59281 (10-2016)

|   |                                    |
|---|------------------------------------|
| Serial Number<br><b>80-106687</b>   | Instrument Location<br><b>PLPD</b> |
| Reason for Install/Repair<br><input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change<br><input type="checkbox"/> Other (Specify) _____ |                                    |

### Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P). **JSL**
  - B. ACA Test (Level 1, Function C). **JSL**
  - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test). **JSL**
- 7. Repair and/or Maintenance Performed (if Any): **Install**
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

CEE  
7.24.19  
RECEIVED

|                               |                        |
|-------------------------------|------------------------|
| Field Inspector Signature<br> | Date<br><b>7/12/19</b> |
|-------------------------------|------------------------|

### Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

|                           |                                  |
|---------------------------|----------------------------------|
| Reviewed/Certified by<br> | Certified Date<br><b>7/24/19</b> |
|---------------------------|----------------------------------|

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006687  
Location = PLPD                8164.14.00 09/16  
07/12/2019                      21:37

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=-|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=-|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            34418080A2  
Standard Cyl #?            43  
Standard Expiration?       02/05/2021  
Oper No?                    120804

Flow Cal. Date:            10/05/2018  
Slope                      696  
Intercept                  -609977

IR Calibration Date:        10/05/2018  
                                    3um                      9um

|                    | 3um    | 9um    |
|--------------------|--------|--------|
| 0th Coef(*100):    | -11478 | -18771 |
| 1st Coef(*100):    | 255256 | 132135 |
| 2nd Coef(*100):    | 4432   | 1676   |
| H2O adj(mg/l*10k): | 604    | 575    |

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
JASON CARTIER

Remarks:


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006687  
Location = PLPD      8164.14.00 09/16  
07/12/2019      21:38

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 21:39 |
| 02 Std. Gas | 0.080 | 21:40 |
| 03 Room Air | 0.000 | 21:40 |
| 04 Std. Gas | 0.080 | 21:40 |
| 05 Room Air | 0.000 | 21:41 |
| 06 Std. Gas | 0.080 | 21:41 |
| 07 Room Air | 0.000 | 21:42 |

Lot No = 34418080A2  
Cyl No = 43  
Exp Date = 02/05/2021  
County = 07      Oper No. = 120804

  
\_\_\_\_\_  
Operator Signature  
JASON CARTIER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006687  
Location = PLPD      8164.14.00 09/16  
07/12/2019      21:43

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 21:44 |
| 02 Std. Gas | 0.080 | 21:44 |
| 03 Room Air | RFI*  | 21:44 |
| 04 Room Air | RFI*  | 21:45 |

\*Invalid Test  
Inhibited - RFI

Lot No = 34418080A2  
Cyl No = 43  
Exp Date = 02/05/2021  
County = 07      Oper No. = 120804



Operator Signature  
JASON CARTIER

Remarks:

Form 106-I8000