

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-007087

A. GENERAL Setup and Checks:

- Power On
 5 min. countdown during initial warm-up
 1 min. countdown to begin test after standby mode
- Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:
 Level 1 (Operator) Level 2 (Field Inspector) Level 3 (Lab Personnel)
- Printer Setup (Level 3,S,P)
- Forms Load (CMI, Inc., CoBrA Software)

B. LEVEL 1 Setup and Checks:

- Internal Printer Test and Paper Feed (Level 1,P). Attach record.
- External Printer Test (Level 1,P). Attach record.
3. Barcode Reader (Level 1,B):
 Operator ID ND License R key: Refusal. Attach record.

C. LEVEL 2 Setup and Checks:

- Continuous Air Blank (Level 2,A)
- Diagnostic (Level 2,D). Attach record.
- Prelim Data Entry (Level 2,E)
4. Ser Number/Version Info (Level 2,V):
Inst Ser#: 80-007087 DSP Ver: 1362.14
Core Ver: 8164.16 CPLD Ver: 7

D. LEVEL 3 Setup and Checks:

- DVM Monitor (Level 3,D,D): 3: 12638 s: 1 C: 47.00
9: 12909 s: 1 B: 45 0
(RFI indicator with keyed radio): 1
2. Function (F) Keys: : F5 – Tone F9 - Fan
 F6 – Green LED F10 – Pump
 F7 – Red LED F11 – Solenoid
3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 951 CF: 1.06
4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 185 F: 0.000 L/s

- 5. Temperature Monitor (Level 3,D,T): C: 47.0 A: 45.5 I: 0.0
B: 45.0 V: 0.0
- 6. Tank Monitor (Level 3,D,G): P (display): 421 psi P (gauge): 425 psi
- 7. General Setup (Level 3,S,G)
- 8. Configure Start Test (Level 3,S,T)
- 9. Network (Level 3,S,N)

E. TESTS (attach all test records):

Configure wet bath simulator standards before running tests (Level 1,S):

- 1. Zero AC (Level 1,C):
Sim ser #: DR3841 Lot #: MILLI Q H₂O Exp date: NA
Known value: 0.000 AC
Results: 0.000 AC 0.000 AC 0.000 AC
- 2. Low AC (Level 1,C):
Sim ser #: DR3378 Lot #: 18020 Exp date: 1.9.20
Known value: 0.020 AC
Results: 0.019 AC 0.019 AC 0.018 AC
- 3. High AC (Level 1,C):
Sim ser #: DR7350 Lot #: 2018034 Exp date: 3.22.20
Known value: 0.300 AC
Results: 0.294 AC 0.294 AC 0.295 AC

Configure dry gas standard before running tests (Level 1,S):

- 4. Ethanol + 0.05% Acetone (Level 1,B):
Sim ser #: DR3451 Lot #: ICS5 Exp date: NA
Known value: 0.10 AC + 0.05% ACETONE
 Display reads "Interferent Detect"
- 5. Acetone (Level 1,B):
Sim ser #: DR5188 Lot #: ACE1 Exp date: NA
Known value: 0.05% ACETONE
 Display reads "Interferent Detect"
- 6. Methanol (Level 1,B):
Sim ser #: DR3445 Lot #: MET2 Exp date: NA
Known value: 0.10% METHANOL
 Display reads "Interferent Detect"
- 7. Isopropanol (Level 1,B):
Sim ser #: DR5189 Lot #: IPA2 Exp date: NA
Known value: 0.10% ISOPROPANOL
 Display reads "Interferent Detect"

- 8. Ambient Fail (Level 1,B): Display and test record read correct.
- 9. Invalid Sample (Level 1,B): Display and test record read correct.
- 10. Sequence Aborted (Level 1,B): Display and test record read correct.
- 11. ABA (Level 1,B): Test record reads correct.
- 12. Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
- 13. RFI Detected (B,C, or CMS mode): Display and test record read correct.
- 14. Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 15. Deficient Samples (CMS mode): Display and test record read correct.
- 16. Improper Sample (CMS mode): Display and test record read correct.
- 17. Difference Too Great (CMS mode): Display and test record read correct.
- 18. CMS (CMS mode): Test record reads correct.
- 19. Stability Test (Level 3,T,S):

Ave: 0.0797 Std Dev: 0.0005 Rel Std Dev: 0.59%

Maintenance/Repairs/Notes: INSPECTION AFTER RECEIVING
FROM CMI FACTORY AFTER WARRANTY REPAIR.

Instrument is acceptable to be used in the field.

Charles E. Ed
Field Inspector's Signature

7/25/19
Date

NA
Reviewed by

NA
Date

Intoxilyzer Test Record and Checklist
NDDAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 14:14

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstuvw
ABCDEFGHIJKLMNOPQRSTUV

32 characters per line
abcdefghijklmnopqrstuvwxy
ABCDEFGHIJKLMNOPQRSTUVWXYZ

42 characters per line
abcdefghijklmnopqrstuvwxy
ABCDEFGHIJKLMNOPQRSTUVWXYZ1234567890

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On MemFull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 34917080A3
Standard Cyl #? 28
Standard Expiration? 02/05/2020
Oper No? 666666

Flow Cal. Date: 07/02/2019
Slope 721
Intercept -792109

IR Calibration Date: 07/02/2019
3um 9um

0th Coef(*100):	-45706	-18865
1st Coef(*100):	253208	131287
2nd Coef(*100):	2874	1457
H2O adj(mg/l*10k):	1171	483

***** Test End *****



Operator Signature
CHARLES EDER

Remarks: INTERNAL

PRINT TEST

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 14:16

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 34917080A3
Standard Cyl #? 28
Standard Expiration? 02/05/2020
Oper No? 666666

Flow Cal. Date: 07/02/2019
Slope 721
Intercept -792109

IR Calibration Date: 07/02/2019

	3um	9um
0th Coef(*100):	-45706	-18865
1st Coef(*100):	253208	131287
2nd Coef(*100):	2874	1457
H2O adj(mg/l*10k):	1171	483

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks: EXTERNAL PRINT

TEST

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 14:39

Test	AC	Time
01 Room Air	0.000	14:40
02 *Subject Test	REF*	14:40
03 Room Air	0.000	14:40

*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = R KEY REFUSAL
Dr. Lic. = ND/DIS921456
Lot No = 34917080A3
Cyl No = 28
Expiration Date = 02/05/2020
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: R KEY REFUSAL

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 14:42

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

DIAGNOSTIC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 14:56

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:57
02 Std. Sol.	0.000	14:58
03 Room Air	0.000	14:58
04 Std. Sol.	0.000	14:59
05 Room Air	0.000	14:59
06 Std. Sol.	0.000	15:00
07 Room Air	0.000	15:01

08 Sim Temp = 34.0°C

Simul Ser No = DR3841
Std Sol No = MILLIQ H2O
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

ZERO AC

0.000 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 15:03

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:03
02 Std. Sol.	0.019	15:04
03 Room Air	0.000	15:05
04 Std. Sol.	0.019	15:05
05 Room Air	0.000	15:06
06 Std. Sol.	0.018	15:06
07 Room Air	0.000	15:07

08 Sim Temp = 34.0°C

Simul Ser No = DR3378
Std Sol No = 18020
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Low AC*
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 15:15

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:15
02 Std. Sol.	0.294	15:16
03 Room Air	0.000	15:17
04 Std. Sol.	0.294	15:18
05 Room Air	0.000	15:18
06 Std. Sol.	0.295	15:19
07 Room Air	0.000	15:19

08 Sim Temp = 34.0°C

Simul Ser No = DR7350
Std Sol No = 201803H
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

HIGH AC
0.300 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 15:22

Test	AC	Time
01 Room Air	0.000	15:22
02 *Subject Test	INT*	15:23
03 Room Air	0.000	15:23

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS921456
Lot No = 34418080A2
Cyl No = 40
Expiration Date = 02/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent Check

CEE ~~0.10ACA~~ 0.10AC + 0.05% ACETONE

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 15:25

Test	AC	Time
01 Room Air	0.000	15:26
02 *Subject Test	INT*	15:26
03 Room Air	0.000	15:27

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS921456
Lot No = 34418080A2
Cyl No = 40
Expiration Date = 02/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent Check
0.05% ACETONE

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 15:28

Test	AC	Time
01 Room Air	0.000	15:29
02 *Subject Test	INT*	15:30
03 Room Air	0.000	15:30

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS921456
Lot No = 34418080A2
Cyl No = 40
Expiration Date = 02/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Interferent Check*
0.10% METHANOL

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 15:32

Test	AC	Time
01 Room Air	0.000	15:33
02 *Subject Test	INT*	15:33
03 Room Air	0.000	15:34

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS921456
Lot No = 34418080A2
Cyl No = 40
Expiration Date = 02/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Interferent Check*

CEE 0.10A 0.10% ISOPROPAWOL

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 15:39

Test	AC	Time
01 Room Air	AMB*	15:40
02 Room Air	0.000	15:40

*Invalid Test
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = AMBIENT FAIL
Dr. Lic. = ND/DIS921456
Lot No = 34418080A2
Cyl No = 40
Expiration Date = 02/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *AMBIENT FAIL*

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 15:41

Test	AC	Time
01 Room Air	0.000	15:42
02 Invalid Sample	X.XXX	15:42
03 Room Air	0.000	15:43

*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = INVALID SAMPLE
Dr. Lic. = ND/DIS921456
Lot No = 34418080A2
Cyl No = 40
Expiration Date = 02/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

INVALID SAMPLE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 15:44

Test	AC	Time
01 Room Air	0.000	15:45
02 *Subject Test	ABT*	15:45
03 Room Air	0.000	15:45

*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = SEQUENCE ABORT
Dr. Lic. = ND/DIS921456
Lot No = 34418080A2
Cyl No = 40
Expiration Date = 02/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: SEQUENCE ABORTED

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 15:46

Test	AC	Time
01 Room Air	0.000	15:47
02 Subject Test 1	0.082	15:47
03 Room Air	0.000	15:48
04 Reported AC	0.082	15:47

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = ABA MODE
Dr. Lic. = ND/DIS921456
Lot No = 34418080A2
Cyl No = 40
Expiration Date = 02/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

ABA MODE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 15:49

Test	AC	Time
01 Diagnostic	OK	15:50
02 Room Air	0.000	15:50
03 Subject Test 1	0.081	15:51
04 Room Air	0.000	15:53
05 Std. Gas	0.088*	15:54
06 Room Air	0.000	15:55

*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = CAL CHK OUT TOL
Dr. Lic. = ND/DIS921456
Lot No = 34418080A2
Cyl No = 40
Expiration Date = 02/05/2021
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
OUT OF TOLERANCE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 15:56

Test	AC	Time
01 Diagnostic	OK	15:57
02 Room Air	0.000	15:58
03 Subject Test 1	0.081	15:58
04 Room Air	0.000	16:01
05 Std. Gas	0.080	16:02
06 Room Air	0.000	16:03
07 *Subject Test	RFI*	16:03
08 Room Air	0.000	16:04

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI DETECT
Dr. Lic. = ND/DIS921456
Lot No = 34418080A2
Cyl No = 40
Expiration Date = 02/05/2021
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *RFI DETECT*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 16:05

Test	AC	Time
01 Diagnostic	OK	16:06
02 Room Air	0.000	16:06
03 *Subject Test	0.000*	16:09
04 Room Air	0.000	16:12
05 Std. Gas	0.080	16:13
06 Room Air	0.000	16:14
07 *Subject Test	0.000*	16:17
08 Room Air	0.000	16:18

No RFI Detected

*Deficient Sample - Value Printed was
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = DEFICIENT SAMP
Dr. Lic. = ND/DIS921456
Lot No = 34418080A2
Cyl No = 40
Expiration Date = 02/05/2021
County = 08 Oper No. = 666666

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

DEFICIENT SAMPLES

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 16:19

Test	AC	Time
01 Diagnostic	OK	16:19
02 Room Air	0.000	16:20
03 Reference	IPS*	16:20
04 Room Air	0.000	16:21

*Improper Sample

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = IMPROPER SAMPLE
Dr. Lic. = ND/DIS921456
Lot No = 34418080A2
Cyl No = 40
Expiration Date = 02/05/2021
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *IMPROPER SAMPLE*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 16:21

Test	AC	Time
01 Diagnostic	OK	16:22
02 Room Air	0.000	16:23
03 Subject Test 1	0.151*	16:23
04 Room Air	0.000	16:25
05 Std. Gas	0.081	16:26
06 Room Air	0.000	16:28
07 Subject Test 2	0.080*	16:28
08 Room Air	0.000	16:29

*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = DIFF TOO GREAT
Dr. Lic. = ND/DIS921456
Lot No = 34418080A2
Cyl No = 40
Expiration Date = 02/05/2021
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

DIFFERENCE
TOO GREAT

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 16:30

Test	AC	Time
01 Diagnostic	OK	16:31
02 Room Air	0.000	16:31
03 Subject Test 1	0.081	16:32
04 Room Air	0.000	16:34
05 Std. Gas	0.080	16:35
06 Room Air	0.000	16:36
07 Subject Test 2	0.080	16:37
08 Room Air	0.000	16:38
09 Reported AC	0.080	16:37

Difference OK
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = CMS MODE
Dr. Lic. = ND/DIS921456
Lot No = 34418080A2
Cyl No = 40
Expiration Date = 02/05/2021
County = 08 Oper No. = 666666

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

CMS MODE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = TOXL 8164.16.00 09/18
07/25/2019 16:39

Test	AC	Time
01 Room Air	0.000	16:40
02 Std. Gas	0.080	16:40
03 Room Air	0.000	16:41
04 Std. Gas	0.080	16:41
05 Room Air	0.000	16:41
06 Std. Gas	0.080	16:42
07 Room Air	0.000	16:42
08 Std. Gas	0.080	16:43
09 Room Air	0.000	16:43
10 Std. Gas	0.079	16:44
11 Room Air	0.000	16:44
12 Std. Gas	0.080	16:44
13 Room Air	0.000	16:45
14 Std. Gas	0.080	16:45
15 Room Air	0.000	16:46
16 Std. Gas	0.080	16:46
17 Room Air	0.000	16:47
18 Std. Gas	0.080	16:47
19 Room Air	0.000	16:48
20 Std. Gas	0.079	16:48
21 Room Air	0.000	16:48
22 Std. Gas	0.080	16:49
23 Room Air	0.000	16:49
24 Std. Gas	0.080	16:50
25 Room Air	0.000	16:50
26 Std. Gas	0.080	16:50
27 Room Air	0.000	16:51
28 Std. Gas	0.079	16:51
29 Room Air	0.000	16:52
30 Std. Gas	0.080	16:52
31 Room Air	0.000	16:53
32 Std. Gas	0.079	16:53
33 Room Air	0.000	16:54
34 Std. Gas	0.079	16:54
35 Room Air	0.000	16:54
36 Std. Gas	0.079	16:55
37 Room Air	0.000	16:55
38 Std. Gas	0.080	16:56
39 Room Air	0.000	16:56
40 Std. Gas	0.080	16:57
41 Room Air	0.000	16:57
42 Cal Check Stats:		
43 Average	0.0797	
44 Std Dev	0.0005	
45 Rel Std Dev	0.59%	



Operator Signature
CHARLES EDER

Remarks: STABILITY TEST
0-080 AC

Form 106-I8000

LOT: 3441808042

EXP: 2/5/21

CYL: 040