



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Jeffrey D. Fleck</i>		
Location <i>LaMoure Co Courthouse</i>		Intoxilyzer® Serial Number <i>80-006680</i>
Gas Lot Number <i>34418080A2</i>	Gas Cylinder Number <i>7</i>	Gas Expiration Date <i>2-5-2021</i>

Check When Done:

*CEE
9-16-19*

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.013 AC
 - B. 0.017 AC
 - C. 0.017 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Jeffrey D. Fleck</i>	Date <i>9/10/2019</i>
Reviewed By (Crime Laboratory Use Only) <i>Charles E. Egan</i>	Date <i>9-16-19</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006680
Location = LAMO 8164.14.00 09/16
09/10/2019 15:10

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:10
02 Std. Gas	0.078	15:11
03 Room Air	0.000	15:11
04 Std. Gas	0.077	15:11
05 Room Air	0.000	15:12
06 Std. Gas	0.077	15:12
07 Room Air	0.000	15:13

Lot No = 34418080A2
Cyl No = 7
Exp Date = 02/05/2021
County = 23 Oper No. = 010605


Operator Signature
JEFFREY D. FLECK

Remarks:

Form 106-I8000