

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

0.44		
Serial Number 80 - 00 70 87	Instrument Location	Police Department
Reason for Install/Repair	OCTORICIT	Torre Derai Man
Install After Receiving From Crime Laboratory Install	After Location Change	
Other (Specify)		
Check When Done:	3	
Surge Protector Installed/Property Grounded.		
2., Telephone Line Connected to Intoxilyzer® 8000.	/A)	
☑ 3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Lo	cation; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function 5)		
☐ 6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode	or Level 1. Function B or C: Key Radio F	uring Teet)
7. Repair and/or Maintenance Performed (if any):	/A	diling rest).
8. Complete the Top Portion of the Intoxilyzer® Record (SF	N50496, Form 120-G) and Place it by the	Intoviluzer® for I lea
9. File Previous Intoxilyzer® Record (SFN504096, Form 12	0-G) at the Intoxilvzer® Location at the Ac	ency
10. Send the Following to the Crime Laboratory: Completed 104-G), Print Test, ACA Test, and REF Test.	Intoxilyzer® 8000 Installation and Repair	Checkout (SFN59281, Form
Field Inspector Signature		Date 11-12-19
Crime Laboratory Use Only		1777
This installation has been reviewed and the instrument is approved to the date the Field Inspector performed the installation. This record on County of Burleigh, North Dakota, is certified to be a true and correct	to be used for the analysis of breath to deter in file at the Office of Attorney General, Cri copy of the documents received.	ermine alcohol concentration from me Laboratory Division, in the
Reviewed/Certified By	0	Certified Date
Charles En	L	11.15.2019

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = WISH 8164.16.00 09/18
11/12/2019 12:11

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Yes Display Third Digit? Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34418080A2 Standard Cyl #? Standard Expiration? 02/05/2021 Oper No? 131099

Flow Cal. Date: 07/02/2019 Slope 721

Intercept -792109

IR Calibration Date: 07/02/2019
3um 9um

Oth Coef(*100): -45706 -18865
1st Coef(*100): 253208 131287
2nd Coef(*100): 2874 1457
H2O adj(mg/l*10k): 1171 483

******* Printer Test End ********

Operator Signature LUCAS KUNTZ

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WISH
11/12/2019

Alcohol Analyzer SN 80-007087 8164.16.00 09/18 12:12

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:23
02 Std. Gas	0.080	12:23
03 Room Air	0.000	12:24
04 Std. Gas	0.081	12:24
05 Room Air	0.000	12:24
06 Std. Gas	0.081	12:25
07 Room Air	0.000	12:25

Lot No = 34418080A2

Cyl No = 12

Exp Date = 02/05/2021

Operator Signature LUCAS KUNTZ

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007087
Location = WISH 8164.16.00 09/18
11/12/2019 12:30

Test	AC	Time
01 Room Air	0.000	12:32
02 *Subject Test	RFI*	12:33
03 Room Air	0.000	12:33

*Invalid Test Inhibited - RFI

Sub Name = RFI TEST, RFI TEST TEST

Sub DOB = 01/01/1990

Sub Sex = Male Weight = 180
Test = OTH Cit = NONE

Dr. Lic. = ND/TES901234

Lot No = 34418080A2

Cyl No = 12

Expiration Date = 02/05/2021

County = 26 Oper No. = 131099

Operator Signature LUCAS KUNTZ

Remarks:

Form 106-I8000