

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 8:0-006513	Instrument Location Bowman County Sh	C.S.C		
Reason for Install/Repair				
Install After Receiving From Crime Laboratory				
Cher (Specify)				
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
6. Run Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
7. Repair and/or Maintenance Performed (if any): CEE NA CEE 1/27/20				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.				
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.				
Field Inspector Signature		Date 1/24/20		

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

	1	
Reviewed/Certified By	Charlacter	Certified Date

l

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006513 8164.13.00 06/09 Location = BOWM10:02 01/24/2020 \*\*\*\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*\*\*\*\*\*\*\*\* abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Yes Display Prelim Rslt? Yes Display Third Digit? Inhib Printer(Y/N)? No No Display Volume? Yes Disable On Memfull? # of Print Copies?
Select Std (D/W/I)? 1 Dry 0.080 Standard Value? 34418080A2 Standard Lot #? Standard Cyl #? Standard Cyl #? Standard Expiration? 02/05/2021 010703 41 Oper No? 08/20/2015 Flow Cal. Date: 680 Slope -696037 Intercept IR Calibration Date: 08/21/2015 3um 9u 9um \_\_\_\_\_ Oth Coef(\*100):-14436-136101st Coef(\*100):2654901361472nd Coef(\*100):44221710 448 H2O adj(mg/l\*10k): 672

\*\*\*\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*\*\*\*

Operator Signature JASON BULLIS

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer<br/>North Dakota Model 8000Alcohol Analyzer<br/>SN 80-006513Location = BOWM<br/>01/24/20208164.13.00 06/09<br/>10:02

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.082 0.000 0.082 0.000 0.082 0.000	10:03 10:03 10:04 10:04 10:05 10:05 10:05

Lot No = 34418080A2 Cyl No = 41 Exp Date = 02/05/2021 County = 06 Oper No. = 010703

Operator Signature JASON BULLIS

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 Alcohol Analyzer CMI, Inc. Intoxilyzer SN 80-006513 North Dakota Model 8000 8164.13.00 06/09 Location = BOWM01/24/2020 10:12 Time AC Test 01 Room Air RFI\* 10:14 10:14 RFI\* 02 Room Air \*Invalid Test Inhibited - RFI Sub Name = DOE, JOHN NA Sub DOB = 01/24/2000Sub Sex = Unknown Weight = NA Cit = NATest = OTHDr. Lic. = ND/NA Lot No = 34418080A2Cyl No = 41Expiration Date = 02/05/2021 Oper No. = 010703County = 06

Operator Signature JASON BULLIS

Remarks:

Form 106-18000