

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-006674	Instrument Location WASH			
Reason for Install/Repair				
🛒 Install After Receiving From Crime Laboratory 🛛 🗌 Install After Location Change				
Other (Specify)				
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
6. Run Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
7. Repair and/or Maintenance Performed (if any): None				
2. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.				
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.				
Field Inspector Signature				
hichod phuson	Date 5/6/2020			
Crime Laboratory Use Only				
This is the first first first sector of the				

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	Chalatel	Certified Date 5 • 26 • 2020

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006674 Location = WASH 8164.14.00 09/16 05/06/2020 14:40 abcdefghijklmnopgrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry 0.080 Standard Lot #? Standard Cyl #? Standard T 24119080A1 22 Standard Expiration? 11/05/2021 Oper No? 118707 Flow Cal. Date: 06/28/2017 707 Slope Intercept -708478 IR Calibration Date: 06/28/2017 3um 9um ______ 0th Coef(*100):-39151st Coef(*100):2565202nd Coef(*100):3843 -18046 127506 1644 H2O adj(mg/l*10k): 428568

Operator Signature RICHARD JOHNSON

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006674Location = WASH8164.14.00 09/1605/06/202014:40

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.080 0.000 0.081 0.000 0.079 0.000	14:41 14:41 14:42 14:42 14:42 14:43 14:43

Lot No = 24119080A1 Cyl No = 22 Exp Date = 11/05/2021 County = 28

Oper No. = 118707

Operator Signature RICHARD JOHNSON

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006674Location = WASH8164.14.00 09/1605/06/202014:50

DRY CAL	CHECK	
Test	AC	Time
01 Room Air 02 Room Air	RFI* 0.000	14:50 14:50
*Invalid Test Inhibited - RFI		
Lot No = 24119080A1 Cyl No = 22 Exp Date = 11/05/2021 County = 28	Oper No.	= 118707

Operator Signature RICHARD JOHNSON

Remarks:

Form 106-18000