# OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

## **INTOXILYZER® 8000 INSPECTION**

Intoxi	lyzer® 8000 Serial Number: <u>80</u>	<u>-006508</u>	Location: 10xL
Α.	General Setup and Checks:  1.	n code (Level 2,E).  R Sign and attach test r	e-set if necessary. ecord. egulator within 50 psi.
B.	Sim. Ser #: <u><i>Mf 3</i></u> b.	ollowing tests (Level 1 1,C): Known Value ≤ 061 Lot #: 2018101 1,C): Known Value ≥	,S). 0.03 AC: 0・0 2 0 AC Exp. Date: 10・24・20 0.25 AC: 0.300 AC Exp. Date: 2・20・21
	Sim. Ser #: <u>DE3 451</u> Legar Sim. Ser #: <u>DE3 451 Legar Sim. Ser #: DE3 451 Legar Sim. Ser #: <u>DE3 451 Legar Sim. S</u></u>	rel 1,B): Known Value of #: \( \frac{TCS}{6} \) Exp. rent Detect"  e) \( \sqrt{Display re} \) (Level 1,C): Known \( \frac{80A_1}{1} \) Cyl #: \( \frac{47}{1} \) Test 4 \( \frac{0.08_1}{0.08_1} \) Test 6 \( \frac{0.08_1}{0.08_1} \)	e: <u>0./0</u> AC + 0.05% Acetone Date: <u>// A</u> ads "RFI Detect" /alue <u>0.080</u> AC _Exp. Date: <u>4.5.22</u> Test 7 <u>0.081</u> Test 8 <u>0.081</u>
C.	Remarks/Maintenance: An REPLACES BREATH H		2C770N
,	trument is acceptable to be used	in the field.	5/27/20
Breath	n Analyst Signature	_	5/27/20 Date WA-
Revie	wed by		Date

02/18

Issuing Authority - CEE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006508
Location = TOXL 8164.13.00 06/09
05/27/2020 11:37

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 05620080A1 Standard Cyl #? 47 04/05/2022 Standard Expiration? Oper No? 666666

Flow Cal. Date: 08/20/2015 Slope 644 Intercept -568181

IR Calibration Date: 08/20/2015
3um 9um

Oth Coef(\*100): -22365 -18627
1st Coef(\*100): 256969 136095

 0th Coef(\*100):
 -22365
 -1862

 1st Coef(\*100):
 256969
 136095

 2nd Coef(\*100):
 3168
 1444

 H2O adj(mg/l\*10k):
 693
 453

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature CHARLES EDER

PRINT TEST

Remarks:

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 05/27/2020 SN 80-006508 8164.13.00 06/09 11:38

## WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:39
02 Std. Sol.	0.018	11:39
03 Room Air	0.000	11:40
04 Std. Sol.	0.018	11:41
05 Room Air	0.000	11:41
06 Std. Sol.	0.018	11:42
07 Room Air	0.000	11:43

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = MP3061 Std Sol No = 201810D

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

LOW AC 0.020 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
05/27/2020

SN 80-006508 8164.13.00 06/09 11:44

#### WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:44
02 Std. Sol.	0.300	11:45
03 Room Air	0.000	11:46
04 Std. Sol.	0.303	11:46
05 Room Air	0.000	11:47
06 Std. Sol.	0.302	11:48
07 Room Air	0.000	11:48

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = MP3069 Std Sol No = 201902D

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

HIGH AC 0.300 AC Form 106-18000

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006508
Location = TOXL 8164.13.00 06/09
05/27/2020 11:51

Test	AC	Time
01 Room Air	0.000	11:52
02 *Subject Test	INT*	11:52
03 Room Air	0.000	11:53

\*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT

Dr. Lic. = ND/DIS921456

Lot No = 05620080A1

Cyl No = 47

Expiration Date = 04/05/2022

County = 08 Oper No. = 666666

Operator Signature CHARLES EDER

Remarks: 7

INTERFERENT CHECK 0.10 AR + 0.05% ACETONE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006508
Location = TOXL 8164.13.00 06/09
05/27/2020 11:54

Test		AC	Time
01	Diagnostic	OK	11:54
02	Room Air	0.000	11:55
03	*Subject Test	RFI*	11:55
	Room Air	0.000	11:56

\*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS921456

Lot No = 05620080A1

Cyl No = 47

Expiration Date = 04/05/2022

County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
05/27/2020

SN 80-006508 8164.13.00 06/09 12:02

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:03
02 Std. Gas	0.081	12:03
03 Room Air	0.000	12:04
04 Std. Gas	0.081	12:04
05 Room Air	0.000	12:05
06 Std. Gas	0.081	12:05
07 Room Air	0.000	12:06

Lot No = 05620080A1

Cyl No = 47

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

ALIBRATION CHECK 0.080 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
05/27/2020

SN 80-006508 8164.13.00 06/09 12:06

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:07
02 Std. Gas	0.082	12:07
03 Room Air	0.000	12:08
04 Std. Gas	0.081	12:08
05 Room Air	0.000	12:09
06 Std. Gas	0.081	12:09
07 Room Air	0.000	12:09

Lot No = 05620080A1

Cyl No = 47

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

0.080AC

PALIBRATION CHECK

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
05/27/2020

SN 80-006508 8164.13.00 06/09 12:10

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:10
02 Std. Gas	0.081	12:11
03 Room Air	0.000	12:11
04 Std. Gas	0.081	12:12
05 Room Air	0.000	12:12
06 Std. Gas	0.081	12:12
07 Room Air	0.000	12:13

Lot No = 05620080A1

Cyl No = 47

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature
CHARLES EDER

ALIBRATION CHECK

Domarke.

Form 106-I8000

0.080 AC