Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8	8000 Serial	Number: 8	0-000	676
----------------	-------------	-----------	-------	-----

Location: ToxL

- Α. General Setup and Checks:
 - Diagnostic Tests Pass and Instrument in Ready Mode 1.
 - 2. Breath tube heated
 - 3 Date, time and location code (Level 2,E). Re-set if necessary.
 - Print test (Level 1,P). Sign and attach test record. 4.
 - Tank monitor (Level 3,D,G). Display and Regulator within 50 psi. 5.
 - Display: 347 psi Regulator: 375 psi
- Β. Tests (Sign and attach test records):

 \square Configure simulator for the following tests (Level 1,S).

- Wet Calibration Check 1.
 - \square Low AC (Level 1,C): Known Value ≤ 0.03 AC: 0.020 AC a. Sim. Ser #: MP3064 Lot #: 20070 Exp. Date: 2.13.22
 - High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC b. Sim. Ser #: MP 3069 Lot #: 201902 D Exp. Date: 2.20.21

Configure dry gas standard for the following tests (Level 1,S).

- 2. Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone Display reads "Interferent Detect"
- 3.
- 4. $\begin{array}{c} \text{Gas Cyl Lot } \# : \underline{O5} \underbrace{O20080 \text{ Al Cyl } \# : \underline{36} \quad \text{Exp. Date: } \underline{4 \cdot 5 \cdot 26} \\ \text{Test 1 } \underbrace{O \cdot 081} \quad \text{Test 4 } \underbrace{O \cdot 081} \quad \text{Test 7 } \underbrace{O \cdot 081} \\ \text{Test 2 } \underbrace{O \cdot 081} \quad \text{Test 5 } \underbrace{O \cdot 081} \quad \text{Test 8 } \underbrace{O \cdot 082} \\ \text{Test 3 } \underbrace{O \cdot 081} \quad \text{Test 6 } \underbrace{O \cdot 081} \quad \text{Test 9 } \underbrace{O \cdot 081} \\ \end{array}$
- Remarks/Maintenance: <u>ANNUAL INSPECTION</u> C.

Average 0.081

A Instrument is acceptable to be used in the field.

Breath Analyst Signature

Reviewed by

02/18 Issuing Authority - CEE

SN 80-006676

 $\frac{5 \cdot 27 \cdot 20}{\text{Date}}$

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006676 Location = TOXL 8164.14.00 09/16 05/27/2020 16:29 abcdefghijklmnopgrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopgrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies?
Select Std (D/W/I)? 1 Dry Standard Value? 0.080 Standard Lot #? 05620080A1 Standard Cyl #? 36 Standard Expiration? 04/05/2022 Oper No? 666666 Flow Cal. Date: 08/22/2016 Slope 710 Intercept -713681 IR Calibration Date: 08/22/2016 3um 9um Oth Coef(*100):-166671st Coef(*100):2591902nd Coef(*100):3393 -22327 132291 1374 H2O adj(mg/l*10k): 611546

Charla El

Operator Signature CHARLES EDER

Remarks:

INT TEST

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006676Location = TOXL8164.14.00 09/1605/27/202016:30

	WET CAL CHECK]
Test	AC	Time
01 Room Air 02 Std. Sol. 03 Room Air 04 Std. Sol. 05 Room Air 06 Std. Sol.	0.000 0.017 0.000 0.017 0.000 0.017	16:31 16:32 16:32 16:33 16:33 16:33
07 Room Air	0.000	16:35

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = MP3064 Std Sol No = 20070 County = 08

Oper No. = 666666

IN L

Operator Signature CHARLES EDER

Low AC

0.020AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006676Location = TOXL8164.14.00 09/1605/27/202016:41

	WET CAL CHECK	
Test	AC	Time
01 Room Air	0.000	16:41
02 Std. Sol.	0.298	16:42
03 Room Air	0.000	16:42
04 Std. Sol.	0.298	16:43
05 Room Air	0.000	16:44
06 Std. Sol.	0.298	16:44
07 Room Air	0.000	16:45

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = MP3069 Std Sol No = 201902D County = 08 Oper No. = 666666

la

Operator Signature CHARLES EDER

HIGH AC 0.300AC Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006676 Location = TOXL 8164.14.00 09/16 05/27/2020 16:46 Test AC Time 01 Room Air 16:47 0.000 02 *Subject Test INT* 16:47 03 Room Air 0.000 16:48 *Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = INTERFERENT Dr. Lic. = ND/DIS921456Lot No = 05620080A1Cyl No = 36Expiration Date = 04/05/2022County = 08Oper No. = 666666

Operator Signature CHARLES EDER

Interferent Check 0.10 AC + 0.05%. ACETONE

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer SN 80-006676 North Dakota Model 8000 Location = TOXL 8164.14.00 09/16 05/27/2020 16:48 AC Time Test 01 Diagnostic 16:49 OK 16:49 02 Room Air 0.000 03 *Subject Test 16:49 RFI* 04 Room Air 0.000 16:50 *Invalid Test Inhibited - RFI Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = RFI CHECK Dr. Lic. = ND/DIS921456Lot No = 05620080A1Cyl No = 36Expiration Date = 04/05/2022County = 08Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006676 Location = TOXL 8164.14.00 09/16 05/27/2020 16:50

[DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.081 0.000 0.081 0.000	16:51 16:51 16:52 16:52 16:53 16:53 16:54

Lot No = 05620080A1Cyl No = 36Exp Date = 04/05/2022County = 08

Oper No. = 666666

1 hi

Operator Signature CHARLES EDER ALIBRATION CHECK

0.080 AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006676 Location = TOXL 8164.14.00 09/16 05/27/2020 16:54

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.081 0.000 0.081 0.000	16:54 16:55 16:55 16:56 16:56 16:57 16:57

Lot No = 05620080A1Cyl No = 36Exp Date = 04/05/2022County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Operator Signature CHARLES EDER ALIBRATION CHECK

0.0804

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000Alcohol Analyzer
SN 80-006676Location = TOXL
05/27/20208164.14.00 09/16
16:57

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.082 0.000 0.081 0.000	16:58 16:59 16:59 17:00 17:00 17:00

Lot No = 05620080A1 Cyl No = 36 Exp Date = 04/05/2022 County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

ALIBRATION CHECK

0.080 AC Form 106-18000