OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

Intoxi	lyzer® 8000 Serial Number: 80-006680 Loc	ation: 70×L
Α.	General Setup and Checks: 1. Diagnostic Tests Pass and Instrument in Ready 2. Breath tube heated 3. Date, time and location code (Level 2,E). Re-set 4. Print test (Level 1,P). Sign and attach test record 5. Tank monitor (Level 3,D,G). Display and Regular a. Display: 999 psi Regulator: 1000 psi	t if necessary.
B.	Tests (Sign and attach test records): Configure simulator for the following tests (Level 1,S). Wet Calibration Check a.	p. Date: <u>/o· 2<i>4</i>· 2</u> o s AC: <u>o. 3oo</u> AC
	2. Interferent Check (Level 1,B): Known Value: O.D. Sim. Ser #: DR3451 Lot #: FCS 6 Exp. Date Display reads "Interferent Detect" 3. XRFI Check (CMS Mode) A Display reads "Interferent Detect" 4. Dry Calibration Check (Level 1,C): Known Value Gas Cyl Lot #: 05620 080 AICyl #: 34 Exp. Date Test 1 0.083 Test 4 0.082 Test Test 2 0.082 Test 5 0.082 Test 3 0.082 Test 6 0.082 Test Average 0.082	OAC + 0.05% Acetone ERFI Detect" O.080 AC Date: 4.5-22 t 7 0.082 t 8 0.082
C.	Remarks/Maintenance: ANNUAL INSPEC	BATTELY
	trument is acceptable to be used in the field. Analyst Signature NA	6-4-20 Date NA
Revie	wed by	Date

02/18

Issuing Authority - CEE

SN 80-006680

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006680
Location = TOXL 8164.14.00 09/16
06/04/2020 11:40

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 05620080A1 34 Standard Expiration? 04/05/2022 Oper No? 666666

Flow Cal. Date: 09/06/2016 Slope 681

Intercept -620591

IR Calibration Date: 09/06/2016 3um 9um

Oth Coef(*100): -14265 -21271
1st Coef(*100): 260743 129687
2nd Coef(*100): 3906 1392
H2O adj(mg/l*10k): 742 684

****** Printer Test End *******

Operator Signature CHARLES EDER

emarks: VO 10/2 To

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/04/2020

Alcohol Analyzer SN 80-006680 8164.14.00 09/16 11:41

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:42
02 Std. Sol.	0.018	11:42
03 Room Air	0.000	11:43
04 Std. Sol.	0.018	11:44
05 Room Air	0.000	11:44
06 Std. Sol.	0.019	11:45
07 Room Air	0.000	11:45

 $08 \text{ Sim Temp} = 34.0 ^{\circ}\text{C}$

Simul Ser No = MP3061 Std Sol No = 201810D

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Low AC

0.018 AC CEE 0.020 AC

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 06/04/2020 Alcohol Analyzer SN 80-006680 8164.14.00 09/16 11:51

WET CAL CHECK

Test		AC	Time
01 Room	Air	0.000	11:52
02 Std.	Sol.	0.294	11:53
03 Room	Air	0.000	11:53
04 Std.	Sol.	0.296	11:54
05 Room	Air	0.000	11:55
06 Std.	Sol.	0.296	11:55
07 Room	Air	0.000	11:56

08 Sim Temp = 34.0°C

Simul Ser No = MP3069 Std Sol No = 201902D

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

416H A

Remarks:

0.300 AC

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006680
Location = TOXL 8164.14.00 09/16
06/04/2020 11:58

Test	AC	Time
01 Room Air	0.000	11:59
02 *Subject Test	INT*	11:59
03 Room Air	0.000	12:00

*Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT

Dr. Lic. = ND/DIS921456

Lot No = 05620080A1

Cyl No = 34

Expiration Date = 04/05/2022

County = 08 Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Interferent Check 0.10AC + 0.05% ACETONE FORM 106-18000

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006680
Location = TOXL 8164.14.00 09/16
06/04/2020 12:01

Test	AC	Time
01 Diagnostic	OK	12:01
02 Room Air	0.000	12:02
03 *Subject Test	RFI*	12:02
04 Room Air	0.000	12:02

*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150 Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS921456

Lot No = 05620080A1

Cyl No = 34

Expiration Date = 04/05/2022

County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/04/2020

Alcohol Analyzer SN 80-006680 8164.14.00 09/16 12:03

DRY CAL CHECK

Test	AC Ti	Time
01 Room Air	0.000	12:03
02 Std. Gas	0.083	12:04
03 Room Air	0.000	12:04
04 Std. Gas	0.082	12:05
05 Room Air	0.000	12:05
06 Std. Gas	0.082	12:05
07 Room Air	0.000	12:06

Lot No = 05620080A1

Cyl No = 34

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

ALIBRATION CHECK

O.080 AC
FORM 106-18000

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/04/2020

Alcohol Analyzer SN 80-006680 8164.14.00 09/16 12:06

DRY CAL CHECK

Test		AC	Time
01 Ro	om Air	0.000	12:07
02 St	d. Gas	0.082	12:07
03 Ro	om Air	0.000	12:08
04 St	d. Gas	0.082	12:08
05 Ro	om Air	0.000	12:08
06 St	d. Gas	0.082	12:09
07 Ro	om Air	0.000	12:09

Lot No = 05620080A1

Cyl No = 34

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks: CAUBRATION CHECK

0.080 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/04/2020

SN 80-006680 8164.14.00 09/16 12:10

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:10
02 Std. Gas	0.082	12:11
03 Room Air	0.000	12:11
04 Std. Gas	0.082	12:12
05 Room Air	0.000	12:12
06 Std. Gas	0.082	12:12
07 Room Air	0.000	12:13

Lot No = 05620080A1

Cyl No = 34

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

SALIBRATION CAECK 0.080AC