OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

Intoxil	yzer® 8000 Serial Number: <u></u> <u></u> <u></u> 80 - 00 6688	Location: ToxL
Α.	General Setup and Checks: 1. Diagnostic Tests Pass and Instrument in Re 2. Diagnostic Tests Pass and Instrument in Re 3. Diagnostic Tests Pass and Instrument in Re 4. Diagnostic Tests Pass and Instrument in Re 5. Diagnostic Tests Pass and Instrument in Re 6. Diagnostic Tests Pass and Instrument i	e-set if necessary.
В.	Tests (Sign and attach test records): Configure simulator for the following tests (Level 1 Wet Calibration Check a. X Low AC (Level 1,C): Known Value ≤ Sim. Ser #: MP3061 Lot #: 201810 D b. X High AC (Level 1,C): Known Value ≥ Sim. Ser #: MP3069 Lot #: 201902 D	0.03 AC: <u>0.020</u> AC Exp. Date: <u>/0.24.2</u> 0 0.25 AC: 0.300 AC
	Configure dry gas standard for the following tests (2.	ads "RFI Detect" alue 0.080 AC Exp. Date: 4.5.22 Test 7 0.082
C.	Remarks/Maintenance: ANNUAL INSPE	EC770N
	rument is acceptable to be used in the field.	1
Breath	Analyst Signature NA	$\frac{6.5.20}{\text{Date}}$
Reviev	ved by	Date

02/18

Issuing Authority - CEE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006688
Location = TOXL 8164.14.00 09/16
06/05/2020 13:50

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 05620080A1 Standard Cyl #? 34 Standard Expiration? 04/05/2022 Oper No? 666666

Flow Cal. Date: 09/07/2016 Slope 698 Intercept -642832

IR Calibration Date: 09/07/2016
3um 9um

Oth Coef(*100): -10822 -23152 1st Coef(*100): 262468 134492 2nd Coef(*100): 3796 1436 H2O adj(mg/l*10k): 665 568

******* Printer Test End *******

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/05/2020

Alcohol Analyzer SN 80-006688 8164.14.00 09/16 13:53

WET CAL CHECK

Test		AC	AC Time	
01	Room Air	0.000	13:53	
02	Std. Sol.	0.019	13:54	
03	Room Air	0.000	13:55	
04	Std. Sol.	0.019	13:55	
05	Room Air	0.000	13:56	
06	Std. Sol.	0.019	13:56	
07	Room Air	0.000	13:57	

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = MP3061 Std Sol No = 201810D

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

LOW AC

Remarks:

0.020 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/05/2020

Alcohol Analyzer SN 80-006688 8164.14.00 09/16 13:58

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:59
02 Std. Sol.	0.295	13:59
03 Room Air	0.000	14:00
04 Std. Sol.	0.298	14:01
05 Room Air	0.000	14:01
06 Std. Sol.	0.298	14:02
07 Room Air	0.000	14:02

08 Sim Temp = 34.0°C

Simul Ser No = MP3069 Std Sol No = 201902D

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

416H M

Remarks:

0.300AC

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006688
Location = TOXL 8164.14.00 09/16
06/05/2020 14:04

Test	AC Time	
01 Room Air	0.000	14:05
02 *Subject Test	INT*	14:05
03 Room Air	0.000	14:06

*Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT

Dr. Lic. = ND/DIS921456

Lot No = 05620080A1

Cyl No = 34

Expiration Date = 04/05/2022

County = 08 Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006688 Location = TOXL 8164.14.00 09/16 06/05/2020

Test	AC	Time
01 Diagnostic	OK	14:09
02 Room Air	0.000	14:10
03 *Subject Test	RFI*	14:10
04 Room Air	0.000	14:10

14:08

*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = FemaleWeight = 150Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS921456

Lot No = 05620080A1

Cyl No = 34

Expiration Date = 04/05/2022

County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

> Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/05/2020

Alcohol Analyzer SN 80-006688 8164.14.00 09/16 14:11

DRY CAL CHECK

Test	AC Time	
01 Room Air	0.000	14:11
02 Std. Gas	0.082	14:12
03 Room Air	0.000	14:12
04 Std. Gas	0.082	14:13
05 Room Air	0.000	14:13
06 Std. Gas	0.082	14:13
07 Room Air	0.000	14:14

Lot No = 05620080A1

Cyl No = 34

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

ALIBRATION CHECK

0.080 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/05/2020

Alcohol Analyzer SN 80-006688 8164.14.00 09/16 14:14

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:15
02 Std. Gas	0.082	14:15
03 Room Air	0.000	14:16
04 Std. Gas	0.082	14:16
05 Room Air	0.000	14:17
06 Std. Gas	0.082	14:17
07 Room Air	0.000	14:17

Lot No = 05620080A1

Cyl No = 34

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks: CALIBRATION CHECK

Form 106-I8000

0.080AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/05/2020

Alcohol Analyzer SN 80-006688 8164.14.00 09/16 14:18

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	14:18
02	Std.	Gas	0.082	14:19
03	Room	Air	0.000	14:19
04	Std.	Gas	0.082	14:20
05	Room	Air	0.000	14:20
06	Std.	Gas	0.082	14:20
07	Room	Air	0.000	14:21

Lot No = 05620080A1

Cyl No = 34

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

CALIBRATION CHECK