

FLD 1079



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
SFN 59281 (06/2018)

Serial Number 80-006689	Instrument Location EMMONS COUNTY SD
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests: **CEE RECEIVED ALL TEST RECORDS. 6.12.20**
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature 	Date 5/29/20
-------------------------------	------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By 	Certified Date 6.12.2020
---------------------------	------------------------------------

NOTE: INSTALLATION DATE FOR THIS INTOXILYZER INSTRUMENT WILL BE 6.8.2020. THIS IS DUE TO A PROPER PRINT TEST PERFORMED ON 6.8.2020 AND SIGNED BY THE FIELD INSPECTOR.

6.12.2020


Intoxilyzer Test Record and Checklist
 NDDAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
 North Dakota Model 8000 SN 80-006689
 Location = LINT 8164.14.00 09/16
 05/29/2020 15:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:51
02 Std. Gas	0.077	15:51
03 Room Air	0.000	15:52
04 Std. Gas	0.077	15:52
05 Room Air	0.000	15:53
06 Std. Gas	0.077	15:53
07 Room Air	0.000	15:53

Lot No = 34418080A2
 Cyl No = 2
 Exp Date = 02/05/2021
 County = 15 Oper No. = 130117


 Operator Signature
 MARK PEARSON

Remarks:

Intoxilyzer Test Record and Checklist
 NDDAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
 North Dakota Model 8000 SN 80-006689
 Location = LINT 8164.14.00 09/16
 05/29/2020 15:54

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	15:55
02 Room Air	0.000	15:55

*Invalid Test
 Inhibited - RFI

Lot No = 34418080A2
 Cyl No = 2
 Exp Date = 02/05/2021
 County = 15 Oper No. = 130117


 Operator Signature
 MARK PEARSON

Remarks:

Form 106-18000

Form 106-18000

Intoxilyzer Test Record and Checklist
NDDAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006689
Location = LINT 8164.14.00 09/16
06/08/2020 13:10

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstuwx
ABCDEFGHIJKLMNQRSTUWX

32 characters per line
abcdefghijklmnopqrstuwx123456
ABCDEFGHIJKLMNQRSTUWXYZ123456

42 characters per line
abcdefghijklmnopqrstuwx1234567890abcdef
ABCDEFGHIJKLMNQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Menfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 34418080A2
Standard Cyl #? 2
Standard Expiration? 02/05/2021
Oper No? 130117

Flow Cal. Date: 05/29/2020
Slope 643
Intercept -576193

IR Calibration Date: 05/29/2020

	3um	9um
0th Coef(*100):	-25870	-28277
1st Coef(*100):	287704	139899
2nd Coef(*100):	26	744
H2O adj(mg/l*10k):	293	328

***** Test End *****


Operator Signature
MARK PEARSON
SN 80-006689