

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 006 5 1 2	Instrument Location Carring P.	D. (CARR)		
Reason for Install/Repair				
☐ Install After Receiving From Crime Laboratory ☐ Install After	er Location Change			
Other (Specify)				
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
☐ 3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
7 6. Run Tests:				
A. Print Test (Level1, Function P).				
☑ B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Dur	ring Test).		
7. Repair and/or Maintenance Performed (if any):				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50	0496, Form 120-G) and Place it by the In	toxilyzer® for Use.		
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Ager	ncy.		
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair Ch	neckout (SFN59281, Form		
Field Inspector Signature Mala-Obs- 4381		Date 6-10-2020		
Crime Laboratory Use Only				
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on fill County of Burleigh, North Dakota, is certified to be a true and correct correct.	e at the Office of Attorney General, Crime			
Reviewed/Certified By	El I	Certified Date 6 - 11-2020		

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006512
Location = CARR 8164.14.00 09/16
06/10/2020 17:35

********* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 2411908CA1 13 Standard Expiration? 11/05/2021 Oper No? 131838

Flow Cal. Date: 08/20/2015 Slope 680 Intercept -660856

IR Calibration Date: 08/21/2015

3um 9um

Oth Coef(*100): -23051 -19261

1st Coef(*100): 266836 137135

2nd Coef(*100): 3974 1466

H2O adj(mg/l*10k): 785 487

******* Printer Test End *******

Operator Signature NICKOLAS OSTER

Remarks: Point Test

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = CARR
06/10/2020

Alcohol Analyzer SN 80-006512 8164.14.00 09/16 17:35

DRY CAL CHECK

Test			AC	Time
01	Rocm	Air	0.000	17:36
02	Std.	Gas	0.078	17:36
03	Room	Air	0.000	17:37
04	Std.	Gas	0.077	17:37
05	Room	Air	0.000	17:38
06	Std.	Gas	0.077	17:38
07	Room .	Air	0.000	17:39

Lot No = 24119080A1

Cyl No = 13

Exp Date = 11/05/2021

County = 16

Oper No. = 131838

Operator Signature NICKOLAS OSTER

Remarks:

ANA Test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = CARR 06/10/2020

Alcohol Analyzer SN 80-006512 8164.14.00 09/16 17:40

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	17:41
02 Room Air	RFI*	17:41

*Invalid Test Inhibited - RFI

Lot No = 24119080A1

Cyl No = 13

Exp Date = 11/05/2021

County = 16

Oper No. = 131838

Operator Signature NICKOLAS OSTER

Remarks:

RFI Test

Form 106-I8C00