Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

Intoxilyzer[®] 8000 Serial Number: 80-006686

Location: ToXL

BrW-004

- General Setup and Checks: Α.
 - X Diagnostic Tests Pass and Instrument in Ready Mode 1.
 - 2. Breath tube heated
 - Date, time and location code (Level 2,E). Re-set if necessary. 3.
 - Print test (Level 1,P). Sign and attach test record. 4.
 - Tank monitor (Level 3,D,G). Display and Regulator within 50 psi. 5.
 - Display:864 psi Regulator:875 psi
- Tests (Sign and attach test records): Β.

 \mathbf{X} Configure simulator for the following tests (Level 1,S).

- Wet Calibration Check
 - X Low AC (Level 1,C): Known Value ≤ 0.03 AC: 0,030 AC a. Sim. Ser #: MP3063 Lot #: 20130 Exp. Date: 3.10.22
 - High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC b. Sim. Ser #: MP3069 Lot #:201902 D Exp. Date: 2.20.21

 \heartsuit Configure dry gas standard for the following tests (Level 1,S).

- Interferent Check (Level 1,B): Known Value: D-10 AC + 0.05% Acetone 2. Sim. Ser #: <u>PR3451</u> Lot #: <u>FCS6</u> Exp. Date: <u>NA</u> Display reads "Interferent Detect"
- 3.
- 4. Gas Cyl Lot #: 05620080 A1 Cyl #: 34 Exp. Date: 4.5-22 Test 1 $\underline{0.080}$ Test 4 $\underline{0.081}$ Test 7 $\underline{0.081}$ Test 2 $\underline{0.080}$ Test 5 $\underline{0.081}$ Test 8 $\underline{0.081}$

	.000	10310	0 -0/	10010	0 01	
Test 3 0	.080	Test 6	0.091	Test 9	0.080	
Average	0.081	_				

C.

Remarks/Maintenance: <u>ANNUAL FNSPECTION</u> TANK SENSOR CALIBRA

Instrument is acceptable to be used in the field.

Breath Analyst Signature

Reviewed by

<u>19/2026</u>

Date

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006686Location = TOXL8164.14.00 09/1606/29/202013:50

abcdefghijklmnopqrstuvwxyz1234567890-= | ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

<pre># of Print Copies? Select Std (D/W/I)? Standard Value? Standard Lot #? Standard Cyl #? Standard Expiration? Oper No? Flow Cal. Date: Slope Intercept</pre>	Enabled DABACABA Yes Yes No No Yes 1 Dry 0.080 05620080A1 34 04/05/2022 666666 09/07/2016 770 -749410
IR Calibration Date: 3um	09/06/2016 9um
Oth Coef(*100): -145 1st Coef(*100): 2522 2nd Coef(*100): 3527 H2O adj(mg/l*10k): 522	

********** Printer Test End **********

call

Operator Signature CHARLES EDER

RINT

Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006686Location = TOXL8164.14.00 09/1606/29/202013:54

	WET CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Sol. 03 Room Air 04 Std. Sol. 05 Room Air 06 Std. Sol. 07 Room Air	0.000 0.029 0.000 0.029 0.000 0.029 0.000	13:55 13:55 13:56 13:57 13:57 13:58 13:58

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = MP3063 Std Sol No = 20130 County = 08 Oper 3

Oper No. = 666666

la Ke

Operator Signature CHARLES EDER

LOW AC 0.030 Ac

Form 106-I8000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006686Location = TOXL8164.14.00 09/1606/29/202014:00

	WET CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Sol. 03 Room Air 04 Std. Sol. 05 Room Air 06 Std. Sol. 07 Room Air	0.000 0.295 0.000 0.296 0.000 0.296 0.000	14:00 14:01 14:02 14:02 14:03 14:04 14:04

08 Sim Temp = 34.0°C

Simul Ser No = MP3069 Std Sol No = 201902D County = 08

Oper No. = 666666

Al 10 00

Operator Signature CHARLES EDER

16H0.300 4

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer SN 80-006686 North Dakota Model 8000 Location = TOXL 8164.14.00 09/16 06/29/2020 14:05 Time Test AC 14:06 01 Room Air 0.000 14:06 02 *Subject Test INT* 14:07 03 Room Air 0.000 *Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = INTERFERENT Dr. Lic. = ND/DIS921456Lot No = 05620080A1Cyl No = 34Expiration Date = 04/05/2022County = 08Oper No. = 666666

Operator Signature CHARLES EDER

Interferent Check 0.10AC +0.05% ACETONE

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 Alcohol Analyzer CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-006686 8164.14.00 09/16 Location = TOXL 06/29/2020 14:08 AC Time Test 01 Diagnostic OK 14:09 14:09 02 Room Air 0.000 03 *Subject Test 14:10 RFI* 04 Room Air 0.000 14:10 *Invalid Test Inhibited - RFI Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = RFI CHECK Dr. Lic. = ND/DIS921456Lot No = 05620080A1Cyl No = 34Expiration Date = 04/05/2022County = 08Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

Form 106-18000

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006686 8164.14.00 09/16 Location = TOXL 14:10 06/29/2020

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.080 0.000 0.080 0.000 0.080 0.080 0.000	14:11 14:12 14:12 14:12 14:13 14:13 14:13

Lot No = 05620080A1Cyl No = 34Exp Date = 04/05/2022County = 08

Oper No. = 666666

in

Operator Signature CHARLES EDER ALIBRATION CHECK

0.080 AC

Form 106-18000

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006686 Location = TOXL 8164.14.00 09/16 06/29/2020 14:14

	DRY CAL CHECK	
Test	AC	Time
01 Room Ai 02 Std. Ga 03 Room Ai 04 Std. Ga 05 Room Ai 06 Std. Ga 07 Room Ai	s 0.081 r 0.000 s 0.081 r 0.000 s 0.081	14:14 14:14 14:15 14:15 14:16 14:16 14:17

Lot No = 05620080A1Cyl No = 34Exp Date = 04/05/2022County = 08

Oper No. = 666666

Operator Signature

ALIBRATION CHECK

0.080 AC

Form 106-I8000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006686Location = TOXL8164.14.00 09/1606/29/202014:17

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.081 0.000 0.080 0.000	14:18 14:19 14:19 14:19 14:20 14:20

Lot No = 05620080A1Cyl No = 34 Exp Date = 04/05/2022County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

NIECK 0 ALIBRATIC 0.080AC

Form 106-18000