

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location Oakes PD				
Reason for Install/Repair					
Install After Receiving From Crime Laboratory					
Conter (Specify)					
Oberly Mitter Deser					
Check When Done:					
1. Surge Protector Installed/Property Grounded.					
2. Telephone Line Connected to Intoxilyzer® 8000.					
☑ 3. Breath Tube Heated.					
🗹 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).					
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).					
6. Run Tests:					
A. Print Test (Level1, Function P).					
B. ACA Test (Level 1, Function C).					
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).					
7. Repair and/or Maintenance Performed (if any):					
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.					
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.					
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.					

Field Inspector Signature	Date	1
Must Chow	8/3	1/20
		/

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	Charla Etch	Certified Date
		-

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006680 Location = OKPD 8164.14.00 09/16 08/31/2020 16:19 abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34418080A2 Standard Cyl #? 50 Standard Expiration? 02/05/2021 Oper No? 130690 Flow Cal. Date: 07/20/2020 Slope 659 Intercept -627239 IR Calibration Date: 07/20/2020 3um 9um 0th Coef(\*100): -19714 -24116 lst Coef(\*100): 276675 2nd Coef(\*100): 822 136189 883 H2O adj(mg/l\*10k): 491 529 \*\*\*\*\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*\*\*\*

Operator Signature MATTHEW O'BRIEN

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006680Location = OKPD8164.14.00 09/1608/31/202016:21

			DRY CAL CHECK	
Te	est		AC	Time
01	Room	Air	0.000	16:22
02	Std.	Gas	0.081	16:22
03	Room	Air	0.000	16:23
04	Std.	Gas	0.081	16:23
05	Room	Air	0.000	16:24
06	Std.	Gas	0.080	16:24
07	Room	Air	0.000	16:24

Operator Signature MATTHEW O'BRIEN

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006680 Location = OKPD 8164.14.00 09/16 08/31/2020 16:26 Test AC Time 01 Room Air RFI\* 16:29 02 Room Air 0.000 16:29 \*Invalid Test Inhibited - RFI Sub Name = RFI TEST, INSTALL NA Sub DOB = 01/01/2001Sub Sex = Unknown Weight = NATest = OTHCit = NA Dr. Lic. = ND/NALot No = 34418080A2Cyl No = 50Expiration Date = 02/05/2021County = 11Oper No. = 130690

Operator Signature MATTHEW O'BRIEN

Remarks:

Form 106-18000