

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-006671

Location: TOXL

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 951 psi Regulator: 950 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check

- a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.015 AC  
Sim. Ser #: MP5319 Lot #: 201911C Exp. Date: 11-19-21
- b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.900 AC  
Sim. Ser #: MP3062 Lot #: 202103E Exp. Date: 3-24-23

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR3451 Lot #: ICS 6 Exp. Date: NA

Display reads "Interferent Detect"

3.  RFI Check (CMS Mode)  Display reads "RFI Detect"

4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC

Gas Cyl Lot #: 24119080A1 Cyl #: 5 Exp. Date: 11-5-21

Test 1 0.083 Test 4 0.083 Test 7 0.082

Test 2 0.083 Test 5 0.082 Test 8 0.082

Test 3 0.082 Test 6 0.082 Test 9 0.082

Average 0.082

C. Remarks/Maintenance: ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

Charles E. Ed  
Breath Analyst Signature

4.14.21  
Date

Kali G. Diet  
Reviewed by

4-19-21  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000       SN 80-006671  
Location = TOXL               8164.14.00 09/16  
04/14/2021                       10:27

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:           Enabled  
Start Test Sequence:       DABACABA  
Display Prelim Rslt?       Yes  
Display Third Digit?       Yes  
Inhib Printer(Y/N)?       No  
Display Volume?            No  
Disable On Memfull?       Yes  
# of Print Copies?         1  
Select Std (D/W/I)?       Dry  
Standard Value?            0.080  
Standard Lot #?            02621080A1  
Standard Cyl #?            29  
Standard Expiration?       03/05/2023  
Oper No?                   666666

Flow Cal. Date:            08/22/2016  
Slope                       703  
Intercept                   -691868

IR Calibration Date:       08/22/2016

	3um	9um
0th Coef(*100):	-17135	-24345
1st Coef(*100):	254798	132856
2nd Coef(*100):	3386	1354
H2O adj(mg/l*10k):	610	539

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006671  
Location = TOXL      8164.14.00 09/16  
04/14/2021      10:29

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:29
02 Std. Sol.	0.012	10:30
03 Room Air	0.000	10:30
04 Std. Sol.	0.014	10:31
05 Room Air	0.000	10:32
06 Std. Sol.	0.014	10:32
07 Room Air	0.000	10:33

08 Sim Temp = 34.0°C

Simul Ser No = MP5319  
Std Sol No = 201911C  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Low AC*  
*0.015 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006671  
Location = TOXL                    8164.14.00 09/16  
04/14/2021                            10:38

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:39
02 Std. Sol.	0.393	10:39
03 Room Air	0.000	10:40
04 Std. Sol.	0.396	10:41
05 Room Air	0.000	10:41
06 Std. Sol.	0.397	10:42
07 Room Air	0.000	10:43

08 Sim Temp = 34.0°C

Simul Ser No = MP3062  
Std Sol No = 202103E  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *HIGH AC*  
*0.400 AC*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006671  
Location = TOXL                8164.14.00 09/16  
04/14/2021                        10:43

Test	AC	Time
01 Room Air	0.000	10:44
02 *Subject Test	INT*	10:44
03 Room Air	0.000	10:45

\*Invalid Test  
Interferent Detected

Sub Name = BILLS, THE BUFFALO  
Sub DOB = 06/14/1993  
Sub Sex = Male                    Weight = NA  
Test = DUI                        Cit = INTERFERENT  
Dr. Lic. = ND/BIL933277  
Lot No = 24119080A1  
Cyl No = 5  
Expiration Date = 11/05/2021  
County = 08                        Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Interferent Check*  
*0.10 AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-006671  
Location = TOXL                              8164.14.00 09/16  
04/14/2021    10:49

Test	AC	Time
01 Diagnostic	OK	10:50
02 Room Air	0.000	10:50
03 *Subject Test	RFI*	10:50
04 Room Air	0.000	10:51

\*Invalid Test  
Inhibited - RFI

Sub Name = BILLS, THE BUFFALO  
Sub DOB = 06/14/1993  
Sub Sex = Male                              Weight = NA  
Test = DUI                                      Cit = RFI CHECK  
Dr. Lic. = ND/BIL933277  
Lot No = 24119080A1  
Cyl No = 5  
Expiration Date = 11/05/2021  
County = 08                                      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *RFI CHECK*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006671  
Location = TOXL      8164.14.00 09/16  
04/14/2021      10:53

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:53
02 Std. Gas	0.083	10:54
03 Room Air	0.000	10:54
04 Std. Gas	0.083	10:55
05 Room Air	0.000	10:55
06 Std. Gas	0.082	10:55
07 Room Air	0.000	10:56

Lot No = 24119080A1  
Cyl No = 5  
Exp Date = 11/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006671  
Location = TOXL      8164.14.00 09/16  
04/14/2021      10:56

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:57
02 Std. Gas	0.083	10:57
03 Room Air	0.000	10:58
04 Std. Gas	0.082	10:58
05 Room Air	0.000	10:58
06 Std. Gas	0.082	10:59
07 Room Air	0.000	10:59

Lot No = 24119080A1  
Cyl No = 5  
Exp Date = 11/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006671  
Location = TOXL      8164.14.00 09/16  
04/14/2021      10:59

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:00
02 Std. Gas	0.082	11:00
03 Room Air	0.000	11:01
04 Std. Gas	0.082	11:01
05 Room Air	0.000	11:02
06 Std. Gas	0.082	11:02
07 Room Air	0.000	11:03

Lot No = 24119080A1  
Cyl No = 5  
Exp Date = 11/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000