



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <i>80-006689</i>	Instrument Location <i>Linton</i>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i>	Date <i>4/15/21</i>
---	------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <i>[Signature]</i>	Certified Date <i>4/23/2021</i>
---	------------------------------------

**Intoxilizer Test Record and Checklist**  
 ND06G Crime Lab. Div., Bismarck, ND 58501

CH1, Inc. Intoxilizer Alcohol Analyzer  
 North Dakota Model 8000 SN 80-006689  
 Location = LINT 8164.14.00 09/16  
 04/15/2021 13:14

\*\*\*\*\* Printer Test \*\*\*\*\*

12 char/line  
 abcdefghijkl  
 ABCDEFGHIJKL

16 char per line  
 abcdefghijklmnop  
 ABCDEFGHIJKLMNOP

21 chars per line  
 abcdefghijklmnopqrstu  
 ABCDEFGHIJKLMNOPQRSTU

24 characters per line  
 abcdefghijklmnopqrstuvw  
 ABCDEFGHIJKLMNOPQRSTUWVX

32 characters per line  
 abcdefghijklmnopqrstuvwxyz123456  
 ABCDEFGHIJKLMNOPQRSTUVWXYZ123456

42 characters per line  
 abcdefghijklmnopqrstuvwxyz1234567890abcdef  
 ABCDEFGHIJKLMNOPQRSTUVWXYZ1234567890BCDEF

Current Instrument Setup Enabled  
 Data Entry Mode: 09884C994  
 Start Test Sequence: Yes  
 Display Prelim Result? Yes  
 Display Third Digit? Yes  
 Inhib Printer(Y/N)? No  
 Display Volume? No  
 Disable On Refill? Yes  
 # of Print Copies? 1  
 Select Std (D/W/L)? Dry  
 Standard Value? 0.060  
 Standard Lot #? 056208091  
 Standard Cyl #? 29  
 Standard Expiration? 04/05/2022  
 Oper No? 131117

Flow Cal. Date: 05/29/2020  
 Slope: 643  
 Intercept: -576193

IR Calibration Date: 05/29/2020  
 Sun Sun

0th Coef(\*100): -25870 -28277  
 1st Coef(\*100): 287704 139899  
 2nd Coef(\*100): 26 744  
 H2O adj(tng/\*100): 293 328

\*\*\*\*\* Test End \*\*\*\*\*

  
 Operator Signature  
 MARK PEARSON

Remarks:

Form 106-18000

**Intoxilizer Test Record and Checklist**  
 ND06G Crime Lab. Div., Bismarck, ND 58501

CH1, Inc. Intoxilizer Alcohol Analyzer  
 North Dakota Model 8000 SN 80-006689  
 Location = LINT 8164.14.00 09/16  
 04/15/2021 13:15

DRY CAL CHECK

Test RC Time

01 Room Air	0.000	13:15
02 Std. Gas	0.078	13:16
03 Room Air	0.000	13:16
04 Std. Gas	0.077	13:16
05 Room Air	0.000	13:17
06 Std. Gas	0.078	13:17
07 Room Air	0.000	13:18

Lot No = 056208091  
 Cyl No = 29  
 Exp Date = 04/05/2022  
 County = 15  
 Oper No. = 130117

  
 Operator Signature  
 MARK PEARSON

Remarks:

Form 106-18000

**Intoxilizer Test Record and Checklist**  
 ND06G Crime Lab. Div., Bismarck, ND 58501

CH1, Inc. Intoxilizer Alcohol Analyzer  
 North Dakota Model 8000 SN 80-006689  
 Location = LINT 8164.14.00 09/16  
 04/15/2021 13:19

Test RC Time

01 Room Air	RF1*	13:21
02 Room Air	RF1*	13:21

\*Invalid Test  
 Inhibited - RF1

Sub Name = PEARSON, MARK ALAN  
 Sub DOB = 02/28/1962 Weight = 280  
 Sub Sex = Male Cit = NH  
 Test = OTH  
 Dr. Lic. = ND/PE629285  
 Lot No = 056208091

Cyl No = 29  
 Expiration Date = 04/05/2022  
 County = 15  
 Oper No. = 130117

  
 Operator Signature  
 MARK PEARSON

Remarks:

Form 106-18000