



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <b>80-007087</b>	Instrument Location <b>Wishek PD</b>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000. - *connected when crime lab requests*
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i> <b>Lucas Kuntz 1664</b>	Date <b>4-23-21</b>
----------------------------------------------------------------------------	------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <i>[Signature]</i>	Certified Date <b>5.4.2021</b>
---------------------------------------------	-----------------------------------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-007087  
Location = WISH                8164.16.00 09/18  
04/23/2021                      14:11

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup


Data Entry Mode:                Enabled  
Start Test Sequence:            DABACABA  
Display Prelim Rslt?            Yes  
Display Third Digit?            Yes  
Inhib Printer(Y/N)?            No  
Display Volume?                No  
Disable On Memfull?            Yes  
# of Print Copies?              1  
Select Std (D/W/I)?            Dry  
Standard Value?                0.080  
Standard Lot #?                02621080A1  
Standard Cyl #?                21  
Standard Expiration?           03/05/2023  
Oper No?                        131099

Flow Cal. Date:                07/02/2019  
Slope                            721  
Intercept                       -792109

IR Calibration Date:            07/02/2019  
                                  3um                            9um

	3um	9um
0th Coef(*100):	-45706	-18865
1st Coef(*100):	253208	131287
2nd Coef(*100):	2874	1457
H2O adj(mg/l*10k):	1171	483

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
LUCAS KUNTZ

Remarks: *Print test*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-007087  
Location = WISH      8164.16.00 09/18  
04/23/2021      14:11

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:11
02 Std. Gas	0.080	14:12
03 Room Air	0.000	14:12
04 Std. Gas	0.080	14:13
05 Room Air	0.000	14:13
06 Std. Gas	0.080	14:13
07 Room Air	0.000	14:14

Lot No = 02621080A1  
Cyl No = 21  
Exp Date = 03/05/2023  
County = 26      Oper No. = 131099

  
Operator Signature  
LUCAS KUNTZ

Remarks: *reinstall cal check*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-007087  
Location = WISH      8164.16.00 09/18  
04/23/2021      14:15

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	14:16
02 Room Air	0.000	14:16

\*Invalid Test  
Inhibited - RFI

Lot No = 02621080A1  
Cyl No = 21  
Exp Date = 03/05/2023  
County = 26      Oper No. = 131099

  
Operator Signature  
LUCAS KUNTZ

Remarks:

RFI test

Form 106-I8000