

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number In	strument Location Lake V	PEVL
Reason for Install/Repair		
☑ Install After Receiving From Crime Laboratory ☐ Install After L	ocation Change	
Other (Specify)		
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
✓ 2. Telephone Line Connected to Intoxilyzer® 8000.		
✓ 3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).		
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
☑ A. Print Test (Level1, Function P).		
☑ B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode or Le	vel 1, Function B or C; Key Radio During	Test).
7. Repair and/or Maintenance Performed (if any): none		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5049)	96, Form 120-G) and Place it by the Intoxil	lyzer® for Use.
☑ 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.		
10. Send the Following to the Crime Laboratory: Completed Intoxil 104-G), Print Test, ACA Test, and RFI Test.	lyzer® 8000 Installation and Repair Check	kout (SFN59281, Form
Field Inspector Signature	Daţ	e
The state of the s	C	28/2021
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to be u	•	

the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	Certified Date 5/21/21	
	7 7	

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = DEVL 8164.14.00 09/16
04/15/2021 19:59

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 02621080A1 Standard Cyl #? 29 Standard Expiration? 03/05/2023 Oper No? 132008

Flow Cal. Date: 08/22/2016 Slope 703 Intercept -691868

H2O adj(mg/l*10k): 610

IR Calibration Date: 08/22/2016
3um 9um

Oth Coef(*100): -17135 -24345
1st Coef(*100): 254798 132856
2nd Coef(*100): 3386 1354

******* Printer Test End ********

Operator Signature

WILLIAM EARL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = DEVL 04/15/2021

Alcohol Analyzer SN 80-006671 8164.14.00 09/16 19:53

DRY CAL CHECK

Test		AC	Time
02 St	om Air 1. Gas	0.000 0.081 0.000 0.081 0.000 0.081 0.000	19:54 19:54 19:55 19:55 19:56 19:56

Lot No = 02621080A1

Cyl No = 29

Exp Date = 03/05/2023

County = 36

Oper No. = 132008

Operator Signature WILLIAM EARL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = DEVL 04/15/2021

Alcohol Analyzer SN 80-006671 8164.14.00 09/16 19:50

DRY CAL CHECK

AC	Time
0.000 0.081 0.000 RFI*	19:51 19:51 19:52 19:52
	0.000 0.081 0.000

*Invalid Test Inhibited - RFI

Lot No = 02621080A1

Cyl No = 29

Exp Date = 03/05/2023

County = 36

Oper No. = 132008

Operator Signature WILLIAM EARL

Remarks:

Form 106-I8000



Intoxilyzer® Serial No 80-006671	umber		Location	
			Devils Lake	EVL
Gas Lot Number 02621080A1		Gas Cylinder Number 29		Gas Expiration Date 03/05/2023
Test Date	Chemical Test Operator Number		Operator's (PRINT Last Na	s Name ame, First, MI)
			988.4	

FORM 120-G

AGENCY RETAIN AT INTOXILYZER® LOCATION