

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 406673	Instrument Location	BURKE	County So Shap
Reason for Install/Repair			1 7
Install After Receiving From Crime Laboratory Install After	er Location Change		
Other (Specify)			
Check When Done:			
1. Surge Protector Installed/Property Grounded.			
2. Telephone Line Connected to Intoxilyzer® 8000.			
3. Breath Tube Heated.			
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	ion; Level 2, Function E)).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
6. Run Tests:			
A. Print Test (Level1, Function P).			
B. ACA Test (Level 1, Function C).			
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or 0	C; Key Radio Du	ring Test).
7. Repair and/or Maintenance Performed (if any):			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5)	0496, Form 120-G) and	Place it by the Ir	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G	6) at the Intoxilyzer® Loc	cation at the Age	ncy.
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installati	on and Repair C	heckout (SFN59281, Form
Field Inspector Signature			Date 5/6/21
Crime Laboratory Use Only			/ (
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on fi County of Burleigh, North Dakota, is certified to be a true and correct co	le at the Office of Attorn	ey General, Crim	
Reviewed/Certified By	50		Certified Date

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006673
Location = BURK 8164.14.00 09/16
05/05/2021 16:59

******* Printer Test *******

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Display Prelim Rslt? Yes
Display Thing Thi Display Third Digit? No Inhib Printer(Y/N)? Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 05620080A1 Standard Cyl #? 33 Standard Expiration? 04/05/2022 Oper No? 120804

Flow Cal. Date: 03/30/2017 Slope 697 Intercept -667725

IR Calibration Date: 03/30/2017

3um 9um

Oth Coef(*100): -9136 -17613
1st Coef(*100): 267126 133333

1st Coef(*100): 267126 13333 2nd Coef(*100): 5438 1698 H2O adj(mg/l*10k): 684 589

****** Printer Test End *******

Operator Signature
JASON CARTIER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BURK
05/05/2021

Alcohol Analyzer SN 80-006673 8164.14.00 09/16 16:59

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:00
02 Std. Gas	0.082	17:01
03 Room Air	0.000	17:01
04 Std. Gas	0.082	17:01
05 Room Air	0.000	17:02
06 Std. Gas	0.082	17:02
07 Room Air	0.000	17:03

Lot No = 05620080A1

Cyl No = 33

Exp Date = 04/05/2022

County = 07

Remarks:

Oper No. = 120804

JASON CARTIER

Form 106-I8000

Operator Signature

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BURK
05/05/2021

Alcohol Analyzer SN 80-006673 8164.14.00 09/16 17:04

DRY CAL CHECK

T^{ϵ}	est		AC	Time
01	Room	Air	RFI*	17:04
02	Room	Air	RFI*	17:05

*Invalid Test Inhibited - RFI

Lot No = 05620080A1

Cyl No = 33

Exp Date = 04/05/2022

County = 07

Oper No. = 120804

Operator Signature
JASON CARTIER

Remarks:

Form 106-I8000