

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 30 - 00 4513 Instrument Location Bowns (	County Sheriff ofc.
Reason for Install/Repair	/
Install After Receiving From Crime Laboratory  Install After Location Change	
Other (Specify)	
Check When Done:	
1. Surge Protector Installed/Property Grounded.	
2. Telephone Line Connected to Intoxilyzer® 8000.	
3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).	·
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
(A. Run Tests:	
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key F	Radio During Test).
7. Repair and/or Maintenance Performed (If any):	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it	by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at	the Agency.
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and I 104-G), Print Test, ACA Test, and RFI Test.	Repair Checkout (SFN59281, Form
Field Inspector Signature	Date 5-17-7021
Crime Laboratory Use Only	
This installation has been reviewed and the instrument is approved to be used for the analysis of breath the date the Field Inspector performed the installation. This record on file at the Office of Attorney Gene County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.	
Reviewed/Certified By	Certified Date

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BOWM
05/17/2021

Alcohol Analyzer
SN 80-006513
8164.13.00 06/09
13:59

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup Data Entry Mode: Start Test Sequence: Display Prelim Rslt? Display Third Digit? Inhib Printer(Y/N)? Display Volume? Disable On Memfull? # of Print Copies? Select Std (D/W/I)? Standard Value? Standard Lot #? Standard Cyl #? Standard Expiration? Oper No?	Enabled DABACABA Yes Yes No No Yes 1 Dry 0.080 05620080A1 8 04/05/2022 130015
Flow Cal. Date:	08/20/2015
Slope	680
Intercept	-696037

IR Calibration Date	9: 08/21/20 3um	015 9um
Oth Coef(*100): 1st Coef(*100): 2nd Coef(*100): H2O adj(mg/l*10k):	-14436 265490 4422 672	-13610 136147 1710 448

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature FRANK EBERLE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BOWM
05/17/2021

SN 80-006513 8164.13.00 06/09 13:46

## DRY CAL CHECK

Test		AC	Time
01 Room 02 Std. 03 Room 04 Std. 05 Room 06 Std. 07 Room	Gas Air Gas Air Gas	0.000 0.080 0.000 0.081 0.000 0.081	13:46 13:47 13:47 13:48 13:49 13:49

Lot No = 05620080A1

Cyl No = 8

Exp Date = 04/05/2022

County = 06

Oper No. = 130015

Operator Signature FRANK EBERLE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000	Alcohol Analyzer SN 80-006513
Location = BOWM	8164.13.00 06/09
05/17/2021	14:05

Test	AC	Time
01 Room Air	0.000	14:07
02 *Subject Test	RFI*	14:07
03 Room Air	0.000	14:08

\*Invalid Test Inhibited - RFI

Sub Name = SAM621234, SUBJECT A

Sub DOB = 03/18/1962

Sub Sex = Male Weight = 220
Test = OTH Cit = ABC123

Dr. Lic. = ND/SAM621234

Lot No = 05620080A1

Cyl No = 8

Expiration Date = 04/05/2022

County = 06

Oper No. = 130015

Operator Signature FRANK EBERLE

Remarks:

Form 106-I8000