Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

Location: 70×L Intoxilyzer[®] 8000 Serial Number: 80-006672

- A. General Setup and Checks:
 - 1. Diagnostic Tests Pass and Instrument in Ready Mode
 - 2. Breath tube heated
 - 3. A Date, time and location code (Level 2,E). Re-set if necessary.
 - 4. Print test (Level 1,P). Sign and attach test record.
 - 5. De Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
 - a. Display: 932 psi Regulator: 950 psi
- B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

- 1. Wet Calibration Check
 - a. \cancel{A} Low AC (Level 1,C): Known Value ≤ 0.03 AC: $\cancel{O} \cdot 020$ AC Sim. Ser #: $\cancel{MP} \cdot 3064$ Lot #: $\cancel{20070}$ Exp. Date: $\cancel{213} \cdot \cancel{22}$
 - b. We High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC Sim. Ser #: MP53/7 Lot #: 2020/2A Exp. Date: 12-8-22

Configure dry gas standard for the following tests (Level 1,S).

- 2. X Interferent Check (Level 1,B): Known Value: $\frac{0 \cdot l^0}{AC} + 0.05\%$ Acetone Sim. Ser # $\frac{0.83847}{D}$ Lot #: $\frac{2C56}{D}$ Exp. Date: $\frac{10}{AC}$
- 3. KRFI Check (CMS Mode) Cisplay reads "RFI Detect"

| 10312 01070 | | 10000 | 0 01 |
|----------------|--------------|--------|------|
| Test 3 0.078 | Test 6 0.079 | Test 9 | 0.07 |
| Average 0-0.79 | | | |

C. Remarks/Maintenance: ANNUAL FUSPECTON FLOW SENSOR ISSUE, HARD TO SUPPLY SAMPLE, FLOW SENSOR READING 56

CLASSROOM ONLY INSTRUMENT UNTIL REPAIRED.

instrument is acceptable to be used in the field. $\mathcal{N}ar{D}$

Analyst Signatøre

Date

02/18 Issuing Authority - CEE

Page 1 of 1

BrW-004

Reviewed by

R

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006672 Location = TOXL8164.14.00 09/16 05/18/2021 11:06 abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 24119080A1 Standard Cyl #? 10 11/05/2021 Standard Expiration? Oper No? 666666 Flow Cal. Date: 07/05/2017 705 Slope -687158 Intercept IR Calibration Date: 07/05/2017 3um 9um _____ 0th Coef(*100): -22113 -17955
 Oth Coef(*100):
 -22113

 1st Coef(*100):
 271833

 2nd Coef(*100):
 2129
130091 1363 H2O adj(mq/l*10k): 691568

Operator Signature CHARLES EDER

Form 106-I8000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006672Location = TOXL8164.14.00 09/1605/18/202111:11

| | | 1 |
|--------------|---------------|-------|
| | WET CAL CHECK | |
| L | | |
| | | |
| Test | AC | Time |
| 01 D | 0.000 | 11 10 |
| 01 Room Air | 0.000 | 11:12 |
| 02 Std. Sol. | 0.018 | 11:12 |
| 03 Room Air | 0.000 | 11:13 |
| 04 Std. Sol. | 0.019 | 11:14 |
| 05 Room Air | 0.000 | 11:14 |
| 06 Std. Sol. | 0.019 | 11:15 |
| 07 Room Air | 0.000 | 11:15 |
| | | |

08 Sim Temp = 34.0°C

Simul Ser No = MP3064 Std Sol No = 20070 County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Low F 0.020

Form 106-18000

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006672 Location = TOXL 8164.14.00 09/16 05/18/2021 11:22

| | WET CAL CHECK | |
|---------------|---------------|-------|
| Test | AC | Time |
| 01 Room Air | 0.000 | 11:23 |
| 02 Std. Sol. | 0.298 | 11:24 |
| 03 Room Air | 0.000 | 11:24 |
| 04 Std. Sol. | 0.298 | 11:25 |
| 05 Room Air | 0.000 | 11:26 |
| 06 Std. Sol. | 0.297 | 11:26 |
| 07 Room Air | 0.000 | 11:27 |
| 08 Sim Temp = | 34.0°C | |

Simul Ser No = MP5317 Std Sol No = 202012ACounty = 08

Oper No. = 666666

Operator Signature CHARLES EDER

HIGH AC D. 300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006672 Location = TOXL8164.14.00 09/16 05/18/2021 11:35 Test AC Time 11:36 01 Room Air 0.000 02 *Subject Test INT* 11:36 03 Room Air 0.000 11:37 *Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 01/08/1987Weight = 150Sub Sex = Male Test = DUI Cit = INTERFERNT Dr. Lic. = ND/DIS871456Lot No = 02621080A1Cyl No = 25Expiration Date = 03/05/2023County = 08Oper No. = 666666

Operator Signature CHARLES EDER Remarks: Interferent Check 0./0AC + 0.05% ACETONE

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006672 8164.14.00 09/16 Location = TOXL 05/18/2021 11:40 Time Test AC 01 Diagnostic 11:41 OK 11:42 02 Room Air 0.000 03 *Subject Test RFI* 11:42 11:42 04 Room Air 0.000 *Invalid Test Inhibited - RFI Sub Name = DISCOVER, THE SPIRIT Sub DOB = 01/08/1987Sub Sex = Male Weight = 150Cit = RFI CHECK Test = DUI Dr. Lic. = ND/DIS871456Lot No = 02621080A1Cyl No = 25Expiration Date = 03/05/2023Oper No. = 666666County = 08

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

Form 106-18000

| CMI, Inc. Intoxilyzer | Alcohol Analyzer |
|-------------------------|------------------|
| North Dakota Model 8000 | SN 80-006672 |
| Location = TOXL | 8164.14.00 09/16 |
| 05/18/2021 | 12:12 |

| | DRY CAL CHECK | |
|---|---|---|
| Test | AC | Time |
| 01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air | 0.000 0.079 0.000 0.078 0.000 0.078 0.000 | $12:13 \\ 12:13 \\ 12:14 \\ 12:14 \\ 12:15 \\ 12:1$ |

Lot No = 02621080A1 Cyl No = 25 Exp Date = 03/05/2023 County = 08

Oper No. = 666666

CAU

Operator Signature CHARLES EDER

0.080 14

Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006672Location = TOXL8164.14.00 09/1605/18/202112:26

| | DRY CAL CHECK | |
|-------------|---------------|-------|
| Test | AC | Time |
| 01 Room Air | 0.000 | 12:26 |
| 02 Std. Gas | 0.079 | 12:27 |
| 03 Room Air | 0.000 | 12:27 |
| 04 Std. Gas | 0.079 | 12:27 |
| 05 Room Air | 0.000 | 12:28 |
| 06 Std. Gas | 0.079 | 12:28 |
| 07 Room Air | 0.000 | 12:29 |

Lot No = 02621080A1 Cyl No = 25 Exp Date = 03/05/2023 County = 08

Oper No. = 666666

CAECK

Operator Signature CHARLES EDER

Remarks:

Form 106-18000

0.080AC

AUBRAT

Alcohol Analyzer CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-006672 8164.14.00 09/16 Location = TOXL 05/18/2021 12:29

| | DRY CAL CHECK | |
|---|---|--|
| Test | AC | Time |
| 01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air | 0.000 0.080 0.000 0.079 0.000 0.078 0.000 | 12:30 12:30 12:31 12:31 12:32 12:32 |

Lot No = 02621080A1Cyl No = 25Exp Date = 03/05/2023County = 08

Oper No. = 666666

al

Operator Signature CHARLES EDER

HECLE 0.080 AC

Form 106-18000