

ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

Chemical Test Operator Name (Print) MOEN NATHAN	3		
Location GFCC C		Intoxilyzer® Serial Nu	
Gas Lot Number 0 2 6 2) 6 8 0 A 1	Gas Cylinder Number		Gas Expiration Date

1	Ch	ack	\A/h	on 1	Done	

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.08 AC
 - B. 0.079 AC
 - C. 0.080 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature	Date ,
nother	6/4/21
Reviewed By (Crime Laboratory Use Only)	Date 6.4.2021
(sum to the	0,000

FORM 105-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006507 Location = GFCC 8164.14.00 09/16 06/04/2021

09:59

DRY CAL CHECK

$T\epsilon$	est		AC	Time
01	Room	Air	0.000	10:00
02	Std.	Gas	0.081	10:00
03	Room	Air	0.000	10:00
04	Std.	Gas	0.079	10:01
05	Room	Air	0.000	10:01
06	Std.	Gas	0.080	10:02
07	Room	Air	0.000	10:02

Lot No = 02621080A1

Cyl No = 16

Exp Date = 03/05/2023

County = 18

Oper No. = 100452

Operator Signature NATHAN MOEN

Remarks:

NO PRI DETECTED

Form 106-I8000