



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <i>80-007086</i>	Instrument Location <i>Lisbon</i>
Reason for Install/Repair	
<input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i>	Date <i>6/14/21</i>
---	------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <i>[Signature]</i>	Certified Date <i>6.17.21</i>
---	----------------------------------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-007086  
Location = LISB                    8164.16.00 09/18  
06/14/2021                            17:43

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890+|=|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            05620080A1  
Standard Cyl #?            44  
Standard Expiration?        04/05/2022  
Oper No?                    020607

Flow Cal. Date:            10/16/2018  
Slope                        688  
Intercept                    -653049

IR Calibration Date:        10/16/2018  
                                  3um                    9um

-----  
0th Coef(\*100):            -47466                    -23263  
1st Coef(\*100):            233222                    130299  
2nd Coef(\*100):            3639                      1618  
H2O adj(mg/l\*10k):        1280                      608

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
TIMOTHY JORDAN

Remarks:

Install


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-007086  
Location = LISB      8164.16.00 09/18  
06/14/2021      17:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:50
02 Std. Gas	0.082	17:51
03 Room Air	0.000	17:51
04 Std. Gas	0.082	17:52
05 Room Air	0.000	17:52
06 Std. Gas	0.082	17:53
07 Room Air	0.000	17:53

Lot No = 05620080A1  
Cyl No = 44  
Exp Date = 04/05/2022  
County = 37      Oper No. = 020607

  
\_\_\_\_\_  
Operator Signature  
TIMOTHY JORDAN

Remarks:

*I n stall*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-007086  
Location = LISB      8164.16.00 09/18  
06/14/2021      18:09

Test	AC	Time
01 Room Air	RFI*	18:13
02 Room Air	RFI*	18:14

\*Invalid Test  
Inhibited - RFI




Operator Signature  
TIMOTHY JORDAN

Remarks:

Install

Form 106-I8000

NOTE: SUBJECT AND CYLINDER INFORMATION  
DO NOT APPEAR ON THIS TEST RECORD.  
THIS IS NORMAL AND HOW CMI PROGRAMS  
THE INTOX 8000. THE FIELD INSPECTOR  
DID NOT ENTER A SUBJECT LAST AND/OR FIRST  
NAME DURING DATA ENTRY. THIS PRODUCES  
A RECORD SUCH AS THIS. THIS IS  
A VALID TEST RECORD DEMONSTRATING  
AN RFI.

  
6.17.21