

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number Inst 80 - 0366.76	KillOCCI Police Department		
Reason for Install/Repair			
Install After Receiving From Crime Laboratory 🔲 Install After Location Change			
☐ Other (Specify)			
Check When Done:			
1. Surge Protector Installed/Property Grounded.			
X 2. Telephone Line Connected to Intoxilyzer® 8000.			
3. Breath Tube Heated.			
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; L	Level 2, Function E).		
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
6. Run Tests:			
A. Print Test (Level1, Function P).			
B. ACA Test (Level 1, Function C).			
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level	vel 1, Function B or C; Key Radio During Test).		
7. Repair and/or Maintenance Performed (if any):	2		
Second (SFN50496,	6, Form 120-G) and Place it by the Intoxilyzer® for Use.		
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at t	the Intoxilyzer® Location at the Agency.		
10. Send the Following to the Crime Laboratory: Completed Intoxilyz 104-G), Print Test, ACA Test, and RFI Test.	vzer® 8000 Installation and Repair Checkout (SFN59281, Form		
Field Inspector Signature	Date 6-15-21		

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Certified Date 6.16.21 Reviewed/Certified By inher

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer
 North Dakota Model 8000
 SN 80-006676

 Location = KDPD
 8164.14.00 09/16
 11:22 06/15/2021 abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No DABACABA Display Volume? No Disable On Memfull? Yes 1 # of Print Copies? # OI Print Copies: 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 24119080A1 14 Standard Expiration? 11/05/2021 100444 Oper No? 08/22/2016 Flow Cal. Date: Slope 710 -713681 Intercept IR Calibration Date: 08/22/2016 3um 9um _____

 0th Coef(*100):
 -16667

 1st Coef(*100):
 259190

 2nd Coef(*100):
 3393

-22327 132291 1374 H2O adj(mg/l*10k): 611546

Operator Signature ERIC BRAATHEN

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006676Location = KDPD8164.14.00 09/1606/15/202111:23

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:24
02 Std. Gas	0.082	11:24
03 Room Air	0.000	11:24
04 Std. Gas	0.081	11:25
05 Room Air	0.000	11:25
06 Std. Gas	0.081	11:26
07 Room Air	0.000	11:26

11 Breethin

Operator Signature ERIC BRAATHEN

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006676Location = KDPD8164.14.00 09/1606/15/202111:28

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air	0.000 RFI* 0.000	11:29 11:30 11:30

*Invalid Test Inhibited - RFI

Lot No = 24119080A1 Cyl No = 14 Exp Date = 11/05/2021 County = 13 Or

Oper No. = 100444

Enc Bratha

Operator Signature ERIC BRAATHEN

Remarks:

Form 106-18000