

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 006684	Instrument Location West Fargo Police Dept.
Reason for Install/Repair	
☐ Install After Receiving From Crime Laboratory ☐ Install After	r Location Change
Other (Specify)	
Check When Done:	
1. Surge Protector Installed/Property Grounded.	
2. Telephone Line Connected to Intoxilyzer® 8000.	
☑ 3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	on; Level 2, Function E).
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
☑ 6. Run Tests;	
A. Print Test (Level1, Function P).	
☑ B. ACA Test (Level 1, Function C).	
☑ C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio During Test).
√ √ 7. Repair and/or Maintenance Performed (if any): NA	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5)	0496, Form 120-G) and Place it by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxllyzer® 8000 Installation and Repair Checkout (SFN59281, Form
Field Inspector Signature	Date
Preid Inspector digitature	06/25/2021
Crime Laboratory Use Only	
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on fi County of Burleigh, North Dakota, is certified to be a true and correct co	e at the Office of Attorney General, Crime Laboratory Division, in the
Reviewed/Certified By	Certified Date (28/202

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer	
North Dakota Model 8000	SN 80-006684	
Location = WFPD 06/25/2021	8164.14.00 09/16 16:34	

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPORSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 05620080A1 Standard Cyl #? Standard Expiration? 04/05/2022

Flow Cal. Date: 09/06/2016 Slope 709 Intercept -725441

Oper No?

IR Calibration Date: 09/06/2016 3um 0th Coef(*100): -16013 -23467 1st Coef(*100): 256431 2nd Coef(*100): 2904 135731 2904 1268 H2O adj(mg/l*10k): 553

130781

514

******** Printer Test End ********

Operator Signature KYLE JOHNSON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WFPD
06/25/2021

Alcohol Analyzer SN 80-006684 8164.14.00 09/16 16:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:36
02 Std. Gas	0.080	16:36
03 Room Air	0.000	16:37
04 Std. Gas	0.081	16:37
05 Room Air	0.000	16:37
06 Std. Gas	0.081	16:38
07 Room Air	0.000	16:38

Lot No = 05620080A1

Cyl No = 38

Exp Date = 04/05/2022

County = 09

Oper No. = 130781

Operator Signature KYLE JOHNSON

Remarks: ACA Test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WFPD
06/25/2021

Alcohol Analyzer SN 80-006684 8164.14.00 09/16 16:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	16:40
02 Room Air	0.000	16:40

*Invalid Test Inhibited - RFI

Lot No = 05620080A1

Cyl No = 38

Exp Date = 04/05/2022

County = 09

Oper No. = 130781

Operator Signature KYLE JOHNSON

Remarks:

RFI Test

Form 106-I8000