

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 006 682	Instrument Location ろここ				
Reason for Install/Repair					
🔀 Install After Receiving From Crime Laboratory 🛛 Install After Location Change					
Other (Specify)					

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- X 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).

🕅 6. Run Tests:

- A. Print Test (Level1, Function P).
- B. ACA Test (Level 1, Function C).
- 🔀 C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any):
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature	ZSM/	Date 6/25/2021

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By		Certified Date
	lind El	1 20 2071
	(Martin Carca	6.20.202

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer
 North Dakota Model 8000
 SN 80-006682

 Location = BCCC
 8164.14.00 09/16
Location = BCCC 16:23 06/25/2021 abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+? Current Instrument Setup Enabled Data Entry Mode: Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No No Display Volume? Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry 0.080 Standard Value? Standard Lot #? 05620080A1 Standard Cyl #? 45 Standard Expiration? 04/05/2022 100012 Oper No? 09/06/2016 Flow Cal. Date: 701 Slope -706193 Intercept IR Calibration Date: 09/06/2016 3um 9um -----Oth Coef(*100):-171501st Coef(*100):2717262nd Coef(*100):831 -23247 134617 1215 518 H2O adj(mg/l*10k): 393

Operator Signature DANA RUSTEBAKKE

Remarks:

Print Test

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer SN 80-006682 North Dakota Model 8000 8164.14.00 09/16 Location = BCCC 16:23 06/25/2021

		DRY CAL CHECK	
Tes	t	AC	Time
02 S 03 R 04 S 05 R	oom Air td. Gas oom Air td. Gas oom Air td. Gas	0.000 0.080 0.000 0.080 0.000 0.000 0.080	16:24 16:24 16:25 16:25 16:26 16:26
	oom Air	0.000	16:26

Lot No = 05620080A1Cyl No = 45Exp Date = 04/05/2022Oper No. = 100012County = 02

Operator Signature DANA RUSTEBAKKE

Remarks:

ALA Test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 Alcohol Analyzer CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-006682 8164.14.00 09/16 Location = BCCC16:28 06/25/2021 Time AC Test 16:30 01 Room Air 0.000 02 *Subject Test RFI* 16:30 03 Room Air 0.000 16:30 *Invalid Test Inhibited - RFI Sub Name = TEST, TEST TEST Sub DOB = 01/01/2000Weight = 175Sub Sex = Male Cit = 000000Test = DUI Dr. Lic. = ND/TESTLot No = 05620080A1Cyl No = 45Expiration Date = 04/05/2022County = 02Oper No. = 100012

Operator Signature DANA RUSTEBAKKE

Remarks:

RFI Test

Form 106-18000