

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number SCO 6691	Instrument Location Kenmove Police	ce Dest.		
Reason for Install/Repair				
🕍 Install After Receiving From Crime Laboratory 📗 Install After Location Change				
Other (Specify)				
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
☑ 2. Telephone Line Connected to intoxilyzer® 8000.				
☑ 3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
☑ 6. Run Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
7. Repair and/or Maintenance Performed (if any):				
8, Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.				
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Int 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form		
Field Inspector Signature		Date 7/29/2021		
Crime Laboratory Use Only				
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on f County of Burleigh, North Dakota, is certified to be a true and correct co	ile at the Office of Attorney General, Crin			

FORM 104-G

Reviewed/Certified By

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006691
Location = KMPD 8164.14.00 09/16
07/29/2021 17:24

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 02621080A1 Standard Cyl #? 30 Standard Expiration? 03/05/2023 Oper No? 132664

Flow Cal. Date: 08/28/2018 Slope 699 Intercept -660610

IR Calibration Date: 08/28/2018

3um 9um

Oth Coef(*100): -9426 -20764

1st Coef(*100): 258193 134855 2nd Coef(*100): 3532 1465 H2O adj(mg/l*10k): 408 547

****** Printer Test End ********

Operator Signature CHRISTOPHER ALMLIE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = KMPD
07/29/2021

Alcohol Analyzer SN 80-006691 8164.14.00 09/16 17:24

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:25
02 Std. Gas	0.080	17:25
03 Room Air	0.000	17:26
04 Std. Gas	0.081	17:26
05 Room Air	0.000	17:27
06 Std. Gas	0.080	17:27
07 Room Air	0.000	17:27

Lot No = 02621080A1

Cyl No = 30

Exp Date = 03/05/2023

County = 51

Oper No. = 132664

Operator Signature CHRISTOPHER ALMLIE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006691
Location = KMPD 8164.14.00 09/16
07/29/2021 17:30

Test	AC	Time
01 Room Air	RFI*	17:32
02 Room Air	0.000	17:33

*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST TEST

Sub DOB = 03/29/1991

Sub Sex = Male Weight = NA Test = OTH Cit = NA

Dr. Lic. = ND/NA Lot No = 02621080A1

Cyl No = 30

Expiration Date = 03/05/2023

County = 51

Oper No. = 132664

Operator Signature CHRISTOPHER ALMLIE

Remarks:

Form 106-I8000