

INTOXILYZER® 8000 REPAIR AND MAINTENANCE

Intoxilyzer® 8000 Serial Number: 80-006678 Location: TOXL

1. Reported Issue:

SCANNER NOT RECOGNIZING OPERATOR CARD OR LICENSES
CAN MANUALLY ENTER "INVALID PASSWORD" MESSAGE,
BAD KEYBOARD

2. Found/Determined:

INVALID PASSWORD → EXPIRED OPERATOR,
BAD KEYBOARD CONNECTION - INTERMITTENT FUNCTION
NO ISSUES FOUND WITH 2-D SCANNER
INSTRUMENT ONLY RAN ON BATTERY POWER → BAD FUSE

3. Repair/Maintenance Performed:

REPLACED FUSE → AC POWER CORD → CHARGING BATTERY NOW
REPLACED KEYBOARD - FUNCTIONS PROPERLY
2-D SCANNER HAS NOT FAILED TO READ A BARCODE PROPERLY.

4. Check and reset if needed; time, date, and location. (Level 2, E)

5. Tank monitor display and regulator within 50 psi. (Level 3, D, G)

Display: 215 psi Regulator: 250 psi

If needed, completed tare of tank sensor. (Level 3, M, C, G)

6. Does the flow need to be calibrated?

If Yes, attach paperwork.

Yes or No

7. Does the optical bench need a calibration adjustment?

If Yes, attach test records.

Yes or No

8. Does the instrument require further testing?

If No, sign and date at end of form.

If Yes, state testing required or: Perform all tests for repair.

Yes or No

ALL TESTS WILL BE PERFORMED DUE TO RETURNING
TO THE LAB FROM THE FIELD.

Configure simulator(s) for the following tests: (Level 1, S) Attach test records

9. Wet Calibration Check: (ACA Test)

Low AC (Level 1, C): Known value \leq 0.03 AC. 0.030 AC
Sim SN: MP3059 Lot #: 202201B Exp. Date: 1.11.24

High AC (Level 1, C): Known value \geq 0.25 AC. 0.400 AC
Sim SN: MP5320 Lot #: 202103E Exp. Date: 3.24.23

Wet Calibration Check results are within \pm 0.005 or \pm 5% (whichever is greater) of stated AC value.

Configure dry gas standard for the following tests: (Level 1, S) Attach test records.

10. Print test. (Level 1, P)

11. Interferent Check. 0.10 AC ethanol plus 0.05% Acetone. (Level 1, B)
Sim SN: DR3847 Lot #: FCS 6 Display reads "Interferent Detect"

12. RFI Check. (CMS mode) Display reads "RFI Detect"

13. Dry Calibration Check/Stability Test (Level 3, T, S, D, 15)

Known: 0.080 AC Ethanol Gas Standard.

Lot No. 02621080A1 Cylinder No. 4 Expiration Date 3.5.23

Average AC: 0.0816 Std. Dev. 0.0005 Rel. Std. Dev. 0.63%

Remarks/Notes:

2-D SCANNER READ ALL BARCODES DURING TESTING.
PERFORMED OPERATOR.DAT FILE UPLOAD.

Instrument is acceptable to be used in the field.

Charles E. Ed
Breath Analyst Signature

3/29/22
Date

NA
Reviewed by

NA
Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006678
Location = TOXL 8164.14.00 09/16
03/29/2022 12:33

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:34
02 Std. Sol.	0.029	12:34
03 Room Air	0.000	12:35
04 Std. Sol.	0.029	12:36
05 Room Air	0.000	12:36
06 Std. Sol.	0.029	12:37
07 Room Air	0.000	12:37

08 Sim Temp = 34.0°C

Simul Ser No = MP3059
Std Sol No = 202201B
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

LOW AC
0.030 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006678
Location = TOXL 8164.14.00 09/16
03/29/2022 12:43

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:44
02 Std. Sol.	0.404	12:45
03 Room Air	0.000	12:45
04 Std. Sol.	0.403	12:46
05 Room Air	0.000	12:46
06 Std. Sol.	0.402	12:47
07 Room Air	0.000	12:48

08 Sim Temp = 34.0°C

Simul Ser No = MP5320
Std Sol No = 202103E
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

HIGH AC
0.400 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006678
Location = TOXL 8164.14.00 09/16
03/29/2022 12:49

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=-|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=-|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 02621080A1
Standard Cyl #? 4
Standard Expiration? 03/05/2023
Oper No? 666666

Flow Cal. Date: 08/22/2016
Slope 686
Intercept -611596

IR Calibration Date: 08/22/2016
 3um 9um

0th Coef(*100): -15445 -21650
1st Coef(*100): 251854 131766
2nd Coef(*100): 3094 1378
H2O adj(mg/l*10k): 549 513

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks: PRINT TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006678
Location = TOXL 8164.14.00 09/16
03/29/2022 12:50

Test	AC	Time
01 Room Air	0.000	12:51
02 *Subject Test	INT*	12:51
03 Room Air	0.000	12:52

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS921456
Lot No = 02621080A1
Cyl No = 4
Expiration Date = 03/05/2023
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Interferent Check*
0.10 AC + 0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006678
Location = TOXL 8164.14.00 09/16
03/29/2022 12:53

Test	AC	Time
01 Diagnostic	OK	12:54
02 Room Air	0.000	12:54
03 *Subject Test	RFI*	12:54
04 Room Air	0.000	12:55

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/DIS921456
Lot No = 02621080A1
Cyl No = 4
Expiration Date = 03/05/2023
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

RFI CHECK

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006678
Location = TOXL 8164.14.00 09/16
03/29/2022 12:55

Test	AC	Time
01 Room Air	0.000	12:56
02 Std. Gas	0.082	12:56
03 Room Air	0.000	12:57
04 Std. Gas	0.082	12:57
05 Room Air	0.000	12:58
06 Std. Gas	0.081	12:58
07 Room Air	0.000	12:59
08 Std. Gas	0.081	12:59
09 Room Air	0.000	12:59
10 Std. Gas	0.081	13:00
11 Room Air	0.000	13:00
12 Std. Gas	0.082	13:01
13 Room Air	0.000	13:01
14 Std. Gas	0.082	13:02
15 Room Air	0.000	13:02
16 Std. Gas	0.082	13:02
17 Room Air	0.000	13:03
18 Std. Gas	0.082	13:03
19 Room Air	0.000	13:04
20 Std. Gas	0.082	13:04
21 Room Air	0.000	13:05
22 Std. Gas	0.082	13:05
23 Room Air	0.000	13:06
24 Std. Gas	0.081	13:06
25 Room Air	0.000	13:06
26 Std. Gas	0.082	13:07
27 Room Air	0.000	13:07
28 Std. Gas	0.081	13:08
29 Room Air	0.000	13:08
30 Std. Gas	0.081	13:09
31 Room Air	0.000	13:09
32 Cal Check Stats:		
33 Average	0.0816	
34 Std Dev	0.0005	
35 Rel Std Dev	0.63%	



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK / STABILITY TEST
8.080 AC

Form 106-I8000