

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-006682

Location: TUCL

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 248 psi Regulator: 280 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check

- a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.020 AC  
Sim. Ser #: MP6036 Lot #: 21340 Exp. Date: 8/24/23
- b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.300 AC  
Sim. Ser #: MP3068 Lot #: 202201F Exp. Date: 1/18/24

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR3847 Lot #: 1CS6 Exp. Date: NA  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 0262108A1 Cyl #: 2 Exp. Date: 3/5/2023

Test 1 <u>0.081</u>	Test 4 <u>0.080</u>	Test 7 <u>0.081</u>
Test 2 <u>0.080</u>	Test 5 <u>0.080</u>	Test 8 <u>0.080</u>
Test 3 <u>0.080</u>	Test 6 <u>0.081</u>	Test 9 <u>0.081</u>
Average <u>0.080</u>		

C. Remarks/Maintenance: Annual inspection.

Replace breath tube o-ring + simulator return o-ring.

Instrument is acceptable to be used in the field.

Janella Potrichillo  
Breath Analyst Signature

6/14/22  
Date

[Signature]  
Reviewed by

06/15/22  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000           SN 80-006682  
Location = TOXL                    8164.14.00 09/16  
06/14/2022                           20:04

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?         No  
Display Volume?             No  
Disable On Memfull?         Yes  
# of Print Copies?          1  
Select Std (D/W/I)?         Wet  
Standard Value?             0.020  
Standard Solution #?         21340  
Sim Serial #?                MP6036  
Oper No?                     132128

Flow Cal. Date:             09/06/2016  
Slope                        701  
Intercept                    -706193

IR Calibration Date:        09/06/2016  
                              3um                    9um

-----  
0th Coef(\*100):            -17150                -23247  
1st Coef(\*100):            271726                134617  
2nd Coef(\*100):            831                    1215  
H2O adj(mg/l\*10k):         393                    518

\*\*\*\*\* Printer Test End \*\*\*\*\*

Operator Signature  
JANELLE PORTSCHELLER

Remarks: *Print test*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006682  
Location = TOXL      8164.14.00 09/16  
06/14/2022      20:06

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:07
02 Std. Sol.	0.019	20:07
03 Room Air	0.000	20:08
04 Std. Sol.	0.019	20:09
05 Room Air	0.000	20:09
06 Std. Sol.	0.019	20:10
07 Room Air	0.000	20:10

08 Sim Temp = 34.0°C

Simul Ser No = MP6036

Std Sol No = 21340

County = 08

Oper No. = 132128

*Janelle Portscheller*

Operator Signature  
JANELLE PORTSCHELLER

Remarks: *Low AC Check 0.020*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006682  
Location = TOXL      8164.14.00 09/16  
06/14/2022      20:16

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:16
02 Std. Sol.	0.291	20:18
03 Room Air	0.000	20:19
04 Std. Sol.	0.293	20:19
05 Room Air	0.000	20:20
06 Std. Sol.	0.294	20:21
07 Room Air	0.000	20:21

08 Sim Temp = 34.0°C

Simul Ser No = MP3068  
Std Sol No = 202201F  
County = 08      Oper No. = 132128

  
\_\_\_\_\_  
Operator Signature  
JANELLE PORTSCHELLER

Remarks: *High Calibration Check 0.300*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006682  
Location = TOXL      8164.14.00 09/16  
06/14/2022      20:24

Test	AC	Time
01 Room Air	0.000	20:25
02 *Subject Test	INT*	20:28
03 Room Air	0.000	20:28

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT CK  
Dr. Lic. = ND/DIS921456  
Lot No = 02621080A1  
Cyl No = 2  
Expiration Date = 03/05/2023  
County = 08      Oper No. = 132128



Operator Signature  
JANELLE PORTSCHELLER

Remarks: ~~RFI check~~ <sup>JP</sup> Interferent check

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006682  
Location = TOXL                    8164.14.00 09/16  
06/14/2022                            20:35

Test	AC	Time
01 Diagnostic	OK	20:36
02 Room Air	0.000	20:36
03 *Subject Test	RFI*	20:37
04 Room Air	0.000	20:37

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHK  
Dr. Lic. = ND/DIS921456  
Lot No = 02621080A1  
Cyl No = 2  
Expiration Date = 03/05/2023  
County = 08                            Oper No. = 132128

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
JANELLE PORTSCHELLER

Remarks: *RFI check*





Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006682  
Location = TOXL                    8164.14.00 09/16  
06/14/2022                            20:45

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:45
02 Std. Gas	0.080	20:45
03 Room Air	0.000	20:46
04 Std. Gas	0.080	20:46
05 Room Air	0.000	20:47
06 Std. Gas	0.081	20:47
07 Room Air	0.000	20:48

Lot No = 02621080A1  
Cyl No = 2  
Exp Date = 03/05/2023  
County = 08                            Oper No. = 132128

  
\_\_\_\_\_  
Operator Signature  
JANELLE PORTSCHELLER

Remarks: *Dry Calibration Check # 4-6*




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006682  
Location = TOXL      8164.14.00 09/16  
06/14/2022      20:48

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:49
02 Std. Gas	0.081	20:49
03 Room Air	0.000	20:49
04 Std. Gas	0.080	20:50
05 Room Air	0.000	20:50
06 Std. Gas	0.081	20:51
07 Room Air	0.000	20:51

Lot No = 02621080A1  
Cyl No = 2  
Exp Date = 03/05/2023  
County = 08      Oper No. = 132128

  
\_\_\_\_\_  
Operator Signature  
JANELLE PORTSCHELLER

Remarks: *Dry Calibration check #7-9*