

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-006671 Location: ToxL

A. General Setup and Checks:

1. Diagnostic Tests Pass and Instrument in Ready Mode
2. Breath tube heated
3. Date, time and location code (Level 2,E). Re-set if necessary.
4. Print test (Level 1,P). Sign and attach test record.
5. Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
 - a. Display: 411 psi Regulator: 450 psi

B. Tests (Sign and attach test records):

- Configure simulator for the following tests (Level 1,S).
1. Wet Calibration Check
 - a. Low AC (Level 1,C): Known Value \leq 0.03 AC: 0.020 AC
Sim. Ser #: MP3070 Lot #: 22CP41184 Exp. Date: 3/1/24
 - b. High AC (Level 1,C): Known Value \geq 0.25 AC: 0.400 AC
Sim. Ser #: MP3065 Lot #: 202103E Exp. Date: 3/24/23

Configure dry gas standard for the following tests (Level 1,S).

2. Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone
Sim. Ser #: DR3847 Lot #: 1656 Exp. Date: NA
 Display reads "Interferent Detect"
3. RFI Check (CMS Mode) Display reads "RFI Detect"
4. Dry Calibration Check (Level 1,C): Known Value 0.080 AC
Gas Cyl Lot #: 02621080A1 Cyl #: 16 Exp. Date: 3/5/23

| | | |
|----------------------|---------------------|---------------------|
| Test 1 <u>0.080</u> | Test 4 <u>0.079</u> | Test 7 <u>0.080</u> |
| Test 2 <u>0.080</u> | Test 5 <u>0.079</u> | Test 8 <u>0.080</u> |
| Test 3 <u>0.081</u> | Test 6 <u>0.080</u> | Test 9 <u>0.080</u> |
| Average <u>0.080</u> | | |

C. Remarks/Maintenance: Annual inspection
Replaced breath tube + simulator return o-rings.

Instrument is acceptable to be used in the field.

Janelle Putscheller
Breath Analyst Signature

6/22/22
Date

[Signature]
Reviewed by

062422
Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = ~~DEVL~~ *JP* 8164.14.00 09/16
06/22/2022 *TOXL* 16:03

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 16:03 |
| 02 Std. Sol. | 0.018 | 16:05 |
| 03 Room Air | 0.000 | 16:06 |
| 04 Std. Sol. | 0.019 | 16:07 |
| 05 Room Air | 0.000 | 16:07 |
| 06 Std. Sol. | 0.019 | 16:08 |
| 07 Room Air | 0.000 | 16:08 |

08 Sim Temp = 34.0°C

Simul Ser No = MP3070
Std Sol No = 22CP41184
County = 08 Oper No. = 132128

Janelle Portscheller

Operator Signature
JANELLE PORTSCHELLER

Remarks: *Low Ac check 0.020*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = ~~DEVL~~ ^{JLP} 8164.14.00 09/16
06/22/2022 ~~TRVL~~ 16:18

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 16:19 |
| 02 Std. Sol. | 0.395 | 16:19 |
| 03 Room Air | 0.000 | 16:20 |
| 04 Std. Sol. | 0.396 | 16:21 |
| 05 Room Air | 0.000 | 16:21 |
| 06 Std. Sol. | 0.397 | 16:22 |
| 07 Room Air | 0.000 | 16:23 |

08 Sim Temp = 34.0°C

Simul Ser No = MP3065

Std Sol No = 202103E

County = 08

Oper No. = 132128



Operator Signature
JANELLE PORTSCHELLER

Remarks: High Ac Check 0.300

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = ~~DEVI~~ 8164.14.00 09/16
06/22/2022 *Taxi* 16:27

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 16:28 |
| 02 *Subject Test | INT* | 16:29 |
| 03 Room Air | 0.000 | 16:29 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT CK
Dr. Lic. = ND/DIS821456
Lot No = 02621080A1
Cyl No = 16
Expiration Date = 03/05/2023
County = 08 Oper No. = 132128

Janelle Portscheller

Operator Signature
JANELLE PORTSCHELLER

Remarks: *Interferent Check*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = ~~DEVL~~ 8164.14.00 09/16
06/22/2022 *TuxL Jar* 16:33

| Test | AC | Time |
|------------------|-------|-------|
| 01 Diagnostic | OK | 16:34 |
| 02 Room Air | 0.000 | 16:35 |
| 03 *Subject Test | RFI* | 16:35 |
| 04 Room Air | 0.000 | 16:36 |

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHK
Dr. Lic. = ND/DIS821456
Lot No = 02621080A1
Cyl No = 16
Expiration Date = 03/05/2023
County = 08 Oper No. = 132128

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Janelle Portscheller

Operator Signature
JANELLE PORTSCHELLER

Remarks: *RFI Check*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = ~~DEVL~~ *TOXL dev* 8164.14.00 09/16
06/22/2022 16:48

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 16:49 |
| 02 Std. Gas | 0.079 | 16:50 |
| 03 Room Air | 0.000 | 16:50 |
| 04 Std. Gas | 0.079 | 16:50 |
| 05 Room Air | 0.000 | 16:51 |
| 06 Std. Gas | 0.080 | 16:51 |
| 07 Room Air | 0.000 | 16:52 |

Lot No = 02621080A1
Cyl No = 16
Exp Date = 03/05/2023
County = 08 Oper No. = 132128

Janelle Portscheller

Operator Signature
JANELLE PORTSCHELLER

Remarks: *Dry Calibration Check #4-6*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = ~~DEVL~~ *TOXL* 8164.14.00 09/16
06/22/2022 16:56

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 16:57 |
| 02 Std. Gas | 0.080 | 16:57 |
| 03 Room Air | 0.000 | 16:57 |
| 04 Std. Gas | 0.080 | 16:58 |
| 05 Room Air | 0.000 | 16:58 |
| 06 Std. Gas | 0.080 | 16:59 |
| 07 Room Air | 0.000 | 16:59 |

Lot No = 02621080A1
Cyl No = 16
Exp Date = 03/05/2023
County = 08 Oper No. = 132128

Janelle Portscheller

Operator Signature
JANELLE PORTSCHELLER

Remarks: *Dry Calibration Check #7-9*