



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <i>80-006674</i>	Instrument Location <i>WASH</i>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

- Check When Done:
- 1. Surge Protector Installed/Property Grounded.
  - 2. Telephone Line Connected to Intoxilyzer® 8000.
  - 3. Breath Tube Heated.
  - 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
  - 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
  - 6. Run Tests:
    - A. Print Test (Level 1, Function P).
    - B. ACA Test (Level 1, Function C).
    - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
  - 7. Repair and/or Maintenance Performed (if any): \_\_\_\_\_
  - 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
  - 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
  - 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>Richard A. Johnson</i>	Date <i>6/23/22</i>
--	------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <i>Janelle Pertschiller</i>	Certified Date <i>06/28/2022</i>
--	-------------------------------------

*Note: The location code was not changed back to WASH after the annual inspection. On 6/27/22 I spoke with Richard Johnson who changed the location code back to WASH and verified no subject tests had been run on the instrument. 06/28/22 JPP*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006674  
Location = TOXL                    8164.14.00 09/16  
06/23/2022                    *WASH 06/28/22*                    13:53

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:                    Enabled  
Start Test Sequence:                DABACABA  
Display Prelim Rslt?                Yes  
Display Third Digit?                Yes  
Inhib Printer(Y/N)?                No  
Display Volume?                    No  
Disable On Memfull?                Yes  
# of Print Copies?                1  
Select Std (D/W/I)?                Dry  
Standard Value?                    0.080  
Standard Lot #?                    26021080A1  
Standard Cyl #?                    31  
Standard Expiration?                10/05/2023  
Oper No?                            118707

Flow Cal. Date:                    06/28/2017  
Slope                                707  
Intercept                            -708478

IR Calibration Date:                06/28/2017  
    3um                                9um

-----  
0th Coef(\*100):                    -3915                                -18046  
1st Coef(\*100):                    256520                                127506  
2nd Coef(\*100):                    3843                                 1644  
H2O adj(mg/l\*10k):                428                                 568

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
\_\_\_\_\_  
Operator Signature  
RICHARD JOHNSON

Remarks:


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-006674  
Location = TOXL 8164.14.00 09/16  
06/23/2022 *WASH 06/28/22* 13:49  
*JA*

DRY CAL CHECK
---------------

Test	AC	Time
01 Room Air	0.000	13:50
02 Std. Gas	0.079	13:50
03 Room Air	0.000	13:51
04 Std. Gas	0.079	13:51
05 Room Air	0.000	13:52
06 Std. Gas	0.079	13:52
07 Room Air	0.000	13:53

Lot No = 26021080A1  
Cyl No = 31  
Exp Date = 10/05/2023  
County = 28 Oper No. = 118707

  
\_\_\_\_\_  
Operator Signature  
RICHARD JOHNSON

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006674  
Location = TOXL      8164.14.00 09/16  
06/23/2022      *WASH 06/28/22*      13:54

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	13:55
02 Room Air	0.000	13:55

\*Invalid Test  
Inhibited - RFI

Lot No = 26021080A1  
Cyl No = 31  
Exp Date = 10/05/2023  
County = 28      Oper No. = 118707

*Richard Johnson*  
\_\_\_\_\_  
Operator Signature  
RICHARD JOHNSON

Remarks: