

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-006690

Location: TOXL

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 516 psi Regulator: 500 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check
  - a.  Low AC (Level 1,C): Known Value  $\leq$  0.03 AC: 0.020 AC  
Sim. Ser #: MP5290 Lot #: 21340 Exp. Date: 8/24/23
  - b.  High AC (Level 1,C): Known Value  $\geq$  0.25 AC: 0.300 AC  
Sim. Ser #: MP3067 Lot #: 21420 Exp. Date: 10/20/23

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR 7345 Lot #: 1CS 7 Exp. Date: NA  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 02621080A1 Cyl #: 31 Exp. Date: 3/5/23

Test 1 <u>0.080</u>	Test 4 <u>0.080</u>	Test 7 <u>0.079</u>
Test 2 <u>0.079</u>	Test 5 <u>0.080</u>	Test 8 <u>0.080</u>
Test 3 <u>0.080</u>	Test 6 <u>0.079</u>	Test 9 <u>0.080</u>
Average <u>0.080</u>		

C. Remarks/Maintenance: Annual inspection

Note: Instrument had an Intoximeters gas cylinder installed  
when the instrument arrived at the lab.  
Replaced breath tube o-ring.

Instrument is acceptable to be used in the field.

Janelle Putschiller  
Breath Analyst Signature

06/28/22  
Date

[Signature]  
Reviewed by

062922  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer          Alcohol Analyzer  
North Dakota Model 8000          SN 80-006690  
Location = TOXL                  8164.14.00 09/16  
06/28/2022                                  13:46

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopQRSTUVWXYZ!@#%&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopQRSTUVWXYZ!@#%&\*()\_+?

Current Instrument Setup


Data Entry Mode:                  Enabled  
Start Test Sequence:              DABACABA  
Display Prelim Rslt?              Yes  
Display Third Digit?              Yes  
Inhib Printer(Y/N)?              No  
Display Volume?                  No  
Disable On Memfull?              Yes  
# of Print Copies?                1  
Select Std (D/W/I)?              Dry  
Standard Value?                  0.082  
Standard Lot #?                  AG030202  
Standard Cyl #?                  50  
Standard Expiration?              28/10/2022  
Oper No?                          132128

Flow Cal. Date:                  08/23/2018  
Slope                              677  
Intercept                        -593298

IR Calibration Date:              08/23/2018  
    3um                              9um

	3um	9um
0th Coef(*100):	-8380	-21663
1st Coef(*100):	265142	130429
2nd Coef(*100):	3933	1429
H2O adj(mg/l*10k):	412	539

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
JANELLE PORTSCHELLER

Remarks: *Print test*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006690  
Location = TOXL      8164.14.00 09/16  
06/28/2022      13:52

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:52
02 Std. Sol.	0.018	13:54
03 Room Air	0.000	13:54
04 Std. Sol.	0.018	13:55
05 Room Air	0.000	13:56
06 Std. Sol.	0.018	13:56
07 Room Air	0.000	13:57

08 Sim Temp = 34.0°C

Simul Ser No = MP5290

Std Sol No = 21340

County = 08

Oper No. = 132128



Operator Signature  
JANELLE PORTSCHELLER

Remarks:

*Low AC check 0.020*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006690  
Location = TOXL                    8164.14.00 09/16  
06/28/2022                            14:11

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:12
02 Std. Sol.	0.291	14:13
03 Room Air	0.000	14:13
04 Std. Sol.	0.293	14:14
05 Room Air	0.000	14:15
06 Std. Sol.	0.292	14:15
07 Room Air	0.000	14:16


08 Sim Temp = 34.0°C

Simul Ser No = MP3067

Std Sol No = 21420

County = 08

Oper No. = 132128



Operator Signature  
JANELLE PORTSCHELLER

Remarks:

*High AC check 0.300*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-006690  
Location = TOXL 8164.14.00 09/16  
06/28/2022 14:18

Test	AC	Time
01 Room Air	0.000	14:19
02 *Subject Test	INT*	14:20
03 Room Air	0.000	14:21

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male Weight = NA  
Test = DUI Cit = INTERFERENT CK  
Dr. Lic. = ND/TES989643  
Lot No = 02621080A1  
Cyl No = 31  
Expiration Date = 03/05/2023  
County = 08 Oper No. = 132128

  
Operator Signature  
JANELLE PORTSCHELLER

Remarks: *Interferent check*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-006690  
Location = TOXL                              8164.14.00 09/16  
06/28/2022    14:21

Test	AC	Time
01 Diagnostic	OK	14:22
02 Room Air	0.000	14:22
03 *Subject Test	RFI*	14:23
04 Room Air	0.000	14:24

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                              Weight = NA  
Test = DUI                                      Cit = RFI CHK  
Dr. Lic. = ND/TES989643  
Lot No = 02621080A1  
Cyl No = 31  
Expiration Date = 03/05/2023  
County = 08                              Oper No. = 132128

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
JANELLE PORTSCHELLER

Remarks: *RFI check*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-006690  
Location = TOXL                              8164.14.00 09/16  
06/28/2022    14:24

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:25
02 Std. Gas	0.080	14:25
03 Room Air	0.000	14:25
04 Std. Gas	0.079	14:26
05 Room Air	0.000	14:26
06 Std. Gas	0.080	14:27
07 Room Air	0.000	14:27

Lot No = 02621080A1  
Cyl No = 31  
Exp Date = 03/05/2023  
County = 08                                      Oper No. = 132128

  
Operator Signature  
JANELLE PORTSCHELLER

Remarks: *Dry Calibration Check #1-3*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006690  
Location = TOXL      8164.14.00 09/16  
06/28/2022      14:28

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:29
02 Std. Gas	0.080	14:29
03 Room Air	0.000	14:29
04 Std. Gas	0.080	14:30
05 Room Air	0.000	14:30
06 Std. Gas	0.079	14:31
07 Room Air	0.000	14:31

Lot No = 02621080A1  
Cyl No = 31  
Exp Date = 03/05/2023  
County = 08      Oper No. = 132128

  
\_\_\_\_\_  
Operator Signature  
JANELLE PORTSCHELLER

Remarks: *Dry Calibration Check #4-6*



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006690  
Location = TOXL      8164.14.00 09/16  
06/28/2022      14:32

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:33
02 Std. Gas	0.079	14:33
03 Room Air	0.000	14:34
04 Std. Gas	0.080	14:34
05 Room Air	0.000	14:34
06 Std. Gas	0.080	14:35
07 Room Air	0.000	14:35

Lot No = 02621080A1  
Cyl No = 31  
Exp Date = 03/05/2023  
County = 08      Oper No. = 132128

  
\_\_\_\_\_  
Operator Signature  
JANELLE PORTSCHELLER

Remarks: *Dry Calibration Check #7-9*