

INTOXILYZER 5000 G and 8000 INSPECTION

Serial No.: 80-002667

Location: TOXL

Check When Done

- A. Is the warm-up time less than 20 minutes? ✓
- B. Is a three-pronged grounded outlet used? ✓
- C. Is the breath tube heated? ✓
- D. Is the diagnostic check complete? ✓
- E. Is the time, date, and year correct (re-set if necessary)? ✓
- F. Print test. (Attach test record.) ✓
- G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN DR5144 Lot No. Rep AC 0.019 ✓
- H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN DR5187 Lot No. Rep AC 0.285 ✓
- I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN DR5114 Lot No. Rep AC 0.100 ✓
- J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) ✓
- K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.)

Lot No. <u>659358</u>	Cyl No. <u>19</u>	Exp. Date <u>03 Sept 2011</u>
Test 1 <u>0.079</u>	Test 4 <u>0.080</u>	Test 7 <u>0.080</u>
Test 2 <u>0.080</u>	Test 5 <u>0.080</u>	Test 8 <u>0.080</u>
Test 3 <u>0.080</u>	Test 6 <u>0.080</u>	Test 9 <u>0.080</u>
Average <u>0.080</u>		

L. Remarks/Maintenance record: Software upgrade 8/16/13. Annual inspection

<u>14 Jan 2010</u> Date	<u>Deb Kashner</u> Field Inspector's Signature
<u>14 Jan 2010</u> Date	<u>Mary Pearson</u> Reviewed by

80-002667
53/61 KH

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.13.00 06/09
01/14/2010 14:02

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:03
02 Std. Sol.	0.019	14:04
03 Room Air	0.000	14:04
04 Std. Sol.	0.020	14:05
05 Room Air	0.000	14:05
06 Std. Sol.	0.019	14:06
07 Room Air	0.000	14:07

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Low AC 0.019 AC

Form 106-18000

80-002667
5/4/11 FR

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.13.00 06/09
01/14/2010 11:05

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659358
Standard Cyl #? 19
Standard Expiration? 09/03/2011
Oper No? 777777

Flow Cal. Date: 09/10/2007
Slope 693
Intercept -570251

IR Calibration Date: 09/10/2007
3um 9um

0th Coef(*100): -21400 -13991
1st Coef(*100): 268177 135590
2nd Coef(*100): 1202 658
H2O adj(mg/l*10k): 369 165

***** Printer Test End *****

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Print Test

Form 106-I8000

80-002667
55/61 PM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.13.00 06/09
01/14/2010 14:41

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:42
02 Std. Sol.	0.291	14:43
03 Room Air	0.000	14:43
04 Std. Sol.	0.292	14:44
05 Room Air	0.000	14:44
06 Std. Sol.	0.293	14:45
07 Room Air	0.000	14:46

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Linearity test 0.285 AC

Form 106-I8000

*80-002667
Siga FM*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.13.00 06/09
01/14/2010 15:07

Test	AC	Time
01 Room Air	0.000	15:08
02 *Subject Test	INT*	15:08
03 Room Air	0.000	15:09

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 19
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Interferent test

Form 106-I8000

80-002667
P. 58/61 PM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.13.00 06/09
01/14/2010 15:10

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:11
02 Std. Gas	0.079	15:11
03 Room Air	0.000	15:12
04 Std. Gas	0.080	15:12
05 Room Air	0.000	15:12
06 Std. Gas	0.080	15:13
07 Room Air	0.000	15:13

Lot No = 659358
Cyl No = 19
Exp Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks: Calibration check 0.080 AC

Form 106-I8000

80-002467
59/01 JTW

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.13.00 06/09
01/14/2010 16:03

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:03
02 Std. Gas	0.080	16:04
03 Room Air	0.000	16:04
04 Std. Gas	0.080	16:04
05 Room Air	0.000	16:05
06 Std. Gas	0.080	16:05
07 Room Air	0.000	16:06

Lot No = 659358
Cyl No = 19
Exp Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kashur
Operator Signature
DEB KASHUR

Remarks:
Calibration check 0.080 AC

Form 106-I8000

80-002667
6/16/11 FM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.13.00 06/09
01/14/2010 16:07

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:07
02 Std. Gas	0.080	16:08
03 Room Air	0.000	16:08
04 Std. Gas	0.080	16:08
05 Room Air	0.000	16:09
06 Std. Gas	0.080	16:09
07 Room Air	0.000	16:10

Lot No = 659358
Cyl No = 19
Exp Date = 09/03/2011
County = 08 Oper No. = 777777



Operator Signature
DEB KASHUR

Remarks:

Calibration check

Form 106-I8000

*80-002667
6/16/10*