

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-006489

1) GENERAL Setup and Checks:

- a) Power On
 - 5 min. countdown during initial warm-up
 - 1 min. countdown to begin test after standby mode
- b) Display Intensity (Ctrl + ↑ or ↓)
- c) Password Security:
 - Level 1 (Operator) Level 2 (Field Inspector) Level 3 (Lab Personnel)
- d) Printer Setup (Level 3,S,P)
- e) Forms Load (CMI, Inc., CoBrA Software)

2) LEVEL 1 Setup and Checks:

- a) Internal Printer Test and Paper Feed (Level 1,P). Attach record. ✓
- b) External Printer Test (Level 1,P). Attach record. ✓
- c) Barcode Reader (Level 1,B):
 - Operator ID ND License R key: Refusal. Attach record. ✓

3) LEVEL 2 Setup and Checks:

- a) Continuous Air Blank (Level 2,A)
- b) Diagnostic (Level 2,D). Attach record. ✓
- c) Prelim Data Entry (Level 2,E)
- d) Ser Number/Version Info (Level 2,V):
 - Inst Ser#: 80-006489 DSP Ver: 1362.14
 - Core Ver: 8164.13 CPLD Ver: 7

4) LEVEL 3 Setup and Checks:

- a) DVM Monitor (Level 3,D,D): 3: 12792 s: 1 C: 47.06
9: 13160 s: 1 B: 45 0
(RFI indicator with keyed radio): 1
- b) Function (F) Keys : F5 – Tone F9 - Fan
 F6 – Green LED F10 – Pump
 F7 – Red LED F11 – Solenoid
- c) Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 939 CF: 1.07
- d) Flow Monitor (Level 3,D,F): Flow Sensor, R: 235 F: 0.000 L/s

- e) Temperature Monitor (Level 3,D,T): C: 47.0 A: 48.5 I: 0.0
B: 45.1 V: 0.0
- f) Tank Monitor (Level 3,D,G): P (display): 856 psi P (gauge): -850 psi
- g) General Setup (Level 3,S,G)
- h) Configure Start Test (Level 3,S,T)
- i) Network (Level 3,S,N)

5) **TESTS (attach all test records):**

Configure wet bath simulator standards before running tests (Level 1,S):

- a) Zero AC (Level 1,C):
Sim ser #: MP3067 Lot #: NA ^{CEC, MILLI} AZO Exp date: NA
Known value: 0.000
Results: 0.000 AC 0.000 AC 0.000 AC
- b) Low AC (Level 1,C):
Sim ser #: MP3066 Lot #: 15160 Exp date: 7/4/17
Known value: 0.020
Results: 0.019 AC 0.019 AC 0.019 AC
- c) High AC (Level 1,C):
Sim ser #: MP3068 Lot #: 14102 Exp date: 4/23/16
Known value: 0.400
Results: 0.406 AC 0.403 AC 0.402 AC

Configure dry gas standard before running tests (Level 1,S):

- d) Ethanol + 0.05% Acetone (Level 1,B):
Sim ser #: DR5135 Lot #: ICS 4 Exp date: 9/9/2017
Known value: 0.100 AC + 0.05% ACETONE
 Display reads "Interferent Detect"
- e) Acetone (Level 1,B):
Sim ser #: DR5143 Lot #: ACE 2 Exp date: 9/29/2017
Known value: 0.05% ACETONE
 Display reads "Interferent Detect"
- f) Methanol (Level 1,B):
Sim ser #: DR3377 Lot #: MET 2 Exp date: 9/29/2017
Known value: 0.10% METHANOL
 Display reads "Interferent Detect"
- g) Isopropanol (Level 1,B):
Sim ser #: DR3453 Lot #: IPA 2 Exp date: 9/29/17
Known value: 0.10% ISOPROPRANOL
 Display reads "Interferent Detect"

- h) Ambient Fail (Level 1,B): Display and test record read correct.
- i) Invalid Sample (Level 1,B): Display and test record read correct.
- j) Sequence Aborted (Level 1,B): Display and test record read correct.
- k) ABA (Level 1,B): Test record reads correct.
- l) Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
- m) RFI Detected (B,C, or CMS mode): Display and test record read correct.
- n) Tank Pressure Below Minimum (CMS mode): Display reads correct.
- o) Deficient Samples (CMS mode): Display and test record read correct.
- p) Improper Sample (CMS mode): Display and test record read correct.
- q) Difference Too Great (CMS mode): Display and test record read correct.
- r) CMS (CMS mode): Test record reads correct.
- s) Stability Test (Level 3,T,S):

Ave: 0.0806 Std Dev: 0.0007 Rel Std Dev: 0.85%

Maintenance/Repairs/Notes: _____

Instrument is acceptable to be used in the field.

Charles E. Eder
Field Inspector's Signature

15 JAN 16
Date

Deb Shanaver
Reviewed by

15-Jan-16
Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab, Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 14:53

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

32 characters per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

42 characters per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU


Current Test Results
Data entry mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 19615080A5
Standard Cyl #? 3
Standard Expiration? 09/05/2017
Oper No? 666666

Flow Cal. Date: 08/18/2015
Slope: 663
Intercept: -631787

IR Calibration Date: 08/18/2015
3um 9um

	3um	9um
0th Coef(*100):	-24626	-16551
1st Coef(*100):	269218	130789
2nd Coef(*100):	3390	1656
H2O adj(ng/l*10K):	742	460

***** Test End *****



Operator Signature
CHARLES EDER

Remarks:

PRINT TEST

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 14:55

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 19615080A5
Standard Cyl #? 3
Standard Expiration? 09/05/2017
Oper No? 666666

Flow Cal. Date: 08/18/2015
Slope 663
Intercept -631787

IR Calibration Date: 08/18/2015
 3um 9um

0th Coef(*100):	-24626	-16551
1st Coef(*100):	269218	130789
2nd Coef(*100):	3390	1656
H2O adj(mg/l*10k):	742	460

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 14:58

Test	AC	Time
01 Room Air	0.000	14:59
02 *Subject Test	REF*	14:59
03 Room Air	0.000	15:00

*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = REFUSAL
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

REFUSAL


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 15:01

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

DIAGNOSTIC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 15:11

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:11
02 Std. Sol.	0.000	15:12
03 Room Air	0.000	15:13
04 Std. Sol.	0.000	15:13
05 Room Air	0.000	15:14
06 Std. Sol.	0.000	15:14
07 Room Air	0.000	15:15

08 Sim Temp = 34.0°C

Simul Ser No = MP3067
Std Sol No = MILLIQ H2O
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

ZERO AC
0.000 AC

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 15:17

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:18
02 Std. Sol.	0.019	15:19
03 Room Air	0.000	15:19
04 Std. Sol.	0.019	15:20
05 Room Air	0.000	15:20
06 Std. Sol.	0.019	15:21
07 Room Air	0.000	15:22

08 Sim Temp = 34.0°C

Simul Ser No = MP3066
Std Sol No = 15160
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Low AC
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 15:47

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:47
02 Std. Sol.	0.406	15:48
03 Room Air	0.000	15:49
04 Std. Sol.	0.403	15:50
05 Room Air	0.000	15:50
06 Std. Sol.	0.402	15:51
07 Room Air	0.000	15:52

08 Sim Temp = 34.0°C

Simul Ser No = MP3068
Std Sol No = 14102
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

HIGH AC
0.400 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 15:58

Test	AC	Time
01 Room Air	0.000	16:00
02 *Subject Test	INT*	16:00
03 Room Air	0.000	16:01

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = 0.10AC+0.05ACE
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *INTERFERENT DETECT*
0.10AC + 0.05% ACETONE


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 16:02

Test	AC	Time
01 Room Air	0.000	16:03
02 *Subject Test	INT*	16:04
03 Room Air	0.000	16:05

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = 0.05 ACETONE
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *INTERFERENT DETECT*
0.05% ACETONE

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 16:06

Test	AC	Time
01 Room Air	0.000	16:07
02 *Subject Test	INT*	16:07
03 Room Air	0.000	16:08

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = 0.10 METHANOL
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *INTERFERENT DETECT*
0.10% METHANOL

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 16:09

Test	AC	Time
01 Room Air	0.000	16:10
02 *Subject Test	INT*	16:11
03 Room Air	0.000	16:11

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = 0.10 IPA
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *INTERFERENT DETECT*
0.10% ISOPROPANOL

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 16:12

Test	AC	Time
01 Room Air	AMB*	16:13
02 Room Air	0.000	16:14

*Invalid Test
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = AMBIENT FAIL
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *AMBIENT FAIL*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 16:15

Test	AC	Time
01 Room Air	0.000	16:16
02 Invalid Sample	X.XXX	16:16
03 Room Air	0.000	16:17

*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INVALID SAMPLE
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

INVALID SAMPLE

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 16:19

Test	AC	Time
01 Room Air	0.000	16:20
02 *Subject Test	ABT*	16:20
03 Room Air	0.000	16:21

*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = SEQUENCE ABORT
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: SEQUENCE ABORTED

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 16:21

Test	AC	Time
01 Room Air	0.000	16:22
02 Subject Test 1	0.014	16:23
03 Room Air	0.000	16:23
04 Reported AC	0.014	16:23

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = ABA
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

ABA

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 16:26

Test	AC	Time
01 Diagnostic	OK	16:26
02 Room Air	0.000	16:27
03 Subject Test 1	0.000	16:27
04 Room Air	0.000	16:29
05 Std. Gas	0.086*	16:30
06 Room Air	0.000	16:31

*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = CAL OUT OF TOL
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature

CHARLES EDER *CHECK*

Remarks: *CALIBRATION OUT OF TOLERANCE*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 16:31

Test	AC	Time
01 Diagnostic	OK	16:33
02 Room Air	0.000	16:33
03 *Subject Test	RFI*	16:33
04 Room Air	0.000	16:34

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI DETECTED
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *RFI DETECTED*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 16:36

Test	AC	Time
01 Diagnostic	OK	16:37
02 Room Air	0.000	16:37
03 *Subject Test	0.000*	16:40
04 Room Air	0.000	16:43
05 Std. Gas	0.081	16:44
06 Room Air	0.000	16:45
07 *Subject Test	0.000*	16:49
08 Room Air	0.000	16:49


No RFI Detected

*Deficient Sample - Value Printed was
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = DEFICIENT
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

DEFICIENT
SAMPLES

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/14/2016 16:50

Test	AC	Time
01 Diagnostic	OK	16:51
02 Room Air	IPS*	16:51
03 Room Air	0.000	16:52

*Improper Sample

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = IMPROPER SAMPLE
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

IMPROPER SAMPLE

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/15/2016 09:16

Test	AC	Time
01 Diagnostic	OK	09:17
02 Room Air	0.000	09:17
03 Subject Test 1	0.191*	09:18
04 Room Air	0.000	09:21
05 Std. Gas	0.079	09:22
06 Room Air	0.000	09:23
07 Subject Test 2	0.139*	09:24
08 Room Air	0.000	09:25

*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = DIFF TOO GREAT
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *DIFFERENCE TOO GREAT*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/15/2016 09:26


Test	AC	Time
01 Diagnostic	OK	09:27
02 Room Air	0.000	09:28
03 Subject Test 1	0.139	09:29
04 Room Air	0.000	09:31
05 Std. Gas	0.078	09:32
06 Room Air	0.000	09:34
07 Subject Test 2	0.137	09:35
08 Room Air	0.000	09:36
09 Reported AC	0.137	09:35

Difference OK
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = CMS TEST
Dr. Lic. = ND/DIS821456
Lot No = 19615080A5
Cyl No = 3
Expiration Date = 09/05/2017
County = 08 Oper No. = 666666

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *CMS TEST*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = TOXL 8164.13.00 06/09
01/15/2016 09:58

Test	AC	Time
01 Room Air	0.000	09:58
02 Std. Gas	0.080	09:59
03 Room Air	0.000	09:59
04 Std. Gas	0.080	10:00
05 Room Air	0.000	10:00
06 Std. Gas	0.080	10:01
07 Room Air	0.000	10:01
08 Std. Gas	0.080	10:02
09 Room Air	0.000	10:02
10 Std. Gas	0.080	10:03
11 Room Air	0.000	10:03
12 Std. Gas	0.080	10:04
13 Room Air	0.000	10:04
14 Std. Gas	0.080	10:05
15 Room Air	0.000	10:05
16 Std. Gas	0.080	10:06
17 Room Air	0.000	10:06
18 Std. Gas	0.081	10:07
19 Room Air	0.000	10:07
20 Std. Gas	0.082	10:08
21 Room Air	0.000	10:08
22 Std. Gas	0.080	10:09
23 Room Air	0.000	10:09
24 Std. Gas	0.080	10:10
25 Room Air	0.000	10:10
26 Std. Gas	0.081	10:11
27 Room Air	0.000	10:11
28 Std. Gas	0.081	10:12
29 Room Air	0.000	10:12
30 Std. Gas	0.081	10:13
31 Room Air	0.000	10:13
32 Std. Gas	0.081	10:13
33 Room Air	0.000	10:14
34 Std. Gas	0.081	10:14
35 Room Air	0.000	10:15
36 Std. Gas	0.081	10:15
37 Room Air	0.000	10:16
38 Std. Gas	0.081	10:16
39 Room Air	0.000	10:17
40 Std. Gas	0.082	10:17
41 Room Air	0.000	10:18
42 Cal Check Stats:		
43 Average	0.0806	
44 Std Dev	0.0007	
45 Rel Std Dev	0.85%	


STABILITY TEST
Charlie E. Edle
1-15-16

Operator Signature
CHARLES EDER

Remarks:

STABILITY TEST

Form 106-I8000


1.15.16

LOT: 19615080A5
EXP. DATE: 9/5/2017
CYL: 003