



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program

SFN 59281 (5-2011)

Serial Number 80-002669	Instrument Location Watford City
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>D. J. [Signature]</i>	Date 1-18-13
Reviewed By (Crime Laboratory Use Only)	Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002669
Location = WATF 8164.13.00 06/09
01/18/2013 23:43

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 11611080A2
Standard Cyl #? 114
Standard Expiration? 06/01/2013
Oper No? 030924

Flow Cal. Date: 09/09/2011
Slope 662
Intercept -510765

IR Calibration Date: 09/09/2011
 3um 9um

0th Coef(*100):	-15623	-9015
1st Coef(*100):	263092	133867
2nd Coef(*100):	3130	1285
H2O adj(mg/l*10k):	642	392

***** Printer Test End *****



Operator Signature
DUSTIN B. NEWMAN

Remarks: Good Print Test


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002669
Location = WATF 8164.13.00 06/09
01/18/2013 23:52

Test	AC	Time
01 Room Air	RFI*	23:53
02 Room Air	RFI*	23:53

*Invalid Test
Inhibited - RFI

Sub Name = TEST, TEST TEST
Sub DOB = 03/25/1980
Sub Sex = Male Weight = 215
Test = OTH Cit = 0000000
Dr. Lic. = ND/000000000
Lot No = 11611080A2
Cyl No = 114
Expiration Date = 06/01/2013
County = 27 Oper No. = 030924



Operator Signature
DUSTIN B. NEWMAN

Remarks:

The Intoxilyzer 8000, serial number 80-002669, was not properly installed on 18 Jan 2013 at McKenzie Co SO due to an RFI test done in the wrong mode. I contacted the Field Inspector to ask for a new RFI test. He said he would run it, but I never received one. During that time frame McKenzie Co SO was extremely busy.

The Intoxilyzer 8000, 80-002669 was returned to the crime lab 16 May 2013 by Dustin Newman for the annual inspection.

Deb Kashur
Forensic Scientist
ND Attorney General, Crime Lab Division