



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59281 (03-2014)

Serial Number 80-005942	Instrument Location Turtle Mountain Agency- Belcourt, ND
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if Any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature 	Date 1-27-2015 ✓
Reviewed By (Crime Laboratory Use Only) 	Date 12 Feb 15

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005942
Location = BELC 8164.13.00 06/09
01/27/2015 09:57

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 33913080A4
Standard Cyl #? 37
Standard Expiration? 01/01/2016
Oper No? 119802

Flow Cal. Date: 06/18/2013
Slope 715
Intercept -740740

IR Calibration Date: 04/07/2014
 3um 9um

0th Coef(*100): -28513 -31653
1st Coef(*100): 267070 139194
2nd Coef(*100): 2061 961
H2O adj(mg/l*10k): 821 713

***** Printer Test End *****

 1-27-2015
Operator Signature
BRIAN F. CRISSLER

Remarks: *Print Test*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005942
Location = BELC 8164.13.00 06/09
01/27/2015 10:00

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:01
02 Std. Gas	0.078	10:01
03 Room Air	0.000	10:02
04 Std. Gas	0.080	10:02
05 Room Air	0.000	10:03
06 Std. Gas	0.079	10:03
07 Room Air	0.000	10:03

Lot No = 33913080A4
Cyl No = 37
Exp Date = 01/01/2016
County = 63 Oper No. = 119802


Operator Signature
BRIAN F. CRISSLER

Remarks: *ACA Test*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005942
Location = BELC 8164.13.00 06/09
01/27/2015 10:17

Test	AC	Time
01 Diagnostic	OK	10:18
02 Room Air	RFI*	10:18
03 Room Air	0.000	10:19

*Invalid Test
Inhibited - RFI

Sub Name = RFI TEST, NA NA
Sub DOB = 01/01/1960
Sub Sex = Unknown Weight = NA
Test = DUI Cit = 0000
Dr. Lic. = ND/NA
Lot No = 33913080A4
Cyl No = 37
Expiration Date = 01/01/2016
County = 63 Oper No. = 119802

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.


Operator Signature
BRIAN F. CRISSLER

Remarks: RFI Test