



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION


Toxicology Section/Breath Alcohol Program

SFN 59282 (10-2013)

| | | |
|---|---|-----------------------------------|
| Chemical Test Operator Name (Print) Thompson, Chad A | | |
| Location Cass Jail | Intoxilyzer® Serial Number 80-005939 | |
| Gas Lot Number 14214080A6 | Gas Cylinder Number 49 | Gas Expiration Date 07/05/2016 |

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.083 AC
 - B. 0.084 AC
 - C. 0.084 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

| | |
|--|-------------------|
| Chemical Test Operator Signature  3524 | Date 1/30/15 |
| Reviewed By (Crime Laboratory Use Only) Deb Kashner | Date 04 Feb 15 |


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005939
Location = CASS 8164.13.00 06/09
01/30/2015 03:44

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 03:45 |
| 02 Std. Gas | 0.083 | 03:45 |
| 03 Room Air | 0.000 | 03:46 |
| 04 Std. Gas | 0.084 | 03:46 |
| 05 Room Air | 0.000 | 03:47 |
| 06 Std. Gas | 0.084 | 03:47 |
| 07 Room Air | 0.000 | 03:48 |

Lot No = 14214080A6
Cyl No = 49
Exp Date = 07/05/2016
County = 09 Oper No. = 050910

 3524

Operator Signature
CHAD THOMPSON

Remarks:

ACA Test