

INTOXILYZER 8000
REPAIR AND MAINTENANCE

Serial No.: 80-003067

Location: TOXL

- | | Check When Done |
|--|-----------------|
| 1. Is a three-pronged grounded outlet used? | <u>✓</u> |
| 2. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| 3. Is the breath tube heated? | <u>✓</u> |
| 4. Is the diagnostic check complete? | <u>✓</u> |
| 5. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| 6. Print test. (Attach test record.) | <u>✓</u> |
| 7. Calibration Check. Use Ethanol Breath Standard (EBS). Do three sets in ACA mode. (Attach test records.) | |

Lot No. 23411080 A Cylinder No. 10 Expir. Date 10/01/2013

Test 1 0.078 Test 4 0.081 Test 7 0.080

Test 2 0.078 Test 5 0.079 Test 8 0.079

Test 3 0.078 Test 6 0.080 Test 9 0.081

Average 0.079

- | | |
|--|----------|
| 8. Low AC. Use < 0.03 AC in ACA mode. (Attach test record.) | <u>✓</u> |
| 9. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.) | <u>✓</u> |
| 10. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode. (Attach test record.) | <u>✓</u> |
| 11. RFI check. Run CMS mode. Key radio on first room air. (Attach test record.) | <u>✓</u> |

Electrical Repairs:

1. Was the calibration adjusted/performed? (Attach test record.) yes, @ CMI
2. Was the light source changed? NO
3. Was the sample chamber cleaned? NO

Other Maintenance/Repairs:

1. Was the printer cleaned/repaired? NO
2. Was the solenoid cleaned? Replaced
3. Was new tubing installed? NO
4. Were additional parts replaced/changed/repaired? yes, see below
5. Remarks/Maintenance record: CMI replaced the check valve, battery, and solenoid.

01 Feb 2013
Date

01 FEB 2013
Date

Deb Kashur
Field Inspector's Signature

Charles E. Edin
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = TOXL 8164.13.00 06/09
02/01/2013 11:02

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:02
02 Std. Gas	0.078	11:03
03 Room Air	0.000	11:03
04 Std. Gas	0.078	11:04
05 Room Air	0.000	11:04
06 Std. Gas	0.078	11:05
07 Room Air	0.000	11:05

Lot No = 23411080A1
Cyl No = 10
Exp Date = 10/01/2013
County = 08 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = TOXL 8164.13.00 06/09
02/01/2013 11:08

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:09
02 Std. Gas	0.081	11:09
03 Room Air	0.000	11:10
04 Std. Gas	0.079	11:10
05 Room Air	0.000	11:11
06 Std. Gas	0.080	11:11
07 Room Air	0.000	11:12

Lot No = 23411080A1

Cyl No = 10

Exp Date = 10/01/2013

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = TOXL 8164.13.00 06/09
02/01/2013 11:12

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:12
02 Std. Gas	0.080	11:13
03 Room Air	0.000	11:13
04 Std. Gas	0.079	11:13
05 Room Air	0.000	11:14
06 Std. Gas	0.081	11:14
07 Room Air	0.000	11:15

Lot No = 23411080A1
Cyl No = 10
Exp Date = 10/01/2013
County = 08 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:
Calibration check 0.080 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = TOXL 8164.13.00 06/09
02/01/2013 13:23

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:23
02 Std. Sol.	0.020	13:24
03 Room Air	0.000	13:25
04 Std. Sol.	0.020	13:25
05 Room Air	0.000	13:26
06 Std. Sol.	0.020	13:27
07 Room Air	0.000	13:27

08 Sim Temp = 34.0°C

Simul Ser No = DR5132
Std Sol No = 12210
County = 08 Oper No. = 777777



Operator Signature
DEB KASHUR

Remarks:

Low AC 0.020 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = TOXL 8164.13.00 06/09
02/01/2013 13:56

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:57
02 Std. Sol.	0.297	13:58
03 Room Air	0.000	13:58
04 Std. Sol.	0.298	13:59
05 Room Air	0.000	14:00
06 Std. Sol.	0.298	14:00
07 Room Air	0.000	14:01

08 Sim Temp = 34.0°C

Simul Ser No = DR3803

Std Sol No = 12050

County = 08

Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Linearity test 0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = TOXL 8164.13.00 06/09
02/01/2013 14:02

Test	AC	Time
01 Room Air	0.000	14:03
02 *Subject Test	INT*	14:04
03 Room Air	0.000	14:04

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 23411080A1
Cyl No = 10
Expiration Date = 10/01/2013
County = 08 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Interferent check 0.10 Ac + acetone

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = TOXL 8164.13.00 06/09
02/01/2013 13:42

Test	AC	Time
01 Diagnostic	OK	13:43
02 Room Air	RFI*	13:43
03 Room Air	0.000	13:43

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 23411080A1
Cyl No = 10
Expiration Date = 10/01/2013
County = 08 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
DEB KASHUR

Remarks:

RFI check

Form 106-I8000