



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59281 (5-2011)

Serial Number 80-004953	Instrument Location New Town
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature Darryl W. Abate	Date 02/03/2012
Reviewed By (Crime Laboratory Use Only) Deb Kashner	Date 03 Apr 2012

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004953
Location = NEWT 8164.13.00 06/09
02/03/2012 19:30

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-==|
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-==|
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 23411080A1
Standard Cyl #? 03
Standard Expiration? 10/01/2013
Oper No? 010101

Flow Cal. Date: 01/19/2010
Slope 691
Intercept -628068

IR Calibration Date: 06/08/2010
 3um 9um

0th Coef(*100): -10343 -18245
1st Coef(*100): 249276 131289
2nd Coef(*100): 4014 1439
H2O adj(mg/l*10k): 697 655

***** Printer Test End *****



Operator Signature
DARCY W. ABERLE

Remarks:


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004953
Location = NEWT 8164.13.00 06/09
02/03/2012 19:30

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:31
02 Std. Gas	0.079	19:32
03 Room Air	0.000	19:32
04 Std. Gas	0.080	19:32
05 Room Air	0.000	19:33
06 Std. Gas	0.079	19:33
07 Room Air	0.000	19:34

Lot No = 23411080A1
Cyl No = 03
Exp Date = 10/01/2013
County = 31 Oper No. = 010101


Operator Signature
DARCY W. ABERLE

Remarks:


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004953
Location = NEWT 8164.13.00 06/09
02/03/2012 19:36

Test	AC	Time
01 Room Air	0.000	19:39
02 Reference	RFI*	19:39
03 Room Air	0.000	19:39

*Invalid Test
 Inhibited - RFI

Sub Name = TEST, TEST T
Sub DOB = 12/12/1980
Sub Sex = Male Weight = 200
Test = OTH Cit = NA
Dr. Lic. = ND/NA
Lot No = 23411080A1
Cyl No = 03
Expiration Date = 10/01/2013
County = 31 Oper No. = 010101



Operator Signature
DARCY W. ABERLE

Remarks: