

INTOXILYZER 5000 G and 8000 INSPECTION

Serial No.: 80-004953

Location: TOX LAB

- Check When Done
- A. Is the warm-up time less than 20 minutes? ✓
  - B. Is a three-pronged grounded outlet used? ✓
  - C. Is the breath tube heated? ✓
  - D. Is the diagnostic check complete? ✓
  - E. Is the time, date, and year correct (re-set if necessary)? ✓
  - F. Print test. (Attach test record.) ✓
  - G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)  
Sim SN DR5190 Lot No. — Rep AC 0.021 AC ✓
  - H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)  
Sim SN DR5113 Lot No. — Rep AC 0.285 AC ✓
  - I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.  
(Attach test record.)  
Sim SN DR5135 Lot No. — Rep AC 0.10 AC + 0.05% Acetone ✓
  - J. RFI check. Run CMS mode. Key radio on first room air.  
(Attach test record.) ✓

- K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.  
(Attach test record.)

Lot No. ~~568~~ <sup>MRP</sup> 668180 Cyl No. 94 Exp. Date 9/03/11

Test 1	<u>0.079</u>	Test 4	<u>0.080</u>	Test 7	<u>0.080</u>
Test 2	<u>0.080</u>	Test 5	<u>0.080</u>	Test 8	<u>0.080</u>
Test 3	<u>0.080</u>	Test 6	<u>0.080</u>	Test 9	<u>0.081</u>

Average 0.080 AC

- L. Remarks/Maintenance record: Annual Inspection
- 8 February 2010  
Date DK 12 Feb 2010  
~~He Feb 2010~~  
Date
- Mary Pearson  
Field Inspector's Signature
- Deb Kashner  
Reviewed by

80-004953  
25/33 FM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004953  
Location = TOXL      8164.13.00 06/09  
02/09/2010      14:54

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup


Data Entry Mode:      Enabled  
Start Test Sequence:      DABACABA  
Display Prelim Rslt?      Yes  
Display Third Digit?      Yes  
Inhib Printer(Y/N)?      No  
Display Volume?      No  
Disable On Memfull?      No  
# of Print Copies?      1  
Select Std (D/W/I)?      Wet  
Standard Value?      0.000  
Standard Solution #?      00000  
Sim Serial #?      000000  
Oper No?      107501

Flow Cal. Date:      01/19/2010  
Slope      691  
Intercept      -628068

IR Calibration Date:      02/09/2010  
                                3um                  9um

-----  
0th Coef(\*100):      -10827                  -9420  
1st Coef(\*100):      244917                  123517  
2nd Coef(\*100):      1933                      1313  
H2O adj(mg/l\*10k): 688                      649

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
MARGY PEARSON

Remarks:

Print

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004953  
Location = TOXL      8164.13.00 06/09  
02/09/2010      14:33

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:34
02 Std. Sol.	0.019	14:35
03 Room Air	0.000	14:36
04 Std. Sol.	0.020	14:36
05 Room Air	0.000	14:37
06 Std. Sol.	0.020	14:38
07 Room Air	0.000	14:38


08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 00000

County = 08

Oper No. = 107501



Operator Signature  
MARGY PEARSON

Remarks:

LOW AC 0.021 AC

Form 106-I8000

80-004953  
27/33 KAL

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004953  
Location = TOXL      8164.13.00 06/09  
02/09/2010      14:39

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:40
02 Std. Sol.	0.278	14:42
03 Room Air	0.000	14:42
04 Std. Sol.	0.278	14:43
05 Room Air	0.000	14:43
06 Std. Sol.	0.280	14:44
07 Room Air	0.000	14:45

08 Sim Temp = 34.0°C

Simul Ser No = 000000  
Std Sol No = 00000  
County = 08      Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

*Linearity 0.285 AC*

Form 106-I8000

80-004953  
28/33 PHL


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004953  
Location = TOXL      8164.13.00 06/09  
02/09/2010      14:57

Test	AC	Time
01 Room Air	0.000	14:58
02 *Subject Test	INT*	14:59
03 Room Air	0.000	14:59

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 11/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS821456  
Simul Ser No = 000000  
Std Sol No = 00000  
County = 08      Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

*Interfer - Acetone + Alcohol*

Form 106-18000

80-004953  
29/33 1572

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004953  
Location = TOXL      8164.13.00 06/09  
02/09/2010      15:00

Test	AC	Time
01 Diagnostic	OK	15:01
02 Room Air	RFI*	15:01
03 Room Air	0.000	15:01

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 11/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS821456  
Simul Ser No = 000000  
Std Sol No = 00000  
County = 08      Oper No. = 107501

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator/Signature  
MARGY PEARSON

Remarks:

RFI

Form 106-I8000

80-004953  
30/33 VHL


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004953  
Location = TOXL      8164.13.00 06/09  
02/09/2010      15:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:20
02 Std. Gas	0.079	15:21
03 Room Air	0.000	15:21
04 Std. Gas	0.080	15:22
05 Room Air	0.000	15:22
06 Std. Gas	0.080	15:22
07 Room Air	0.000	15:23

Lot No = 659358  
Cyl No = 94  
Exp Date = 09/03/2011  
County = 08      Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

*ACA*

Form 106-I8000

80-004953  
31/33 YAL

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004953  
Location = TOXL      8164.13.00 06/09  
02/09/2010      15:26

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:26
02 Std. Gas	0.080	15:27
03 Room Air	0.000	15:27
04 Std. Gas	0.080	15:27
05 Room Air	0.000	15:28
06 Std. Gas	0.080	15:28
07 Room Air	0.000	15:29


Lot No = 659358

Cyl No = 94

Exp Date = 09/03/2011

County = 08

Oper No. = 107501



Operator Signature  
MARGY PEARSON

Remarks:

ACA

Form 106-I8000

80-004953  
32/33 KAL




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004953  
Location = TOXL      8164.13.00 06/09  
02/09/2010      15:29

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:30
02 Std. Gas	0.080	15:30
03 Room Air	0.000	15:31
04 Std. Gas	0.080	15:31
05 Room Air	0.000	15:32
06 Std. Gas	0.081	15:32
07 Room Air	0.000	15:33

Lot No = 659358  
Cyl No = 94  
Exp Date = 09/03/2011  
County = 08      Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

ACA

Form 106-I8000

80-004953  
33/33 192