

## Intoxilyzer 8000 Initial Inspection

**Serial Number** 80-004949

- A. Power On 5 min. countdown 1 min.
- B. Security Checks. Level 1 Level 2 Level 3
- C. Gas Cylinder Setup
- D. Printer. Level 1 Func P Paper Feed Internal External
- E. Display is working.
- F. Configure Start Test Switch. Level 3 Func S Sub T
- G. Setup Preliminary Data. Level 2 Func E
- H. Configure General Setup. Level 3 Func S Sub G.
- I. Configure Port. Level 3 Func S Sub N
- J. Check Camera. ABA Mode Oper ID ND Lic. Hot Key "R"
- K. Diagnostic Test. Level 2 Func D
- L. Check Ser Num / Version Level 2 Func V  
80-004949 Serial No. 816413 Version 136214 Display 7 CPLD
- M. Detector Voltages. Level 3 Func D Sub D  
3 12618 s 1 C 4700  
9 13011 s 3 B 44 0 RFI Ind.
- N. Hot Keys  
F5 Tone  
F6 Green LED  
F7 Red LED  
F9 Fan  
F10 Pump  
F11 Sol Output
- O. RFI Dectector
- P. Atmospheric Monitor. Level 3 Func D Sub A. P 1026 CF 0.98

08/09

G:\Crime Lab\Manuals\ISO SOP\Toxicology Section\Breath Alcohol Program\Instrument Certification and Maintenance Unit\Intoxilyzer 8000\Forms-Worksheets\Form BrF-005-8 Intoxilyzer 8000 Initial Inspection 08-09

80-004949  
4/24/12

OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 Instrument Certification and Maintenance Unit

BrF-005-8

Q. Flow Monitor. Level 3 Func D Sub F. R 214 F 0.000 L/s  
 R. Temperature Monitor. Level 3 Func D Sub T  
 C 47.0 A 47.5 I 0.0  
 B 44.9 V 0.0  
 S. Gas Tank Monitor. Level 3 Func D Sub G  
 T. Configure Wet Bath Sim. Level 1 Func S  
 U. Zero Alcohol Test Level 1 Func C  
DR3409 Sim SN        Lot 0.000 Rep Value  
0.000 AC 0.000 AC 0.000 AC  
 V. Low AC Test. Level 1 Func C  
DR5144 Sim SN        Lot 0.019 Rep Value  
0.018 AC 0.018 AC 0.018 AC  
 W. Linearity Test. Level 1 Func C  
DR5187 Sim SN        Lot 0.285 Rep Value  
0.285 AC 0.284 AC 0.284 AC  
 X. Interferent Test.  
 Level 1 Func B  
 Acetone DR3845 Sim SN        Lot 0.05% Conc  
 Display  Printout  
 Methanol 10573 Sim SN        Lot 0.10% Conc  
 Display  Printout  
 IPA DR5191 Sim SN        Lot 0.10% Conc  
 Display  Printout  
 Y. Ambient Test  
 Z. Invalid Sample Test  
 AA. RFI Detector Test. CMS Mode  
 Display  Printout  
 BB. Configure Dry Gas Standard. Level 1 Func S  
 CC. Stability Test at 0.080AC. Level 3 Func T Sub S  
0.0801 Ave. 0.0003 Std Dev 0.39% Rel Std Dev  
 DD. Maintenance/Repairs Recalibrated.

10 Feb 2010  
 Date  
15 Feb 2010  
 Date  
 MAP

Deb Kashur  
 Signature  
Mary Pearson  
 Reviewed by

80-004999  
 7/24/10

08/09

TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-004949  
 02/09/2010 14:20:47

Auto Calibration

pg 1 of 2

	<<<<<	3um	>>>>>	<<<<<	9um	>>>>>
-----						
Solution =	0.000 g/210L	or	0.0000 mg/l,	Samples = 4,	Discarded = 1	
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)	
Sample #1	0.1160	(0.0020)		0.2840	(-0.0080)	
Sample #2	0.1080	(0.0740)		0.2660	(0.0280)	
Sample #3	0.1380	(0.0990)		0.3010	(0.0260)	
Sample #4	0.1330	(0.1440)		0.2840	(0.0640)	
Avg % Abs	0.1263	(0.1057)		0.2837	(0.0393)	
STD DEV	0.0161	(0.0355)		0.0175	(0.0214)	
REL STD DEV	12.723	(33.571)		6.170	(54.370)	
-----						
Solution =	0.040 g/210L	or	0.1905 mg/l,	Samples = 4,	Discarded = 1	
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)	
Sample #1	0.8040	(0.0120)		1.6430	(0.0000)	
Sample #2	0.8180	(0.0220)		1.6680	(-0.0070)	
Sample #3	0.8340	(0.0380)		1.6660	(0.0030)	
Sample #4	0.8150	(0.0490)		1.6520	(0.0170)	
Avg % Abs	0.8223	(0.0363)		1.6620	(0.0043)	
STD DEV	0.0102	(0.0136)		0.0087	(0.0121)	
REL STD DEV	1.242	(37.368)		0.525	(278.202)	
-----						
Solution =	0.107 g/210L	or	0.5095 mg/l,	Samples = 4,	Discarded = 1	
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)	
Sample #1	2.0950	(-0.0240)		4.0890	(-0.0150)	
Sample #2	2.0580	(-0.0030)		4.0550	(0.0070)	
Sample #3	2.0710	(0.0080)		4.0610	(0.0130)	
Sample #4	2.0220	(0.0330)		4.0380	(0.0180)	
Avg % Abs	2.0503	(0.0127)		4.0513	(0.0127)	
STD DEV	0.0254	(0.0184)		0.0119	(0.0055)	
REL STD DEV	1.238	(145.643)		0.294	(43.481)	
-----						
Solution =	0.147 g/210L	or	0.7000 mg/l,	Samples = 4,	Discarded = 1	
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)	
Sample #1	2.7330	(-0.0100)		5.3600	(0.0200)	
Sample #2	2.7170	(0.0230)		5.3570	(0.0460)	
Sample #3	2.7480	(0.0140)		5.3790	(0.0320)	
Sample #4	2.7300	(0.0080)		5.3570	(0.0420)	
Avg % Abs	2.7317	(0.0150)		5.3643	(0.0400)	
STD DEV	0.0156	(0.0075)		0.0127	(0.0072)	
REL STD DEV	0.570	(50.332)		0.237	(18.028)	
-----						
Solution =	0.285 g/210L	or	1.3571 mg/l,	Samples = 4,	Discarded = 1	
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)	
Sample #1	5.2410	(-0.0170)		10.0310	(0.0020)	
Sample #2	5.2070	(0.0220)		10.0290	(0.0270)	
Sample #3	5.2460	(0.0100)		10.0760	(0.0260)	
Sample #4	5.2650	(0.0070)		10.0840	(0.0360)	
Avg % Abs	5.2393	(0.0130)		10.0630	(0.0297)	
STD DEV	0.0296	(0.0079)		0.0297	(0.0055)	
REL STD DEV	0.564	(61.056)		0.295	(18.565)	
-----						

*Calibration  
 Deb Kashan  
 p. 1 of 2*

*80-004949  
 8/24 PM*

TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-004949  
 PPPP 0>3A  
 Zero Order Coef -329.45 -372.95  
 PP

Act (g/210L)	Fit (g/210L)	Residual (g/210L)	Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0002	0.000	0.000	-0.0002
0.040	0.039	0.0006	0.040	0.039	0.0006
0.107	0.108	-0.0012	0.107	0.108	-0.0011
0.147	0.146	0.0008	0.147	0.146	0.0008
0.285	0.285	-0.0001	0.285	0.285	-0.0000

<<<< 3um >>>> <<<< 9um >>>>

Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1

Sample	3um	9um
Sample #1	3208.00	3117.00
Sample #2	3200.00	3156.00
Sample #3	3112.00	3142.00
Sample #4	3213.00	3174.00
Avg	3175.0000	3157.3333
STD DEV	54.9454	16.0416
REL STD DEV	1.731	0.508
H2O adjust (mg/l*10k)	634	652

Atmospheric Pressure = 1026

\*\*\*\*\*CALIBRATION SUCCESSFUL\*\*\*\*\*

*Calibration  
 Deb Kashner  
 p. 2 of 2*

80-004949  
 9/24 FM

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-004949  
Location = TOXL 8164.13.00 06/09  
02/09/2010 15:54

\*\*\*\*\* Printer Test \*\*\*\*\*

12 char/line  
abcdefghijkl  
ABCDEFGHIJKL

16 char per line  
abcdefghijklmnop  
ABCDEFGHIJKLMN

21 chars per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMNQRSTU

24 characters per line  
abcdefghijklmnopqrstuwx  
ABCDEFGHIJKLMNQRSTUWX

32 characters per line  
abcdefghijklmnopqrstuwx123456  
ABCDEFGHIJKLMNQRSTUWXYZ123456

42 characters per line  
abcdefghijklmnopqrstuwx1234567890abcdef  
ABCDEFGHIJKLMNQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup

Data Entry Mode: Enabled  
Start Test Sequence: DABACABA  
Display Prelim Rslt? Yes  
Display Third Digit? Yes  
Inhib Printer(Y/N)? No  
Display Volume? No  
Disable On MemFull? No  
# of Print Copies? 1  
Select Std (D/W/I)? Dry  
Standard Value? 0.080  
Standard Lot #? 668180  
Standard Cyl #? 54  
Standard Expiration? 09/21/2011  
Oper No? 777777

Flow Cal. Date: 01/18/2010  
Slope 771  
Intercept -891468

IR Calibration Date: 02/09/2010

	3um	9um
0th Coef(*100):	-32945	-37295
1st Coef(*100):	268565	134695
2nd Coef(*100):	-610	387
H2O adj(ng/l*10K):	634	652

\*\*\*\*\* Test End \*\*\*\*\*

*Deb Keshur*

Operator Signature  
DEB KASHUR

Remarks:

*Print Test*

Form 106-18000

80-004949  
10/24 RM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000       SN 80-004949  
Location = TOXL                8164.13.00 06/09  
02/09/2010                      15:55

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup  
Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?         Yes  
Inhib Printer(Y/N)?         No  
Display Volume?             No  
Disable On Memfull?         No  
# of Print Copies?          1  
Select Std (D/W/I)?         Dry  
Standard Value?             0.080  
Standard Lot #?             668180  
Standard Cyl #?             54  
Standard Expiration?        09/21/2011  
Oper No?                    777777

Flow Cal. Date:             01/18/2010  
Slope                        771  
Intercept                    -891468

IR Calibration Date:        02/09/2010

	3um	9um
0th Coef(*100):	-32945	-37295
1st Coef(*100):	268565	134695
2nd Coef(*100):	-610	387
H2O adj(mg/l*10k):	634	652

\*\*\*\*\* Printer Test End \*\*\*\*\*

Deb Kasher  
Operator Signature  
DEB KASHUR

Remarks:

*Print Test*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004949  
Location = TOXL                8164.13.00 06/09  
02/09/2010                        16:03

Test	AC	Time
01 Room Air	0.000	16:04
02 *Subject Test	REF*	16:04
03 Room Air	0.000	16:05

\*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 668180  
Cyl No = 54  
Expiration Date = 09/21/2011  
County = 08                          Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*Camera/Refusal check*

Form 106-I8000

80-004949  
10/24 PM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004949  
Location = TOXL                    8164.13.00 06/09  
02/09/2010                            16:05

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 777777

Deb Kashur  
Operator Signature  
DEB KASHUR

Remarks:

*Diagnostic Test*

Form 106-I8000

*80-004949  
13/24 AK*



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004949  
Location = TOXL      8164.13.00 06/09  
02/09/2010      16:12

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:12
02 Std. Sol.	0.000	16:13
03 Room Air	0.000	16:13
04 Std. Sol.	0.000	16:14
05 Room Air	0.000	16:15
06 Std. Sol.	0.000	16:15
07 Room Air	0.000	16:16

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*zero alcohol      0.000 AC*

Form 106-I8000

80-004949  
14/24 RAR

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004949  
Location = TOXL      8164.13.00 06/09  
02/09/2010      16:24

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:25
02 Std. Sol.	0.018	16:25
03 Room Air	0.000	16:26
04 Std. Sol.	0.018	16:27
05 Room Air	0.000	16:27
06 Std. Sol.	0.018	16:28
07 Room Air	0.000	16:28

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*Low AC 0.018 AC*

Form 106-18000

*80-004949  
15/24 FK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004949  
Location = TOXL      8164.13.00 06/09  
02/09/2010      16:30

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:30
02 Std. Sol.	0.285	16:31
03 Room Air	0.000	16:31
04 Std. Sol.	0.284	16:32
05 Room Air	0.000	16:33
06 Std. Sol.	0.284	16:33
07 Room Air	0.000	16:34

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*Linearity test 0.285AC*

Form 106-I8000

80-004949  
10/24 FM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004949  
Location = TOXL                    8164.13.00 06/09  
02/10/2010                            08:17

Test	AC	Time
01 Room Air	0.000	08:18
02 *Subject Test	INT*	08:18
03 Room Air	0.000	08:19

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 668180  
Cyl No = 54  
Expiration Date = 09/21/2011  
County = 08                            Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*acetone interferent*

Form 106-I8000

80-004949  
17/04 YTH

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004949  
Location = TOXL                    8164.13.00 06/09  
02/10/2010                            08:21

Test	AC	Time
01 Room Air	0.000	08:22
02 *Subject Test	INT*	08:22
03 Room Air	0.000	08:23

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 668180  
Cyl No = 54  
Expiration Date = 09/21/2011  
County = 08                            Oper No. = 777777

Deb Kashur  
Operator Signature  
DEB KASHUR

Remarks:

*Methanol interferent*

Form 106-I8000

80-004949  
18/24 KAR

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004949  
Location = TOXL                8164.13.00 06/09  
02/10/2010                        08:26

Test	AC	Time
01 Room Air	0.000	08:26
02 *Subject Test	INT*	08:27
03 Room Air	0.000	08:27

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                Weight = 150  
Test = DUI                        Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 668180  
Cyl No = 54  
Expiration Date = 09/21/2011  
County = 08                        Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*IPA interferent*

Form 106-I8000

80-004949  
19/24 FM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004949  
Location = TOXL      8164.13.00 06/09  
02/09/2010      16:40

Test	AC	Time
01 Room Air	AMB*	16:41
02 Room Air	0.000	16:42

\*Invalid Test  
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS921456  
Simul Ser No = 00000  
Std Sol No = 000  
County = 08      Oper No. = 777777

Deb Kashur  
Operator Signature  
DEB KASHUR

Remarks:

*Ambient test*

Form 106-I8000

80-004949  
20/24 KM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004949  
Location = TOXL                    8164.13.00 06/09  
02/09/2010                            16:44

Test	AC	Time
01 Room Air	0.000	16:44
02 Invalid Sample	X.XXX	16:45
03 Room Air	0.000	16:46

\*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = NA  
Dr. Lic. = ND/DIS921456  
Simul Ser No = 00000  
Std Sol No = 000  
County = 08                            Oper No. = 777777

Deb Kashur  
Operator Signature  
DEB KASHUR

Remarks:

*Invalid sample test*

Form 106-I8000

80-004949  
2/24 PM



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004949  
Location = TOXL                8164.13.00 06/09  
02/10/2010                        08:29

Test	AC	Time
01 Diagnostic	OK	08:30
02 Room Air	RFI*	08:30
03 Room Air	0.000	08:31

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                Weight = 150  
Test = DUI                        Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 668180  
Cyl No = 54  
Expiration Date = 09/21/2011  
County = 08                        Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Kashur  
Operator Signature  
DEB KASHUR

Remarks:

*RFI check*

Form 106-I8000

*80-004949  
22/24 KM*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004949  
Location = TOXL                8164.13.00 06/09  
02/09/2010                      15:09

Test	AC	Time
01 Room Air	0.000	15:10
02 Std. Gas	0.080	15:10
03 Room Air	0.000	15:11
04 Std. Gas	0.080	15:11
05 Room Air	0.000	15:12
06 Std. Gas	0.080	15:12
07 Room Air	0.000	15:12
08 Std. Gas	0.080	15:13
09 Room Air	0.000	15:13
10 Std. Gas	0.080	15:14
11 Room Air	0.000	15:14
12 Std. Gas	0.080	15:15
13 Room Air	0.000	15:15
14 Std. Gas	0.080	15:15
15 Room Air	0.000	15:16
16 Std. Gas	0.080	15:16
17 Room Air	0.000	15:17
18 Std. Gas	0.081	15:17
19 Room Air	0.000	15:18
20 Std. Gas	0.080	15:18
21 Room Air	0.000	15:18
22 Std. Gas	0.080	15:19
23 Room Air	0.000	15:19
24 Std. Gas	0.081	15:20
25 Room Air	0.000	15:20
26 Std. Gas	0.080	15:21
27 Room Air	0.000	15:21
28 Std. Gas	0.080	15:21
29 Room Air	0.000	15:22
30 Std. Gas	0.080	15:22
31 Room Air	0.000	15:23
32 Std. Gas	0.080	15:23
33 Room Air	0.000	15:24
34 Std. Gas	0.080	15:24
35 Room Air	0.000	15:25
36 Std. Gas	0.080	15:25
37 Room Air	0.000	15:25
38 Std. Gas	0.080	15:26
39 Room Air	0.000	15:26
40 Std. Gas	0.080	15:27
41 Room Air	0.000	15:27
42 Cal Check Stats:		
43 Average	0.0801	
44 Std Dev	0.0003	
45 Rel Std Dev	0.39%	

*Stability test  
Deb Kashner  
p. 1 of 2*

*80-004949  
23/24 FN*

Operator Signature  
DEB KASHUR

Remarks:

Lot #: 668180 Cyl.#: 54 Exp. date: 09/21/2011

Form 106-I8000

Stability test  
Deb Kashur  
P. 2012

80-004949  
2/24/12