

**INTOXILYZER 5000 G and 8000 INSPECTION**

Serial No.: 80-004946

Location: TOX LAB

- |   | Check When Done |
|---|-----------------|
| A. Is the warm-up time less than 20 minutes?  | <u>✓</u>        |
| B. Is a three-pronged grounded outlet used?   | <u>✓</u>        |
| C. Is the breath tube heated?   | <u>✓</u>        |
| D. Is the diagnostic check complete?  | <u>✓</u>        |
| E. Is the time, date, and year correct (re-set if necessary)?   | <u>✓</u>        |
| F. Print test. (Attach test record.)  | <u>✓</u>        |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5190</u> Lot No. <u>—</u> Rep AC <u>0.021 AC</u>  | <u>✓</u>        |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5113</u> Lot No. <u>—</u> Rep AC <u>0.285 AC</u>   | <u>✓</u>        |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR5135</u> Lot No. <u>—</u> Rep AC <u>0.10% + 0.05% Acetone</u> | <u>✓</u>        |
| J. RFI check. Run CMS mode. Key radio on first room air.<br>(Attach test record.)   | <u>✓</u>        |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.)   |                 |

Lot No. 659358 Cyl No. 15 Exp. Date 9/3/11

Test 1 0.079 Test 4 0.078 Test 7 0.080

Test 2 0.079 Test 5 0.078 Test 8 0.079

Test 3 0.079 Test 6 0.078 Test 9 0.079

Average 0.079

L. Remarks/Maintenance record: Annual Inspection

11 February 2010

Date DK

16 Feb 2010 12 Feb 2010

Date

Mary Besum

Field Inspector's Signature

Deb Kashur

Reviewed by

80-004946

25/33 JHL

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004946  
Location = TOXL      8164.13.00 06/09  
02/10/2010      10:23

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup


Data Entry Mode:      Enabled  
Start Test Sequence:      DABACABA  
Display Prelim Rslt?      Yes  
Display Third Digit?      Yes  
Inhib Printer(Y/N)?      No  
Display Volume?      No  
Disable On Memfull?      No  
# of Print Copies?      1  
Select Std (D/W/I)?      Dry  
Standard Value?      0.080  
Standard Lot #?      12345678  
Standard Cyl #?      1234567  
Standard Expiration?      11/11/2011  
Oper No?      107501

Flow Cal. Date:      01/18/2010  
Slope      717  
Intercept      -728279

IR Calibration Date:      01/19/2010

	3um	9um
0th Coef(*100):	-14511	-21944
1st Coef(*100):	272115	138613
2nd Coef(*100):	1283	1052
H2O adj(mg/l*10k):	361	409

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
MARGY PEARSON

Remarks:

*Print*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004946  
Location = TOXL      8164.13.00 06/09  
02/10/2010      11:59

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:59
02 Std. Sol.	0.020	12:00
03 Room Air	0.000	12:01
04 Std. Sol.	0.019	12:02
05 Room Air	0.000	12:02
06 Std. Sol.	0.020	12:03
07 Room Air	0.000	12:03


08 Sim Temp = 34.0°C

Simul Ser No = 0000000

Std Sol No = 000000

County = 08

Oper No. = 107501



Operator Signature  
MARGY PEARSON

Remarks:

Low 0.021 AC

Form 106-I8000

80-004946  
27/33 PH

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004946  
Location = TOXL      8164.13.00 06/09  
02/10/2010      12:10

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:11
02 Std. Sol.	0.281	12:12
03 Room Air	0.000	12:13
04 Std. Sol.	0.282	12:13
05 Room Air	0.000	12:14
06 Std. Sol.	0.282	12:15
07 Room Air	0.000	12:15


08 Sim Temp = 34.0°C

Simul Ser No = 0000000

Std Sol No = 000000

County = 08

Oper No. = 107501



Operator Signature

MARGY PEARSON

Remarks:

Linearity 0.285 AC

Form 106-I8000

80-004946  
28/33 194


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004946  
Location = TOXL      8164.13.00 06/09  
02/11/2010      16:08

Test	AC	Time
01 Room Air	0.000	16:09
02 *Subject Test	INT*	16:09
03 Room Air	0.000	16:10

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 11/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS821456  
Lot No = 659358  
Cyl No = 15  
Expiration Date = 09/03/2011  
County = 08      Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

*Inter-Alc + Acetone*

Form 106-I8000

80-004946  
29/33 FM


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004946  
Location = TOXL                8164.13.00 06/09  
02/11/2010                        16:11

Test	AC	Time
01 Room Air	RFI*	16:12
02 Room Air	0.000	16:12

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 11/01/1982  
Sub Sex = Female                Weight = 150  
Test = DUI                        Cit = NA  
Dr. Lic. = ND/DIS821456  
Lot No = 659358  
Cyl No = 15  
Expiration Date = 09/03/2011  
County = 08                        Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

*RFI*

Form 106-I8000

80-004946  
30/33 RFL


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004946  
Location = TOXL      8164.13.00 06/09  
02/11/2010      14:44

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:45
02 Std. Gas	0.079	14:45
03 Room Air	0.000	14:46
04 Std. Gas	0.079	14:46
05 Room Air	0.000	14:47
06 Std. Gas	0.079	14:47
07 Room Air	0.000	14:48

Lot No = 659358  
Cyl No = 15  
Exp Date = 09/03/2011  
County = 08      Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

ACA

Form 106-I8000

80-004946  
31/33 JH


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004946  
Location = TOXL      8164.13.00 06/09  
02/11/2010      14:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:51
02 Std. Gas	0.078	14:51
03 Room Air	0.000	14:52
04 Std. Gas	0.078	14:52
05 Room Air	0.000	14:53
06 Std. Gas	0.078	14:53
07 Room Air	0.000	14:53

Lot No = 659358  
Cyl No = 15  
Exp Date = 09/03/2011  
County = 08      Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

*Cal check - A&A*

Form 106-I8000

80-004946  
32/33 KHL




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004946  
Location = TOXL      8164.13.00 06/09  
02/11/2010      14:59

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:00
02 Std. Gas	0.080	15:00
03 Room Air	0.000	15:01
04 Std. Gas	0.079	15:01
05 Room Air	0.000	15:02
06 Std. Gas	0.079	15:02
07 Room Air	0.000	15:02

Lot No = 659358  
Cyl No = 15  
Exp Date = 09/03/2011  
County = 08      Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

*Cal Check ACA*

Form 106-I8000

80-004946  
3/23 YML