



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59281 (03-2014)

Serial Number 80-005938	Instrument Location MOUNTAIN COUNTY JAIL
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000. *Available*
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if Any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>John S. Darnell 1235</i>	Date <i>2-3-16 ✓</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Shanaver</i>	Date <i>08 Mar 2016</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005938
Location = STAN 8164.13.00 06/09
03/03/2016 19:58

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-|=|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 13314080A2
Standard Cyl #? 054
Standard Expiration? 06/05/2016
Oper No? 031123

Flow Cal. Date: 07/15/2013
Slope 664
Intercept -578364

IR Calibration Date: 07/27/2015
 3um 9um

	3um	9um
0th Coef(*100):	-14232	-19463
1st Coef(*100):	289698	139293
2nd Coef(*100):	2380	1071
H2O adj(mg/l*10k):	253	294

***** Printer Test End *****

 1235
Operator Signature
JOHN TRONSET

Remarks:


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005938
Location = STAN 8164.13.00 06/09
03/03/2016 19:59

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:00
02 Std. Gas	0.081	20:00
03 Room Air	0.000	20:01
04 Std. Gas	0.080	20:01
05 Room Air	0.000	20:02
06 Std. Gas	0.080	20:02
07 Room Air	0.000	20:03

Lot No = 13314080A2
Cyl No = 054
Exp Date = 06/05/2016
County = 31 Oper No. = 031123


Operator Signature
JOHN TRONSET

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005938
Location = STAN 8164.13.00 06/09
03/03/2016 20:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:05
02 Reference	RFI*	20:05
03 Room Air	0.000	20:05

*Invalid Test
Inhibited - RFI

Lot No = 13314080A2
Cyl No = 054
Exp Date = 06/05/2016
County = 31 Oper No. = 031123

John *TRONSET* *1235*
Operator Signature
JOHN TRONSET
TEST

Remarks: