



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59281 (5-2011)

Serial Number <i>80-004200</i>	Instrument Location <i>BSPD</i>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P). ✓
 - B. ACA Test (Level 1, Function C). ✓
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air). ✓
- 7. Repair and/or Maintenance Performed (if Any): *NA*
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>Pat Ray</i>	Date <i>3-15-13</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date <i>03 Apr 13</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = BSPD 8164.13.00 06/09
03/15/2013 13:45

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 16512080A2
Standard Cyl #? 36
Standard Expiration? 08/01/2014
Oper No? 029710

Flow Cal. Date: 12/10/2012
Slope 637
Intercept -382487

IR Calibration Date: 12/10/2012
 3um 9um

0th Coef(*100): -14942 -22738
1st Coef(*100): 256196 136921
2nd Coef(*100): 3535 1396
H2O adj(mg/l*10k): 492 515

***** Printer Test End *****



Operator Signature
PATRICK RENZ

Remarks: New Intoxilyzer Installed

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = BSPD 8164.13.00 06/09
03/15/2013 13:22

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:22
02 Std. Gas	0.082	13:23
03 Room Air	0.000	13:23
04 Std. Gas	0.082	13:23
05 Room Air	0.000	13:24
06 Std. Gas	0.082	13:24
07 Room Air	0.000	13:25


Lot No = 16512080A2

Cyl No = 36

Exp Date = 08/01/2014

County = 08

Oper No. = 029710



Operator Signature
PATRICK RENZ

Remarks: New Intoxilyzer Installed
New Cylinder Installed

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = BSPD 8164.13.00 06/09
03/15/2013 13:46

Test	AC	Time
01 Diagnostic	OK	13:49
02 Room Air	RFI*	13:49
03 Room Air	0.000	13:49

*Invalid Test
 Inhibited - RFI

Sub Name = TEST, NEW N
Sub DOB = 01/20/1980
Sub Sex = Male Weight = 185
Test = DUI Cit = 0000019
Dr. Lic. = ND/000990000
Lot No = 16512080A2
Cyl No = 36
Expiration Date = 08/01/2014
County = 08 Oper No. = 029710

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
PATRICK RENZ

Remarks: *New Intoxilyzer Installed*