



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59281 (03-2014)

Serial Number 80-004937	Instrument Location Cass Co. Jail
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P). ✓
 - B. ACA Test (Level 1, Function C). ✓
 - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test). ✓ - Detected (X2)
- 7. Repair and/or Maintenance Performed (if Any): None
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature 	Date 3/23/2016 ✓
Reviewed By (Crime Laboratory Use Only) Deb Shanaver	Date 28 Mar 2016

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004937
Location = CASS 8164.13.00 06/09
03/23/2016 23:43

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 29614080A2
Standard Cyl #? 36
Standard Expiration? 12/05/2016
Oper No? 050910

Flow Cal. Date: 01/18/2010
Slope 730
Intercept -779991

IR Calibration Date: 04/08/2015
 3um 9um

0th Coef(*100): -10967 -18591
1st Coef(*100): 254229 128431
2nd Coef(*100): 2837 1327
H2O adj(mg/l*10k): 348 392

***** Printer Test End *****



Operator Signature
CHAD THOMPSON

Remarks: *Print test*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004937
Location = CASS 8164.13.00 06/09
03/23/2016 23:43

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	23:44
02 Std. Gas	0.081	23:44
03 Room Air	0.000	23:44
04 Std. Gas	0.081	23:45
05 Room Air	0.000	23:45
06 Std. Gas	0.081	23:46
07 Room Air	0.000	23:46

Lot No = 29614080A2
Cyl No = 36
Exp Date = 12/05/2016
County = 09 Oper No. = 050910



Operator Signature
CHAD THOMPSON

Remarks: ACA Test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004937
Location = CASS 8164.13.00 06/09
03/23/2016 23:49

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	23:49
02 Room Air	RFI*	23:49

*Invalid Test
Inhibited - RFI

Lot No = 29614080A2
Cyl No = 36
Exp Date = 12/05/2016
County = 09 Oper No. = 050910

 3524

Operator Signature
CHAD THOMPSON

Remarks:

RFI Test
Detected (X2)

Form 106-I8000