

**INTOXILYZER® 8000 INITIAL INSPECTION**

Intoxilyzer® 8000 Serial Number: 80-006502

**A. GENERAL Setup and Checks:**

- Power On
  - 5 min. countdown during initial warm-up
  - 1 min. countdown to begin test after standby mode
- Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:
  - Level 1 (Operator)    Level 2 (Field Inspector)    Level 3 (Lab Personnel)
- Printer Setup (Level 3,S,P)
- Forms Load (CMI, Inc., CoBrA Software)

**B. LEVEL 1 Setup and Checks:**

- Internal Printer Test and Paper Feed (Level 1,P). Attach record.
- External Printer Test (Level 1,P). Attach record.
3. Barcode Reader (Level 1,B):
  - Operator ID    ND License    R key: Refusal. Attach record.

**C. LEVEL 2 Setup and Checks:**

- Continuous Air Blank (Level 2,A)
- Diagnostic (Level 2,D). Attach record.
- Prelim Data Entry (Level 2,E)
4. Ser Number/Version Info (Level 2,V):  
Inst Ser#: 80-006502 DSP Ver: 1362.14  
Core Ver: 8164.13 CPLD Ver: 7

**D. LEVEL 3 Setup and Checks:**

- DVM Monitor (Level 3,D,D): 3: 12564 s: 1 C: 47.00  
9: 13391 s: 1 B: 45 0  
(RFI indicator with keyed radio): 1
2. Function (F) Keys: :  F5 – Tone                       F9 - Fan  
 F6 – Green LED                       F10 – Pump  
 F7 – Red LED                          F11 – Solenoid
3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 959 CF: 1.05
4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 252 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T): C: 47.0 A: 48.5 I: 0.0  
B: 45.1 V: 0.0
6. Tank Monitor (Level 3,D,G): P (display): 1006 psi P (gauge): 1000 psi
7.  General Setup (Level 3,S,G)
8.  Configure Start Test (Level 3,S,T)
9.  Network (Level 3,S,N)

**E. TESTS (attach all test records):**

**Configure wet bath simulator standards before running tests (Level 1,S):**

1. Zero AC (Level 1,C):  
Sim ser #: MP3067 Lot #: MILLI QH<sub>2</sub>O Exp date: NA  
Known value: 0.000 AC  
Results: 0.000 AC 0.000 AC 0.000 AC
2. Low AC (Level 1,C):  
Sim ser #: MP3066 Lot #: 15160 Exp date: 7.14.17  
Known value: 0.020 AC 0.017  
Results: 0.017 AC ~~0.018~~ AC 0.018 AC
3. High AC (Level 1,C):  
Sim ser #: MP3068 Lot #: 14102 Exp date: 4.23.16  
Known value: 0.400 AC  
Results: 0.406 AC 0.406 AC 0.406 AC

**Configure dry gas standard before running tests (Level 1,S):**

4. Ethanol + 0.05% Acetone (Level 1,B):  
Sim ser #: DR5135 Lot #: ICS4 Exp date: 9.9.17  
Known value: 0.10 AC + 0.05% ACETONE  
 Display reads "Interferent Detect"
5. Acetone (Level 1,B):  
Sim ser #: DR5143 Lot #: ACE2 Exp date: 9.29.17  
Known value: 0.05% ACETONE  
 Display reads "Interferent Detect"
6. Methanol (Level 1,B):  
Sim ser #: DR3377 Lot #: MET2 Exp date: 9.29.17  
Known value: 0.10% METHANOL  
 Display reads "Interferent Detect"
7. Isopropanol (Level 1,B):  
Sim ser #: DR3453 Lot #: IPA2 Exp date: 9.29.17  
Known value: 0.10% IPA CEE ISOPROPANOL  
 Display reads "Interferent Detect"

- 8.  Ambient Fail (Level 1,B): Display and test record read correct.
- 9.  Invalid Sample (Level 1,B): Display and test record read correct.
- 10.  Sequence Aborted (Level 1,B): Display and test record read correct.
- 11.  ABA (Level 1,B): Test record reads correct.
- 12.  Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
- 13.  RFI Detected (B,C, or CMS mode): Display and test record read correct.
- 14.  Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 15.  Deficient Samples (CMS mode): Display and test record read correct.
- 16.  Improper Sample (CMS mode): Display and test record read correct.
- 17.  Difference Too Great (CMS mode): Display and test record read correct.
- 18.  CMS (CMS mode): Test record reads correct.
- 19. Stability Test (Level 3,T,S):

Ave: 0.0822 Std Dev: 0.0006 Rel Std Dev: 0.75%

Maintenance/Repairs/Notes: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

Instrument is acceptable to be used in the field.

Charles E. Edr  
Field Inspector's Signature

24 MAR 16  
Date

Deb Shanaver  
Reviewed by

24 Mar 16  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006502  
Location = TOXL      8164.13.00 06/09  
03/23/2016      13:02

\*\*\*\*\* Printer Test \*\*\*\*\*

12 char/line  
abcdefghijkl  
ABCDEFGHIJKL

16 char per line  
abcdefghijklmnop  
ABCDEFGHIJKLMN

21 chars per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMNQRSTU

24 characters per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMNQRSTU

32 characters per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMNQRSTU

42 characters per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMNQRSTU

Current Instrument Record

Data Entry Mode: Enabled  
Start Test Sequence: DABACABA  
Display Prelim Rslt? Yes  
Display Third Digit? Yes  
Inhib Printer(Y/N)? No  
Display Volume? No  
Disable On Menfull? No  
# of Print Copies? 1  
Select Std (D/W/I)? Dry  
Standard Value? 0.080  
Standard Lot #? 80  
Standard Cyl #? 125  
Standard Expiration? 11/11/2020  
Oper No? 666666

Flow Cal. Date: 08/19/2015  
Slope 674  
Intercept -654491

IR Calibration Date: 08/19/2015  
3um 9um

0th Coef(*100):	-24446	-18103
1st Coef(*100):	267787	135417
2nd Coef(*100):	3387	1610
H2O adj(mg/l*10k):	740	513

\*\*\*\*\* Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks: **PRINT TEST**

Form 106-18000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006502  
Location = TOXL                    8164.13.00 06/09  
03/23/2016                            13:05

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:                    Enabled  
Start Test Sequence:                DABACABA  
Display Prelim Rslt?                Yes  
Display Third Digit?                Yes  
Inhib Printer(Y/N)?                No  
Display Volume?                    No  
Disable On Memfull?                No  
# of Print Copies?                  1  
Select Std (D/W/I)?                Dry  
Standard Value?                    0.080  
Standard Lot #?                    80  
Standard Cyl #?                    125  
Standard Expiration?                11/11/2020  
Oper No?                            666666

Flow Cal. Date:                    08/19/2015  
Slope                                674  
Intercept                          -654491

IR Calibration Date:                08/19/2015  
    3um                                9um

-----  
0th Coef(\*100):                    -24446                                -18103  
1st Coef(\*100):                    267787                                135417  
2nd Coef(\*100):                    3387                                    1610  
H2O adj (mg/l\*10k):                740                                    513

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006502  
Location = TOXL                    8164.13.00 06/09  
03/23/2016                            13:07

Test	AC	Time
01 Room Air	0.000	13:08
02 *Subject Test	REF*	13:08
03 Room Air	0.000	13:08

\*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = REFUSAL  
Dr. Lic. = ND/DIS821456  
Lot No = 29614080A2  
Cyl No = 6  
Expiration Date = 12/05/2016  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

REFUSAL


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006502  
Location = TOXL      8164.13.00 06/09  
03/23/2016      13:09

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

DIAGNOSTIC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006502  
Location = TOXL      8164.13.00 06/09  
03/23/2016      13:22

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:22
02 Std. Sol.	0.000	13:23
03 Room Air	0.000	13:24
04 Std. Sol.	0.000	13:24
05 Room Air	0.000	13:25
06 Std. Sol.	0.000	13:25
07 Room Air	0.000	13:26

08 Sim Temp = 34.0°C

Simul Ser No = MP3067  
Std Sol No = MILLIQ H2O  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: ZERO AC  
0.000 AC



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006502  
Location = TOXL      8164.13.00 06/09  
03/23/2016      13:36

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:37
02 Std. Sol.	0.017	13:37
03 Room Air	0.000	13:38
04 Std. Sol.	0.017	13:39
05 Room Air	0.000	13:39
06 Std. Sol.	0.018	13:40
07 Room Air	0.000	13:40

08 Sim Temp = 34.0°C

Simul Ser No = MP3066  
Std Sol No = 15160  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Low free AC*  
*0.020 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006502  
Location = TOXL      8164.13.00 06/09  
03/23/2016      13:49

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:50
02 Std. Sol.	0.406	13:51
03 Room Air	0.000	13:51
04 Std. Sol.	0.406	13:52
05 Room Air	0.000	13:53
06 Std. Sol.	0.406	13:53
07 Room Air	0.000	13:54

08 Sim Temp = 34.0°C

Simul Ser No = MP3068

Std Sol No = 14102

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

HIGH AC  
0.400 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006502  
Location = TOXL                    8164.13.00 06/09  
03/23/2016                            15:41

Test	AC	Time
01 Room Air	0.000	15:42
02 *Subject Test	INT*	15:42
03 Room Air	0.000	15:43

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 29614080A2  
Cyl No = 6  
Expiration Date = 12/05/2016  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*INTERFERENT*

*0.10 AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006502  
Location = TOXL      8164.13.00 06/09  
03/23/2016      15:45

Test	AC	Time
01 Room Air	0.000	15:46
02 *Subject Test	INT*	15:46
03 Room Air	0.000	15:47

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 29614080A2  
Cyl No = 6  
Expiration Date = 12/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*INTERFERENT*  
*0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006502  
Location = TOXL      8164.13.00 06/09  
03/23/2016      15:47

Test	AC	Time
01 Room Air	0.000	15:48
02 *Subject Test	INT*	15:48
03 Room Air	0.000	15:49

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 29614080A2  
Cyl No = 6  
Expiration Date = 12/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*INTERFERENT*  
*0.10% METHANOL*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006502  
Location = TOXL      8164.13.00 06/09  
03/23/2016      15:52

Test	AC	Time
01 Room Air	0.000	15:52
02 *Subject Test	INT*	15:53
03 Room Air	0.000	15:54

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 29614080A2  
Cyl No = 6  
Expiration Date = 12/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:      0.10% ISOPROPANOL


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006502  
Location = TOXL      8164.13.00 06/09  
03/23/2016      16:14

Test	AC	Time
01 Room Air	AMB*	16:14
02 Room Air	0.000	16:16

\*Invalid Test  
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = AMBIENT FAIL  
Dr. Lic. = ND/DIS821456  
Lot No = 29614080A2  
Cyl No = 6  
Expiration Date = 12/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: **AMBIENT FAIL**

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006502  
Location = TOXL      8164.13.00 06/09  
03/23/2016      16:16

Test	AC	Time
01 Room Air	0.000	16:17
02 Invalid Sample	X.XXX	16:18
03 Room Air	0.000	16:18

\*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INVALID SAMPLE  
Dr. Lic. = ND/DIS821456  
Lot No = 29614080A2  
Cyl No = 6  
Expiration Date = 12/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

INVALID SAMPLE



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006502  
Location = TOXL                    8164.13.00 06/09  
03/23/2016                            16:19

Test	AC	Time
01 Room Air	0.000	16:20
02 *Subject Test	ABT*	16:20
03 Room Air	0.000	16:21

\*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = SEQ ABORTED  
Dr. Lic. = ND/DIS821456  
Lot No = 29614080A2  
Cyl No = 6  
Expiration Date = 12/05/2016  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: SEQUENCE ABORTED

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000           SN 80-006502  
Location = TOXL                   8164.13.00 06/09  
03/23/2016   16:24

Test	AC	Time
01 Room Air	0.000	16:25
02 Subject Test 1	0.000	16:25
03 Room Air	0.000	16:25
04 Reported AC	0.000	16:25

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                   Weight = 150  
Test = DUI                        Cit = ABA TEST  
Dr. Lic. = ND/DIS821456  
Lot No = 29614080A2  
Cyl No = 6  
Expiration Date = 12/05/2016  
County = 08                   Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

ABA TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006502  
Location = TOXL      8164.13.00 06/09  
03/23/2016      16:26

Test	AC	Time
01 Diagnostic	OK	16:27
02 Room Air	0.000	16:27
03 Subject Test 1	0.028	16:28
04 Room Air	0.000	16:30
05 Std. Gas	0.088*	16:31
06 Room Air	0.000	16:32

\*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = CAL OUT OF TOL  
Dr. Lic. = ND/DIS821456  
Lot No = 29614080A2  
Cyl No = 6  
Expiration Date = 12/05/2016  
County = 08      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
OUT OF TOLERANCE

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006502  
Location = TOXL                    8164.13.00 06/09  
03/23/2016                            16:33

Test	AC	Time
01 Diagnostic	OK	16:34
02 Room Air	0.000	16:34
03 *Subject Test	RFI*	16:35
04 Room Air	0.000	16:35

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI DETECT  
Dr. Lic. = ND/DIS821456  
Lot No = 29614080A2  
Cyl No = 6  
Expiration Date = 12/05/2016  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

RFI DETECTED

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006502  
Location = TOXL                    8164.13.00 06/09  
03/23/2016                            16:36

Test	AC	Time
01 Diagnostic	OK	16:37
02 Room Air	0.000	16:38
03 *Subject Test	0.000*	16:41
04 Room Air	0.000	16:43
05 Std. Gas	0.081	16:44
06 Room Air	0.000	16:46
07 *Subject Test	0.000*	16:49
08 Room Air	0.000	16:53

No RFI Detected

\*Deficient Sample - Value Printed was  
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = DEFICIENT SAMP  
Dr. Lic. = ND/DIS821456  
Lot No = 29614080A2  
Cyl No = 6  
Expiration Date = 12/05/2016  
County = 08                            Oper No. = 666666

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the  
instructions displayed by the Intoxilyzer  
in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

DEFICIENT  
SAMPLES

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006502  
Location = TOXL      8164.13.00 06/09  
03/23/2016      16:54

Test	AC	Time
01 Diagnostic	OK	16:55
02 Room Air	IPS*	16:55
03 Room Air	0.000	16:55

\*Improper Sample

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = IMPROPER SAMP  
Dr. Lic. = ND/DIS821456  
Lot No = 29614080A2  
Cyl No = 6  
Expiration Date = 12/05/2016  
County = 08      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*IMPROPER  
SAMPLE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006502  
Location = TOXL                    8164.13.00 06/09  
03/23/2016                            16:56

Test	AC	Time
01 Diagnostic	OK	16:57
02 Room Air	0.000	16:58
03 Subject Test 1	0.048*	16:58
04 Room Air	0.000	17:01
05 Std. Gas	0.082	17:02
06 Room Air	0.000	17:03
07 Subject Test 2	0.077*	17:04
08 Room Air	0.000	17:04

\*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = DIFF TOO GREAT  
Dr. Lic. = ND/DIS821456  
Lot No = 29614080A2  
Cyl No = 6  
Expiration Date = 12/05/2016  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *DIFFERENCE  
TOO GREAT*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006502  
Location = TOXL                    8164.13.00 06/09  
03/23/2016                            17:07

Test	AC	Time
01 Diagnostic	OK	17:08
02 Room Air	0.000	17:09
03 Subject Test 1	0.016	17:09
04 Room Air	0.000	17:11
05 Std. Gas	0.082	17:12
06 Room Air	0.000	17:14
07 Subject Test 2	0.015	17:14
08 Room Air	0.000	17:15
09 Reported AC	0.015	17:14

Difference OK  
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = CMS TEST  
Dr. Lic. = ND/DIS821456  
Lot No = 29614080A2  
Cyl No = 6  
Expiration Date = 12/05/2016  
County = 08                            Oper No. = 666666

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*CMS TEST*



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006502  
Location = TOXL      8164.13.00 06/09  
03/23/2016      17:16

Test	AC	Time
01 Room Air	0.000	17:17
02 Std. Gas	0.083	17:17
03 Room Air	0.000	17:18
04 Std. Gas	0.082	17:18
05 Room Air	0.000	17:19
06 Std. Gas	0.082	17:19
07 Room Air	0.000	17:19
08 Std. Gas	0.082	17:20
09 Room Air	0.000	17:20
10 Std. Gas	0.082	17:21
11 Room Air	0.000	17:21
12 Std. Gas	0.082	17:22
13 Room Air	0.000	17:22
14 Std. Gas	0.081	17:22
15 Room Air	0.000	17:23
16 Std. Gas	0.081	17:23
17 Room Air	0.000	17:24
18 Std. Gas	0.082	17:24
19 Room Air	0.000	17:25
20 Std. Gas	0.082	17:25
21 Room Air	0.000	17:26
22 Std. Gas	0.083	17:26
23 Room Air	0.000	17:26
24 Std. Gas	0.082	17:27
25 Room Air	0.000	17:27
26 Std. Gas	0.082	17:28
27 Room Air	0.000	17:28
28 Std. Gas	0.082	17:29
29 Room Air	0.000	17:29
30 Std. Gas	0.082	17:29
31 Room Air	0.000	17:30
32 Std. Gas	0.082	17:30
33 Room Air	0.000	17:31
34 Std. Gas	0.083	17:31
35 Room Air	0.000	17:32
36 Std. Gas	0.083	17:32
37 Room Air	0.000	17:33
38 Std. Gas	0.083	17:33
39 Room Air	0.000	17:33
40 Std. Gas	0.083	17:34
41 Room Air	0.000	17:34
42 Cal Check Stats:		
43 Average	0.0822	
44 Std Dev	0.0006	
45 Rel Std Dev	0.75%	

STABILITY  
TEST

*Charles Gold*  
24 MAR 16  
pg 1 of 2

Operator Signature  
CHARLES EDER

Remarks: Lot: 29614080A2  
Exp: 12-15<sup>05</sup> 2016  
CYL: 006 cer  
Form 106-I8000

STABILITY  
TEST

*Charles Eder*  
24 MAR 16