

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-005941

Location: ToxL

- Check When Done
- A. Is the breath tube heated? ✓
- B. Is the diagnostic check complete? ✓
- C. Is the time, date, and year correct (re-set if necessary)? ✓
- D. Print test. (Attach test record.) ✓
- E. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN DR5191 Lot No. 12111 Rep AC 0.020 ✓
- F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN DR5189 Lot No. 13300 Rep AC 0.300 ✓
- G. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN DR5135 Lot No. ICS #2 Rep AC 0.100 AC + 0.05% ACETONE ✓
- H. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) ✓
- I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.)

Lot No. 16512080A2 Cyl No. 034 Exp. Date 08/01/2014

Test 1 <u>0.083</u>	Test 4 <u>0.083</u>	Test 7 <u>0.083</u>
Test 2 <u>0.083</u>	Test 5 <u>0.082</u>	Test 8 <u>0.083</u>
Test 3 <u>0.083</u>	Test 6 <u>0.083</u>	Test 9 <u>0.083</u>

Average 0.083

J. Tank pressure, Level 3 Func D Sub G P: 306 psi Reg 300 psi

K. Remarks/Maintenance record: ANNUAL INSPECTION

01 APR 14
Date

[Signature]
Field Inspector's Signature

01 Apr 2014
Date

[Signature]
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005941
Location = TOXL 8164.13.00 06/09
04/01/2014 12:00

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 16512080A2
Standard Cyl #? 34
Standard Expiration? 08/01/2014
Oper No? 666666

Flow Cal. Date: 06/18/2013
Slope 695
Intercept -675557

IR Calibration Date: 06/18/2013

	3um	9um
0th Coef(*100):	-20306	-25137
1st Coef(*100):	291687	138772
2nd Coef(*100):	2450	1058
H2O adj(mg/l*10k):	562	452

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks: *PRINT TEST*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005941
Location = TOXL 8164.13.00 06/09
04/01/2014 13:28

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:29
02 Std. Sol.	0.021	13:30
03 Room Air	0.000	13:30
04 Std. Sol.	0.022	13:31
05 Room Air	0.000	13:31
06 Std. Sol.	0.022	13:32
07 Room Air	0.000	13:33


08 Sim Temp = 34.0°C

Simul Ser No = DR5191

Std Sol No = 12111

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

0.020 AC

Low AC

Exp: 12/04/14

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005941
Location = TOXL 8164.13.00 06/09
04/01/2014 13:35

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:36
02 Std. Sol.	0.305	13:37
03 Room Air	0.000	13:38
04 Std. Sol.	0.307	13:38
05 Room Air	0.000	13:39
06 Std. Sol.	0.308	13:40
07 Room Air	0.000	13:40

08 Sim Temp = 34.0°C

Simul Ser No = DR5189
Std Sol No = 13300
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *HIGH AC*
 0.300 AC
 Exp: 11/04/15
 Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005941
Location = TOXL 8164.13.00 06/09
04/01/2014 13:44

Test	AC	Time
01 Diagnostic	OK	13:45
02 Room Air	RFI*	13:45
03 Room Air	0.000	13:45

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/DIS821456
Lot No = 16512080A2
Cyl No = 34
Expiration Date = 08/01/2014
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: RFI CHECK


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005941
Location = TOXL 8164.13.00 06/09
04/01/2014 13:45

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:46
02 Std. Gas	0.083	13:46
03 Room Air	0.000	13:47
04 Std. Gas	0.083	13:47
05 Room Air	0.000	13:48
06 Std. Gas	0.083	13:48
07 Room Air	0.000	13:49

Lot No = 16512080A2
Cyl No = 34
Exp Date = 08/01/2014
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: 0.080 AC
CALIBRATION CHECK


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005941
Location = TOXL 8164.13.00 06/09
04/01/2014 13:49

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:50
02 Std. Gas	0.083	13:50
03 Room Air	0.000	13:51
04 Std. Gas	0.082	13:51
05 Room Air	0.000	13:51
06 Std. Gas	0.083	13:52
07 Room Air	0.000	13:52

Lot No = 16512080A2
Cyl No = 34
Exp Date = 08/01/2014
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: 0.080AC
CALIBRATION CHECK

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005941
Location = TOXL 8164.13.00 06/09
04/01/2014 13:53

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:53
02 Std. Gas	0.083	13:54
03 Room Air	0.000	13:54
04 Std. Gas	0.083	13:54
05 Room Air	0.000	13:55
06 Std. Gas	0.083	13:55
07 Room Air	0.000	13:56

Lot No = 16512080A2
Cyl No = 34
Exp Date = 08/01/2014
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: 0.080 AC
 CALIBRATION CHECK